



## FLEXLED EXAMINATION LIGHT

With its compact LED design and flexible arm with unlimited range of movement, FlexLED allows accurate positioning of the lamp head above almost any work surface.

FlexLED also delivers a powerfully efficient 40,000 lux @ 1 metre and includes a hygienically sealed arm and head as well as a touch sensitive membrane switch for controls - making cleaning easy.

FlexLED is the ultimate affordable examination light solution.



### TOUCH MEMBRANE SWITCH

#### STAY SAFE

FlexLED's touch membrane switch gives quick access to ON/OFF and 4 stage dimming control. And with its easy-wipe and IP65-rated head, superior infection control has never been easier.

### DESIGNED AND MADE IN AUSTRALIA



FlexLED is designed and manufactured entirely in Australia and built to conform to and exceed rigorous Australian standards and testing.



TROLLEY OPTION  
AVAILABLE

### PRODUCT HIGHLIGHTS

- High output 40,000 lux @ 1 metre
- Flexible arm with unlimited range of movement
- Low power consumption at only 12 Watts Total
- Low-profile, easy-clean head with ON/OFF and dimming control
- Hygienically sealed arm and head
- Durable all metal construction
- 60601 medical compliant power supply
- Quick assembly and plug & play

### FLEXLED EXAMINATION LIGHT



#### APPLICATIONS

- General Practice
- General Examination
- Nursing and Workstations
- First Aid Stations



WALL



DESK




TROLLEY





## UNLIMITED CONTROL


FlexLED's ultra-flexible arm and compact modern LED design makes perfect positioning a breeze.



See our range of movement datasheet for more information.



WALL MOUNT



PLUG & PLAY  
No electrician required



DESK MOUNT



TROLLEY MOUNT

\*All measurements in millimetres (mm)

## FLEXLED EXAMINATION LIGHT

PLANET LIGHTING  
Made in Australia.  
Exported to the World.



PRODUCT CODE	DESCRIPTION	POWER
WM	Wall mount	12W
TM	Trolley mount	12W
DM	Desk mount	12W