




Everyday Impervious Gown

The PrimeOn Everyday Impervious Gown provides high level protection for healthcare workers against blood and bodily fluids. The gown is designed with comfort and safety in mind, featuring an impervious plastic film and full back for an extra level of protection. Suitable for everyday tasks such as patient contact and general cleaning, PrimeOn's Everyday Impervious Gown reduces the potential for cross-contamination.



PrimeOn Everyday Impervious Gown

TGA Registration	ARTG 303904
Material	Lightweight Spunbond fabric with laminate
Gown Length	110.5cm
Impervious	Yes
Cuff Type	Elastic
Latex Free	Yes
Colour	Yellow

REORDER CODE
AGO40UN

FEATURES

- Neck ties
- Full back
- Universal size

PACKAGING

10 gowns per pack
5 packs per carton

STANDARDS

AATCC 42, AATCC 127

MANUFACTURING ACCREDITATIONS
EN ISO 13485

A brand by



Mun Australia
Toll free: 1800 456 837
munglobal.com.au

