



AAMI Level 2 Procedure Gown

Made with a lightweight material, the PrimeOn AAMI Level 2 procedure gown is ideal for healthcare workers performing routine tasks, such as suturing, with a low risk of exposure to bodily fluids. This impervious gown is appropriate for use in public and private hospitals, day surgeries, general practices, dental practices, veterinaries and pathology labs.



Neck
Ties

Full
Back



Knitted
Cuff

PrimeOn AAMI Level 2 Procedure Gown

| | |
|------------------|---------------------------|
| TGA Registration | ARTG 303904 |
| Material | Spunbond PP with laminate |
| Gown Length | 115.0cm \pm 3.8cm |
| Impervious | Yes |
| Cuff Type | Knitted |
| Latex Free | Yes |
| Colour | Blue |

REORDER CODE
AGO52UN

FEATURES

- Neck ties
- Waist ties
- Light and comfortable feel

PACKAGING

10 gowns per bag
10 bags per carton

STANDARDS

AAMI PB70, AATCC 42, AATCC 127

MANUFACTURING ACCREDITATIONS
EN ISO 13485

A brand by



Mun Australia

Toll free: 1800 456 837
munglobal.com.au



Protection Always On