




AAMI Level 4 Procedure Gown

Designed to provide the highest level of barrier protection, the PrimeOn AAMI Level 4 surgical procedure gown is fully impervious and resistant to fluids, viruses and chemotherapy drugs. Recommended for use by healthcare workers who are at the greatest risk of exposure to blood-borne pathogens, this gown is suitable for use in public and private hospitals, day surgeries, private practice and veterinary practices.



PrimeOn AAMI Level 4 Procedure Gown

TGA Registration	ARTG 303904
Material	PP SMS Non-woven fabric with PE film
Gown Length	120.5cm \pm 2.0cm
Impervious	Yes
Cuff Type	Knitted
Latex Free	Yes
Colour	Blue

REORDER CODE AGO54UN

FEATURES

- Chemotherapy drug resistant
- Velcro strap with ties at neck
- Waist ties
- Light and breathable at back

PACKAGING

10 gowns per bag
10 bags per carton

STANDARDS

AAMI PB70, ASTM F1670, ASTM F1671,
ASTM F739

MANUFACTURING ACCREDITATIONS

EN ISO 13485

A brand by



Mun Australia

Toll free: 1800 456 837

munglobal.com.au



Protection Always On