













Examination Gloves Comparative Chart



Product Description		Features				Measurements			Packaging		Standards				
Glove Name	Material	Colour	Texture	Powder Free	Landfill Biodegradable	Length (mm)	Thickness (mm)		Tensile Strength (Mpa) Before Ageing	Box Qty	Carton Qty	Chemo Tested	Food Safe	Biocompatibility Tested	Exam Glove Standards
							Palm	Finger							
AVALON Biodegradable Gloves 	Nitrile	Violet blue 	Fingertip		Yes - up to 30% in 7 months*	Standard Cuff (≥230)	0.07 ± 0.02	0.10 ± 0.02	≥18	200 (XS-L) 180 (XL)	10				ASTM D6319 EN 455 part 1, 2, 3, & 4
CELESTE 	Nitrile	Violet blue 	Fingertip		No	Standard Cuff (≥230)	0.07 ± 0.02	0.09 ± 0.02	≥18	200 (XS-L) 180 (XL)	10				ASTM D6319 EN 455 part 1, 2, 3, & 4

*Based on ongoing ASTM D5526 testing data, glove biodegradation continues beyond 7 months.