



4-Compartment Procedure Tray



The EnvirOn range of medical hollowware is made from sugarcane bagasse, a readily renewable resource and an environmentally friendly alternative to plastics because unlike plastic which take years to break down in landfill, it is compostable and biodegradable.

The EnvirOn 4-Compartment Procedure Tray is an ideal solution for organising and holding a variety of medical essentials thanks to its versatile layout of one large compartment and three smaller separate compartments, which can accommodate a variety of equipment, fluids, medicine as well as instruments to suit different procedures according to their workflow. This tray's thoughtful design not only enhances medical practitioners' efficiency but also contributes to maintaining a hygienic and organised medical environment.

4-Compartment Procedure Tray

Dimensions	266 x 204 x 30mm
Material	Sugarcane Bagasse
Features	Rigid Base, leak & oil proof, puncture resistant
Colour	Natural
Surface Finish	Smooth
Liquid Capacity	Holds liquids up to 90°C for at least 24 hours
Liquid Compatibility	Not compatible with alcohol based liquids

REORDER CODE
EHD14PT

STANDARDS
EN 13432

RECOMMENDED USE
Multipurpose and versatile, may be used as an anaesthetic tray.

PACKAGING
250 units per carton

MANUFACTURING ACCREDITATIONS
ISO 14001 Environment Management

A brand by



Mun Australia
Toll free: 1800 456 837
munglobal.com.au

