




## Artemis Procedure Face Mask

The PrimeOn Artemis Procedure Face Mask provides excellent protection while remaining light and comfortable for everyday use. Resistant against fluid splashes, this latex-free design provides comfort and confidence for the wearer. Elastic earloops ensure the mask is easy to don and remove, yet will stay firmly in place whilst performing your daily tasks.



### PrimeOn Artemis with Earloops

TGA Registration	ARTG 290961
Fixation	Earloops
Number of Fabric Layers	3
Anti-fog Foam	No
Latex Free	Yes
Colour	White
Fluid Resistance, mmHg	120
Bacterial Filtration Efficiency (BFE), %	> 99
Differential Pressure ( $\Delta P$ ), mm H <sub>2</sub> O/cm <sup>2</sup>	< 5.0

**REORDER CODE**  
FPM21UN

**LEVEL BARRIER CLASSIFICATION**  
Level 2

**RECOMMENDED USE**  
Where minimal blood or bodily fluid droplet exposure may occur

**PACKAGING**  
50 masks per box  
6 boxes per carton

**STANDARDS**  
AS4381-2015, ASTM F2100

**MANUFACTURING ACCREDITATIONS**  
ISO 9001, ISO 13485

A brand by



**Mun Australia**  
Toll free: 1800 456 837  
[munglobal.com.au](http://munglobal.com.au)



Protection Always On