




Athena Surgical Face Mask



The PrimeOn Athena Surgical Face Mask provides protection at a premium level for surgical tasks. Our anti-fog foam technology prevents condensation from escaping out of the top of the mask, stopping eyewear from fogging, and ensuring maximum visibility at all times. The PrimeOn Athena surgical face mask is resistant against fluid splashes and conforms to the Australian Standard. The tie design ensures comfort for different head shapes and sizes.



PrimeOn Athena with Ties

TGA Registration	ARTG 290961
Fixation	Ties
Number of Fabric Layers	4
Anti-fog Foam	Yes
Latex Free	Yes
Colour	White
Fluid Resistance, mmHg	160
Bacterial Filtration Efficiency (BFE), %	> 99
Differential Pressure (ΔP), mm H ₂ O/cm ²	< 5.0

REORDER CODE
FPM33UN

LEVEL BARRIER CLASSIFICATION
Level 3

RECOMMENDED USE
Where the risk of exposure to blood or bodily fluid splash are high

PACKAGING
50 masks per box
6 boxes per carton

STANDARDS
AS4381-2015, ASTM F2100

MANUFACTURING ACCREDITATIONS
ISO 9001, ISO 13485

A brand by



Mun Australia
Toll free: 1800 456 837
munglobal.com.au



Protection Always On