



Athena Procedure Face Mask

The PrimeOn Athena Procedure Face Mask provides protection at a premium level for many procedures. Our anti-fog foam technology prevents condensation from escaping out the top of the mask, preventing eyeglasses from fogging, and ensuring maximum visibility at all times.

The PrimeOn Athena Procedure Face Mask is resistant against fluid splashes and meets the Australian Standard for the fluid resistant requirement. The earloop design ensures this mask can be easily donned and doffed while providing a secure fit for all head shapes and sizes.




Fluid Repellent Filter
Fluid repellent non-woven layer

Outer Layer
Additional fluid repellent material

Inner Layer
Fluid repellent to prevent bacterial growth

Middle Layer
Meltblown filter



PrimeOn Athena with Earloops

TGA Registration	ARTG 290961
Fixation	Earloops
Number of Fabric Layers	4
Anti-fog Foam	Yes
Latex Free	Yes
Colour	Blue
Fluid Resistance, mmHg	160
Bacterial Filtration Efficiency (BFE), %	> 99
Differential Pressure (ΔP), mm H ₂ O/cm ²	< 5.0

REORDER CODE
FPM34UN

LEVEL BARRIER CLASSIFICATION
Level 3

RECOMMENDED USE
Where the risk of exposure to blood or bodily fluid splash are high

PACKAGING
50 masks per box
6 boxes per carton

STANDARDS
AS4381-2015

MANUFACTURING ACCREDITATIONS
ISO 9001, ISO 13485

A brand by



Mun Australia
Toll free: 1800 456 837
munglobal.com.au

