

gloveon Innova



Latex Exam Gloves Powder Free, Standard Cuff

GloveOn® Innova is our light weight powder free latex glove, with a unique fully textured surface for better grip in both wet and dry environments. These versatile gloves come with high tensile strength and superb barrier protection properties, which are suitable for any demanding task in any market segment.



		GloveOn® Innova	
Length (mm)		≥ 230	
Thickness Measurements (mm)			
Palm (centre of Palm)		0.11 ± 0.03	
Finger (13mm ± 3mm from tip)		0.14 ± 0.03	
Physical Properties	Before Ageing	After Ageing	
Tensile Strength (MPa)	≥ 18	≥ 16	
Elongation (%)	≥ 650	≥ 500	
Inspection Levels & AQL	Inspection Level	AQL	
Watertightness	G1	1.5	
Physical Dimensions	S2	4.0	
Tensile Strength	S2	4.0	
Visual Inspection (Major)	S4	2.5	
Visual Inspection (Minor)	S4	4.0	
Particulate Residue	N = 5	≤ 2mg/glove	
Protein Content	N = 3	≤ 50µg/g	

REORDER CODE

INV22XS	X-SMALL
INV22SS	SMALL
INV22MM	MEDIUM
INV22LL	LARGE
INV22XL	X-LARGE

FEATURES

- Fully textured
- Powder free
- Made with natural rubber latex
- Protein level ≤ 50µg/g
- Ambidextrous
- Standard cuff
- Natural colour

PACKAGING

100 gloves per box
10 boxes per carton

REGULATORY COMPLIANCE

ARTG 127794, FDA 510(k),
MDD 93/42/EEC, REACH, EU 10/2011,
EC 1935/2004, EU 2016/425

STANDARDS

AS/NZS 4011-1, ASTM D5151,
ASTM D6124, ASTM D3578, ASTM D5712,
EN ISO 341-1 (Type B), EN 16523-1,
EN 374 part 2, 4 & 5, EN 420, EN 388,
EN 1186, EN 13130, CEN/TS 14234,
ISO 10993 part 5 & 10, HACCP Certified

MANUFACTURING ACCREDITATIONS

ISO 9001, ISO 13485, EN ISO 13485

*Full test reports on Sensitisation Testing, Primary Skin Irritation and Residual Proteins Testing are available upon request.

A brand by



Mun Australia
Toll free: 1800 456 837
munglobal.com.au



Proud Supporter of
Birthing Kit Foundation (Australia)



Protection Always On

MKG2019-051