

gloveon Hammer



Nitrile Exam Gloves Powder Free, Standard Cuff

GloveOn® Hammer is the nitrile glove designed for the tough jobs that demand high tensile strength, elongation and a definitive colour. These gloves are suitable for security guards, tattoo artists, beauty stylists, motor mechanics, cleaners and general maintenance professionals.



GloveOn® Hammer		
Length (mm)	≥ 230	
Thickness Measurements (mm)		
Palm (centre of Palm)	0.07 ± 0.02	
Finger (13mm ± 3mm from tip)	0.11 ± 0.02	
Physical Properties	Before Ageing	After Ageing
Tensile Strength (MPa)	≥ 18	≥ 16
Elongation (%)	≥ 500	≥ 500
Elongation after Isopropanol immersion (%)	≥ 500	N/A
Inspection Levels & AQL	Inspection Level	AQL
Watertightness	G1	1.5
Physical Dimensions	S2	4.0
Tensile Strength	S2	4.0
Visual Inspection (Major)	S4	2.5
Visual Inspection (Minor)	S4	4.0
Particulate Residue	N = 5	≤ 2mg/glove

REORDER CODE

NTR23XS	X-SMALL
NTR23SS	SMALL
NTR23MM	MEDIUM
NTR23LL	LARGE
NTR23XL	X-LARGE

FEATURES

- Fingertip textured
- Powder free
- Not made with natural rubber latex
- UL certified to meet NFPA 1999
- Ambidextrous
- Standard cuff
- Black colour

PACKAGING

100 gloves per box for XS to L
90 gloves per box for XL
10 boxes per carton

REGULATORY COMPLIANCE

ARTG 300229, FDA 510(k),
REACH, ROHS Directive 2002/95/EC,
89/686/EEC

STANDARDS

EN 374 part 1, 2 & 3, EN 420, EN 388,
EN 455 part 1, 2 & 3, EN 16523-1,
HACCP Certified

MANUFACTURING ACCREDITATIONS

ISO 9001, ISO 13485, EN ISO 13485

A brand by



Mun Australia
Toll free: 1800 456 837
munglobal.com.au



Proud Supporter of
Birthing Kit Foundation (Australia)



Protection Always On

MKG2019-064