



gloveon

Victor

Polyisoprene Surgical Gloves Powder Free

Our GloveOn® Victor polyisoprene surgical glove provides surgeons with a fine touch for surgical procedures requiring a high tactile response. Combined with great dexterity, comfort and feel, GloveOn® Victor allows you to concentrate fully on the job at hand. Best used when worn in conjunction with GloveOn® Victor Under Glove for a complete double donning solution that aids in puncture detection.



GloveOn® Victor

Physical Dimensions				
Length (mm)	290 ± 10			
Size	5.5 6.0 6.5 7.0 7.5 8.0 8.5 9			
Palm width ± 4 mm	72 78 84 89 96 103 110 116			
Thickness Measurements (mm)				
Cuff (25 mm from cuff edge)	0.15 ± 0.04			
Palm (centre of palm)	0.19 ± 0.04			
Finger (13 mm ± 3 mm from tip)	0.21 ± 0.04			
Physical Properties	Before Aging	After Aging		
	AS/NZS 4179	Typical Value	AS/NZS 4179	Typical Value
Force at break (N)	≥ 12.5	16	≥ 9.5	14
Elongation at break (%)	≥ 700	1053	≥ 550	1069
Modulus at 300% elongation (N)	≤ 2.0	1.2	-	-
Inspection Levels & AQL		Inspection Level	AQL	
Watertightness		G1	0.65	
Physical dimensions		S2	4.0	
Physical properties		S2	4.0	
Visual inspection (major)		S4	2.5	
Visual inspection (minor)		S4	4.0	
Particulate residue		N = 5	≤ 2 mg/glove	
Seal strength		N = 20	≥ 1.2 N	
Dye penetration test		N = 20	Accept 0 Reject 1	

REORDER CODE

SPI2855	Size 5.5
SPI2860	Size 6.0
SPI2865	Size 6.5
SPI2870	Size 7.0
SPI2875	Size 7.5
SPI2880	Size 8.0
SPI2885	Size 8.5
SPI2890	Size 9.0

FEATURES

- Powder free
- Not made with natural rubber latex
- Textured
- Curved fingers
- Hand specific
- Natural white colour
- Sterilised using gamma irradiation

PACKAGING

1 pair per pouch
50 pairs per box
4 boxes per carton

REGULATORY COMPLIANCE

ARTG 291855, FDA 510 (k),
MDD 93/42/EEC

STANDARDS

AS/NZS 4179, ASTM D3577, ASTM D5151,
ASTM D6124, ASTM D6978, ASTM D7102,
ASTM F1671, EN 455 part 1, 2 & 3, EN 16523-1,
ISO 10282, ISO 10993 part 5 & 10

MANUFACTURING ACCREDITATIONS

ISO 9001, ISO 13485, EN ISO 13485

A brand by



Mun Australia
Toll free: 1800 456 837
munglobal.com.au

