

gloveon VinClear

Vinyl Exam Gloves Powder Free, Standard Cuff

GloveOn® VinClear is a lightweight, powder free vinyl glove made from 100% PVC, eliminating the threat of latex allergy for the user. GloveOn® VinClear gloves are the right fit for keeping your hands clean in a variety of light duty tasks.



	GloveOn VinClear	
Length (mm)		
	≥ 230	
Thickness Measurements (mm)		
Palm (centre of Palm)	0.08 ± 0.03	
Finger (13mm ± 3mm from tip)	0.08 ± 0.03	
Physical Properties	Before Ageing	After Ageing
Tensile Strength (MPa)	≥ 15	≥ 15
Elongation (%)	≥ 400	≥ 400
Inspection Levels & AQL	Inspection Level	AQL
Watertightness	G1	1.50
Physical Dimensions	S2	4.00
Tensile Strength	S2	4.00
Visual Inspection (Major)	G1	2.50
Visual Inspection (Minor)	G1	4.00
Particulate Residue	S2	≤ 2mg/glove

REORDER CODE

VNY11XS	X-SMALL
VNY11SS	SMALL
VNY11MM	MEDIUM
VNY11LL	LARGE
VNY11XL	X-LARGE

FEATURES

- Powder Free
- Not made with natural rubber latex
- Ambidextrous
- Standard cuff
- Transparent

PACKAGING

100 gloves per box
10 boxes per carton

REGULATORY COMPLIANCE

ARTG 100687

STANDARDS

AS/NZS 4011.2, EN 455 part 1, 2 & 3
EN 1186, EN 13130, HACCP Certified

MANUFACTURING ACCREDITATIONS

ISO 9001, ISO 13485, EN ISO 13485

A brand by



Mun Australia
Toll free: 1800 456 837
munglobal.com.au



Proud Supporter of
Birthing Kit Foundation (Australia)



Protection Always On