

NULINE

REFRIGERATION

A SIGN OF QUALITY

HR Pharmacy



Nuline Refrigeration's HR Series Pharmaceutical Refrigerators have been tested and approved by The Pharmacy Guild Of Australia to meet every requirement and to ensure peace of mind. All refrigerators are fitted with a controller that enables you to check your temperature min/max at any stage, simply by pushing a button. There are 3 sizes available from 135 litres to 600 litres capacity.

Features & Benefits

- QCPP approved
- Digital temperature indicator
- Built in alarm
- Door lock

Optional

- Glass/Solid door
- Data logger
- Chart recorder
- Voltage free connection



Model	Dimensions	Capacity	Shelves	Door
HR200	600W x 620D x 850H	135 Litre	3 Plus Base	Solid
HR400	600W x 620D x 1850H	350 Litre	5 Plus Base	Solid
HR600	770W x 750D x 1900H	570 Litre	5 Plus Base	Solid