



20 Liquid Nitrogen Container for Storage

Feature:

- *More than five year's life
- *Made of high-strength aluminum alloy.
- *Adopting advanced powder coating technique
- *A leather (PU) bag is available for models of the capacity less than 50L.

Technical Parameters:

Model		Unit	LNC-20-50	
Capacity of LN2		L	20	
Weight empty		kg	9.5	
Caliber(neck opening of Dewar)		mm	50	
Outer diameter of Dewar		mm	394	
Height of Dewar		mm	672	
Static evaporation rate		L/d	0.12	
Static holding time		days	168	
Canister	Canister outer diameter	mm	38	
	Canister height	mm	120/276	
	Canisters No. (pc)	pc	3~ 6	
	Canisters No. (standard) (pc)	pc	1	
Straws	holding capacity (single 120mm)	0.5ml	pc	792
		0.25ml	pc	1788
	Straw holding capacity (double 276mm)	0.5ml	pc	1284
		0.25ml	pc	2832