



# NEW PRODUCT

## PREMIER NITROSPRAY PLUS CRYOSURGICAL UNIT

Premier Nitrospray Plus Cryosurgical Unit constructed of a double-walled vacuum stainless steel container (canister) to hold liquid nitrogen. The outflow track extension port, located on the top of the instrument has a luer-lock connector to easily engage/disengage spray tips and other attachments as required.

### Intended Use

To provide a means of transporting liquid nitrogen to the patient and to dispense the liquid nitrogen as a cold gas in a controlled manner during general and dermal cryosurgical procedures. The gaseous nitrogen, when dispensed is at a temperature of  $-196^{\circ}\text{C}$ .

**Indicated for treatment of:** Benign and premalignant lesions and some selected malignant lesions

### Warranty:

Nitrospray Unit warranted for three (3) years from date of purchase against defects in workmanship or parts. Warranty is void if the unit is subjected to improper use or damage and unauthorised repairs.

**Package Inclusions:** Flask (300 or 475ml), Set of 5 spray tips and spray tip protector



### Order Codes:

#### Devices

LNSLITEPLUS Nitrospray Lite Plus 300ml

LNSPLUS Nitrospray Plus 475ml

#### Accessories

LNSTSET Tip Set 16, 17, 18, 19 & 20 Gauge

LNST16G Tip 16 Gauge

LNST17G Tip 17 Gauge

LNST18G Tip 18 Gauge

LNST19G Tip 19 Gauge

LNST20G Tip 20 Gauge

LNSTP Spray Tip Protector

LNSABA Accessory Adaptor For Screw Tips

LNSNER Needle Extender Rigid

LNSTEB Spray Tip Extension Bent

### Specifications:

Holding Times	
Static	6 – 8 hours
Intermittent Use	3 – 6 hours of patient treatment time
Capacity	
Nitrospray Lite	300ml
Nitrospray Plus	475ml
Materials	
303 Stainless Steel – ASTM A582 304 Stainless Steel – ASTM A269 ACETAL – ASTM D6778	

## Pacific Medical

Phone: 1300 076 332

Email: [info@pacificmedical.com.au](mailto:info@pacificmedical.com.au)

[www.pacificmedical.com.au](http://www.pacificmedical.com.au)



#### Team Medical Supplies

Phone 1300 22 44 50

Fax 1300 22 44 60

[info@teammed.com.au](mailto:info@teammed.com.au)

[www.teammed.com.au](http://www.teammed.com.au)



#### BP Medical (WA)

Phone 1300 33 55 60

Fax 1300 33 55 70

[info@bpmedical.com.au](mailto:info@bpmedical.com.au)

[www.bpmedical.com.au](http://www.bpmedical.com.au)