



**PLATINUM
HEALTH**
SUPPLIES FOR LIFE



PRODUCT SPECIFICATION SHEET

AUSTRALIAN MADE N95/P2 RESPIRATORS

SKU: SGRM-002

TGA LISTING: 341917



APPLICABLE DIRECTIVE & PRODUCT STANDARDS:

This product is approved for use in Healthcare, Medical, Civil and Industrial environments to provide protection for half the face.

PRODUCT DESCRIPTION & MATERIALS:

White, flat folded particular respirator and surgical mask with ear loops.
The masks are latex-free and disposable.

SPECIFICATIONS:

- EN 149 Certified
- AS1716 Certified
- EN149 FFP2 NR (granted)
- NIOSH 42C FR84 N95 Certified
- BFE 98% and higher
- PFE 95% and higher

PACKING:

MASKS PER BOX - 10
BOXES PER CARTON -
CARTONS PER CONTAINER -

SUGGESTED STORAGE CONDITIONS AND EXPIRY:

Please store in a cool, dry and well-ventilated room. Keep away from fire, dust exposure and direct sunlight.
Shelf life of 2 years after the month of manufacturing.



Team Medical Supplies

Phone 1300 22 44 50

Fax 1300 22 44 60

info@teammed.com.au

www.teammed.com.au



BP Medical (WA)

Phone 1300 33 55 60

Fax 1300 33 55 70

info@bpmedical.com.au

www.bpmedical.com.au

