



AQ Medical 46

VACCINE REFRIGERATORS

KEY FEATURES

- QCPP cold chain accredited
- Digital temperature control system preconfigured to hold temperature from +2°C to +8°C in line with 'Strive for 5' storage guidelines
- Lockable door
- Adjustable shelves
- Low energy consumption, low noise output
- Audible and visual auto resetting alarms
- Dedicated technical support line

OPTIONAL EXTRAS

- UPS & Battery Backup
- Controller Backup Battery (temperature logging)
- AQ Box cloud temperature monitoring system
- Comprehensive programmed maintenance

PRODUCT SPECIFICATIONS

Dimensions	Height 543 mm
	Width 430 mm
	Depth 475 mm
Gross Weight	22 kg
Gross Capacity	41 L
Temperature Range	+2°C to +8°C
Number of shelves	1.5
Power Consumption/24 hrs	0.12 kWh
Refrigerant	R600a
Warranty (parts & labour)	24 months

