

Medisafe Plus 80

VACCINE REFRIGERATORS



KEY FEATURES

- QCPP cold chain accredited
- Dual temperature probes for both cabinet temperature and simulated vaccine temperature at +2°C to +8°C
- Integrated door lock(s)
- Fully integrated AQ Box provides IoT cloud-based monitoring system with customizable email and text alarm notifications
- Audible and visual alarms for: temperature, door open, power failure

PRODUCT SPECIFICATIONS

| | |
|----------------------------|--|
| Dimensions | Height: 830 mm Width: 500 mm Depth: 560 mm |
| Gross Weight: | 43 kg |
| Gross Capacity: | 80 L |
| Doors: | 1 |
| Temperature Range: | +2°C to +8°C |
| Number of shelves: | 3 |
| Power Consumption/24 hrs: | 1.05 kWh |
| Refrigerant: | R600a |
| Warranty (parts & labour): | 5 year |

OPTIONAL EXTRAS

- UPS Battery Backup
- Controller Backup Battery (temperature logging)
- Customisable controller parameters
- Laser printed branding (bulk orders)
- Comprehensive programmed maintenance

