



## Medisafe Plus VACCINE REFRIGERATORS



Medisafe Plus are premium vaccine storage units that come complete with the built in Medisafe Plus cloud monitoring system. Known as the AQ Box, the Medisafe Plus cloud monitoring system can provide access to live data via any web-enabled device such as a mobile phone, tablet, computer or laptop. If there is a temperature breach or power outage, it sends alerts via SMS and email.

Medisafe Plus units are uniquely equipped with two temperature probes for rigorous monitoring of both the cabinet temperature and a simulated vaccine temperature, yielding the most accurate measurements possible. This is a unique feature in the market of vaccine protection units. In addition, Medisafe Plus units are built with a Bluetooth enabled controller that provides additional data collection capacity. Working with the AQ Box portal, users can set up personalised alerts and reporting.

All Medisafe Plus units are purpose-built for vaccine storage in line with the Australian Government Department of Health's 'Strive for 5' vaccine storage guidelines and are QCPP accredited, giving you confidence that your vaccines will be stored securely. Medisafe Plus customers have access to our fully dedicated and local support team and national network of service technicians for all their support needs.