

Medisafe Plus 420

VACCINE REFRIGERATORS



KEY FEATURES

- QCPP cold chain accredited
- Cutting edge new controller with Bluetooth capability, holding over a week of temperature and alarm logging
- Dual temperature probes for both cabinet temperature and simulated vaccine temperature
- Integrated door lock(s)
- 6 full size shelves and 1 half shelf
- Fully integrated AQ Box provides IoT cloud-based monitoring system with customizable email and text alarm notifications
- Audible and visual alarms for: temperature, door open, power failure

PRODUCT SPECIFICATIONS

Dimensions	Height: 2200mm Width: 600 mm Depth: 700 mm
Gross Weight:	92 kg
Gross Capacity:	422 L
Doors:	1
Temperature Range:	+2°C to +8°C
Number of shelves:	6
Power Consumption/24 hrs:	2.46Wh
Refrigerant:	R290
Warranty (parts & labour):	5 year

OPTIONAL EXTRAS

- UPS Battery Backup
- Controller Backup Battery (temperature logging)
- Customisable controller parameters
- Laser printed branding (bulk orders)
- Comprehensive programmed maintenance

