

Medisafe Plus 1000

VACCINE REFRIGERATORS



KEY FEATURES

- QCPP cold chain accredited
- Low energy consumption, low noise output
- Dual temperature probes for both cabinet temperature and simulated vaccine temperature
- Integrated door lock(s)
- Adjustable shelves
- Fully integrated AQ Box provides IoT cloud-based monitoring system with customizable email and text alarm notifications
- Audible and visual alarms for: temperature, door open, power failure
- Dedicated technical support line

OPTIONAL EXTRAS

- UPS Battery Backup
- Controller Backup Battery (temperature logging)
- Customisable controller parameters
- Laser printed branding (bulk orders)
- Comprehensive programmed maintenance

PRODUCT SPECIFICATIONS

Dimensions	Height: 2200mm Width: 1130mm Depth: 745mm
Gross Weight:	230 kg
Gross Capacity:	980L
Doors:	2
Temperature Range:	+2°C to +8°C
Number of shelves:	10
Power Consumption/24 hrs:	3.79kWh
Refrigerant:	R290
Warranty (parts & labour):	5 year

