

Veinlite LED+[®]

The Ultimate Tool for Clinical and Dental Vein Access Procedures

Veinlite LED+ uses clinically proven transillumination technology to reveal underlying venous structures, including veins, valves, and bifurcations up to 6mm deep. Its patented C-shaped design assists in venous access by gently providing traction, acting as a tourniquet, and preventing veins from rolling during vein access.



Why choose Veinlite LED+

- Clinically proven Veinlite technology*
- Fewer IV access failures
- Improved patient satisfaction
- Works on patients of all ages and skin tones
- Meets INS and CDC Guidelines for infection control
- FDA and CE registered



www.veinlite.com
info@veinlite.com
Tel: 281.240.3111

Follow us



Veinlite[®], Veinlite LED+[®], and the Veinlite logo are registered trademarks of TransLite LLC

*Patent Numbers: 7,874,698; 7,841,751; 8,032,205; 8,177,808