

From Australia's leading supplier of magnifiers, medical and specialty lighting...

# Light Up your Surgery

with precision LED technology from  
**Derungs** 

- ✓ Intuitive
- ✓ Reliable
- ✓ High-performance

*NEW  
Ceiling  
Magnifier Lamp*



**Clarity of vision** and **adequate illumination** are crucial for the **day-to-day procedures** in any successful surgery or clinic. The Swiss-based medical lighting specialists **Derungs** have made their mark on the global medical lighting industry by stepping up to the challenge.

This **premier range of medical lights** feature sealed constructions for **easy disinfection**, **intuitive** and **ergonomically-tested** design for daily operation, **generous lighting areas** and **sturdy** material choices.

These **high-performance illuminators** are adaptable to any clinic or surgery, with a **range of sizes**, optional **wall, ceiling or floor mounts**, and a wide variety of **colour temperatures** and **brightness levels** to select from.

Proudly distributed in Australia by

**RedbankGroup**  
For a clearer tomorrow

1300 788 239 [redbank.net.au](http://redbank.net.au)

# Illumination in Detail

## The OpticLux LED Magnifier (with UV option)

The award-winning OpticLux is designed with everything a modern magnifier luminaire is meant to have: a robust 1.9x magnifier, the latest, maintenance-free LED modules, a huge 1200mm reach with ergonomically-tested arm joints and a highly intuitive operation.

Its daylight-quality light field and outstanding colour rendering properties are of great value in dermatology, human medicine, dentistry, veterinary medicine, in laboratories and for monitoring sterilized products.

- Bright, homogeneous light field of 6000 lux at 0.15m
- Large viewing field of Ø153mm
- Colour temperature daylight white 6500 K
- Acrylic lens with 3.5 dioptries provides 1.9x magnification
- Continuously dimmable
- Colour rendering properties CRI → 93, R9 → 90, R13 → 95
- Optional Wood light model at 365 nm



OpticLux Ceiling Magnifier  
#ML 164

UV Model  
#ML 166

OpticLux Wall Magnifier  
#ML 168-W

UV Model  
#ML 167-W

OpticLux Mobile Magnifier  
#ML 165-S

UV Model  
#ML 167-S

# Our Best & Brightest

## The Visiano LED Exam Light

Give your medical staff unparalleled clarity in minor procedures and general examinations with the Visiano LED exam light by Derungs.

It offers the best in LED medical lighting technology, combining low heat output with adjustable colour and brightness to adapt smoothly to any situation.

- Delivers a powerful 60,000 lux at 0.5m
- Large light field diameter of Ø210mm
- Brings out the most subtle details with changeable colour temperature (4500K cool /3500K warm)
- 4 dimming levels
- Total freedom of movement with 360° rotatable light head and integrated handles
- Easily cleanable with disinfectant spray - fully enclosed light head rated to IP20



Visiano Wall Examination Light  
#ML 515-W

Visiano Mobile Examination Light  
#ML 515-S

Visiano Ceiling Examination Light  
#ML 514





Derungs 

## The Everyday Performer

*The Visiano 10 LED Exam Light*

The Visiano 10 is an LED examination light for everyday use, with excellent illumination, a large light field, and a durable, sealed construction.

A sturdy, flexible arm ensures positional accuracy in every direction.

- High illumination strength of 50,000 lux at 0.5m
- Generous light field of Ø180 mm
- Colour temperature of 4400 K
- Unique diamond optics and maintenance-free premium LEDs provide excellent light quality
- 280° rotation of the luminaire head, high-quality aluminium linkage and durable friction-joints
- Slimline design with ergonomic handle
- Precise colour rendering Ra →93



Visiano 10 Mobile Light  
#ML 501-S



Visiano 10 Wall Light  
#ML 501-W

## Versatile and Adaptable

*The Halux LED Exam Light*

An innovation in the field of illumination, the Halux colour-mixes different LEDs to provide a range of colour temperatures.

Along with the 5 dimming levels, the user has the ability – depending on the skin, tissue or wound – to select the optimal inspection light via an intuitive control unit at the lamp head.

- Delivers 40,000 Lux at 0.5m
- Light field of Ø180mm
- 3 adjustable colour temperatures (3300K, 3800K and 4400K)
- 5 dimming levels from 100% to 20%
- 950mm reach with spring-balanced arm system for precise light positioning



Halux 20-3 Mobile Light  
#ML 509-S

Halux 20-3 Wall Light  
#ML 509-W



## Pure Brilliance

### *The Iris LED Treatment Light*

The Derungs Iris LED treatment light provides high-performance illumination in surgeries, clinics and hospitals around Australia.

The Swiss-made Iris is designed to be highly ergonomic, intuitive and adjustable, providing total control and illumination when it's needed most.

- Provides 50,000 lux at 1m
- Generous light field of Ø170mm
- Includes correction filter for 4300K white light
- 16watt LED, with a dual output intensity of -60% or 100%
- Unique reflector design guarantees homogeneous illumination
- Weight-balanced articulation joint system for easy adjusting and stable positioning
- Easily disinfected with smooth surfaces, resistant materials and removable, sterilisable handle

Proudly distributed in Australia by

**RedbankGroup**  
*For a clearer tomorrow*

1300 788 239 [redbank.net.au](http://redbank.net.au)

