

NEW

MaggyScan

Portable Maggylamp



A New Horizon

MaggyScan Portable Maggylamp


We've heard from medical professionals across the country, and have joined forces with Australian designers and manufacturers to bring you the MaggyScan Portable Maggylamp.

With a cordless operation, user-tested intuitive design, and 2.5 hour run time, this is the perfect inspection tool to supplement your everyday medical procedures.


- No more tangled cords
- Recharges as soon as placed in wall dock
- 2.5 hour continuous run time
- Convenient wall charging dock
- Brighter than the current LED Maggylamp, 6000lux @ 25cm
- 125mm glass lens, 2x magnification
- Wall mounting fixtures included

#ML 219 MaggyScan Portable Maggylamp

Well-protected, reusable packaging



Ideal for skin inspection



Manufactured in Australia by

RedbankGroup
For a clearer tomorrow

1300 788 239
redbank.net.au

Supplied by

team medical
supplies

Team Medical | www.teammed.com.au | info@teammed.com.au | phone 1300 22 44 50