


## Treatment luminaire

TRIANGO 100-1 C D15884000-00802695

**Part No. ML542**



<b>fitted with</b>	LED - Energy efficiency category A+
<b>Illuminancefield @ 100 cm</b>	Ø 18 cm
<b>illuminance @ 100 cm</b>	100 000 lx
<b>color rendering index (CRI)</b>	> 95
<b>colour temperatur</b>	4300 K
<b>work equipment</b>	electronic converter
<b>mains lead</b>	direct power supply
<b>luminaire body</b>	polyamide (PA)
<b>weight (net)</b>	15.9 kg
<b>Power consumption</b>	40 W
<b>usage</b>	keypad at luminaire head
<b>balance of articulation joints</b>	spring
<b>class of protection</b>	I
<b>Accessory</b>	sterilizable handle, adapter for suspended ceiling
<b>norms</b>	EN 60601-1, EN 60601-2-41
<b>special features</b>	4 dimming levels, Diamond optics, High-quality arm system, sterilisable grip, maintenance free
<b>Fastening</b>	ceiling mount, when mounting ceiling tube could be cut to length 25-120 cm
<b>System of protection</b>	IP 20 (43)
<b>Part number</b>	ML542

### light distribution

