

# Actimove® Umerus

## Shoulder Immobiliser



- Improves shoulder stabilisation through the waist strap which can be adjusted to suit most patients
- Allows for better hand positioning thanks to the thumb loop which limits the rotation of the forearm while still allowing full finger movement
- Simple to apply and adjust for an individual fit with easy threading through the freely rotating O-rings
- Machine washable at 40°C
- Radiolucent
- Metal-free

Actimove® Umerus is the cost-effective shoulder and arm immobiliser for limited abduction, extension, internal and external rotation of the humerus.

It helps to keep your inventory low, as it fits on both left and right arm, available in four sizes for your patients who will appreciate its easy application.



### Indications

- Post operative support
- Fracture management (e.g. clavicle, scapular, humeral, and greater tuberosity fractures)
- Joint dislocations (e.g. acromioclavicular & sternoclavicular dislocations, shoulder luxation)
- Inflammatory diseases (e.g. calcific tendinitis of the rotator cuff)
- Minor instability of the shoulder
- Support for hemiplegic patients

**LATEX  
FREE  
FORMULA**

### Ordering information Sizing Chart



Ref. No.	Size	Measurements	Items per Unit
72819-65	S	< 37 cm	1 immobiliser
72819-66	M	37 – 40 cm	1 immobiliser
72819-67	L	40 – 43 cm	1 immobiliser
72819-68	XL	> 43 cm	1 immobiliser

Orthopaedic Soft Goods

### Contents

Non-invasive fracture treatment pathway

Primary Immobilisation

Secondary Immobilisation

Undercast Materials

Cast Removal Equipment

Cast Protection

Orthopaedic Soft Goods