



**Aqualife
Products®**

**EXCELLENCE IN WATER TREATMENT
& LOW COST MAINTENANCE SYSTEMS**



**Quality Assurance
Guaranteed**

NEW

AquaClave M900 Series

YOUR PRACTICE NEEDS

Engineered for thermal washers

AND

Steam Sterilisation equipment

The NEW AquaClave Mark 3 900 Series

remains the proven 5 stage Aqualife RO system but with an enhanced production capacity. Sufficient water to feed thermal washer disinfection and steam sterilisation equipment.

The **M900** includes a larger capacity membrane, double ended high flow pump, a 60 litre storage tank or with 80 litre tank option plus all the enhanced features of the standard AquaClave system.

Check your practice equipment needs with your authorised and approved Aqualife agent.



Reverse Osmosis -

The AquaClave M900 series model has a *high capacity* membrane with excellent abrasive resistance. The membrane will remove approximately 95% of all dissolved minerals and salts.

The two stage polishing system feeding the AutoClave will remove remaining organic and non-organic residues thereby finishing the water quality to meet all autoclave operating standards. The water for the thermal washes by-passes this final treatment and is supplied direct from the membrane.

Commercial Advantages -

The availability of this volume of water enhances the economics of your investment and makes your choice of an **AquaClave M900 series** an intelligent decision.

Dimensions:

(System) 46cm W x 54cm H x 20cm D
(60L Tank) Ø40cm x 60cm H

Please note that this system is only suitable for use on disinfected potable water supplies from your local water authority. For areas where tank or other water is required please consult Aqualife. To achieve the water quality indicated above the system must be operated in accordance with the operating and maintenance instructions. On site verification of the water quality is required for non-municipal water supplies.

Aqualife Products (Australia Pty Ltd) Freecall 1800 623 887

**Phone 03 9310 5259 | Facsimile 03 9310 5387 | Email aqualife@aqualife.com.au
Unit 4, 8-12 Butler Way, Tullamarine VIC 3043 www.aqualife.com.au**

COMPONENT LIST AQUACLAVE M900 KIT

1 x AquaClave Cartridge Set (4)	1 x 60 Litre Tank SS316 Connection
1 x M900 Membrane	1 x TDS Meter
3 x USA Cartridge Housings	1 x Low Pressure Switch
1 x 50HF Pump	1 x High Pressure Switch
1 x Auto Control Box	1 x Leak Detector Kit
1 x 24V Power Supply	2 x Date Change Due Labels
1 x Secure Power Supply Lead	1 x Owner Installation & Operation Manual
1 x 60 Litre Tank	1 x Warranty Registration Card
1 x RMC 500 kPA Pressure Limiting Valve	

Please note that Conductivity reading as per AS4187 Table 7.2 Final rinse water for Washer disinfectors complies on the basis of incoming feed water quality. In case the conductivity reading is above the specified limits, additional water treatment has to be provided. Please contact Aqualife for more information.

Aqualife Products (Australia Pty Ltd) Freecall 1800 623 887

Phone 03 9310 5259 | Facsimile 03 9310 5387 | Email aqualife@aqualife.com.au

Unit 4, 8-12 Butler Way, Tullamarine VIC 3043

www.aqualife.com.au