

PENCIL LOAD KIT FOR ULTRASONIC CLEANERS



# STSPRO HEX TEST

TO AS2773.2

**INCLUDES**  
Hex Test Disc  
HB Pencil  
Ring Holder  
Recording Table  
Zip Lock Bag



*The STS Hex Test kit is manufactured to AS2773.2-1999, Ultrasonic Cleaners for health care facilities, Part 2: Benchtop for the Pencil Load listed under 6.2 (b)(ii) & is designed for up to 60 tests. More details & instructions for use found enclosed.*

**[www.stspro.com.au](http://www.stspro.com.au)**