

Introducing CoFlex® TLC:

Two Layer Compression Bandage system designed to meet all your needs.



COMPRESSION BANDAGES

CoFlex® TLC

Two Layer Compression

Non-Medicated

- New itch control first layer. CoFlex® TLC has added a new formula to Layer 1 that now includes aloe to make it more itch-resistant and comfortable on sensitive skin.
- CoFlex® TLC uses an odour control foam 1st layer that absorbs 85% of the odours of standard foam.
- CoFlex® TLC outperforms other 2 layer compression kits from the superior quality of its patented Latex Free cohesive bandage formula and consistent non-constricting compression.

Medicated


- Layer 1 is a Zinc or Calamine impregnated foam roll that rolls on cleanly and easily as a dry foam roll.
- Layer 2 is the only Latex Free compression bandage to stick directly on a moist dressing allowing the use of a low-profile medicated compression.
- Visual indicators have been incorporated into the cohesive bandage to assist with application in all medicated and non - medicated kits.
- Calamine and Zinc Impregnated Foam rolls are available separately.

• For ABI > 0.5

CoFlex Zinc-TN, Calamine-TN, TLC Lite
25-30mmHg light compression

• For ABI > 0.8

CoFlex TLC TN, TLC XL, Calamine -SC, Zinc-SC
35-40mmHg standard compression

Code	Application	Size	
A8840UBC-TN	CoFlex® TLC Calamine Lite Kit	L1 - 10cm x 5.5m unstretched, L2 - 10cm x 4.76m unstretched	
A8840UBC-SC	CoFlex® TLC Calamine STD Kit	L1 - 10cm x 5.5m unstretched, L2 - 10cm x 4.76m unstretched	
AUB840C-01	CoFlex® TLC Calamine Medicated Roll	10cm x 5.5m unstretched	
A8840UBZ-TN	CoFlex® TLC Zinc Lite Kit	L1 - 10cm x 5.5m unstretched, L2 - 10cm x 4.76m unstretched	
A8840UBZ-SC	CoFlex® TLC Zinc STD Kit	L1 - 10cm x 5.5m unstretched, L2 - 10cm x 4.76m unstretched	
AUB840Z-01	CoFlex® TLC Zinc Medicated Roll	10cm x 5.5m unstretched	
A7802TLC-TN	CoFlex® TLC Lite Kit	L1 - 10cm x 3.1m unstretched, L2 - 10cm x 3.5m unstretched	
A7800TLC-TN	CoFlex® TLC STD Kit	L1 - 10cm x 3.1m unstretched, L2 - 10cm x 3.5m unstretched	
A7800TLC-XL	CoFlex® TLC X/Length STD Kit	L1 - 10cm x 4.9m unstretched, L2 - 10cm x 4.8m unstretched	

Sutherland Medical Pty Limited

Unit 2, 195 Chesterville Road Moorabbin Vic 3189
Ph: 1300 664 027
Email: orders@sutherlandmedical.com.au
Web Site: www.sutherlandmedical.com.au

ARTG 305810 (rolls only)
ARTG 218947 (Kits only)

Version 5 / Created 7/10/2020