

Spinal Needles


- Ergonomic design
- Clear polycarbonate hub
- Rapid and easy recognition of cerebrospinal fluid (CSF) flashback

Quincke (K-3 Lancet)

The K-3 Lancet point (modified Quincke), with secondary side bevels, lowers puncture force¹ and minimises dural trauma.



Uniever Quincke Spinal Needles

Product Code	Colour	Gauge	Needle Length	Carton Qty
			mm	Each
1110866		18G	90	25
1120866		19G	90	25
1131166		20G	90	25
1140476		21G	70	25
1140766		21G	90	25
1150376		22G	38	25
1150976		22G	70	25
1151266		22G	90	25
1160676		23G	70	25
1160966		23G	90	25
1182966		25G	90	25

Note: Uniever Spinal Needles are not available in New Zealand