



LATEX FREE



SAFEED[®] Extension Tube Set

The highly flexible and kink resistant tubing minimises disruption to the flow of fluids. The rotating hub prevents the tube from coiling up when connecting between a device and the IV access point. With a luer lock on one end and a slip tip at the other, the extension set can easily be connected to an infusion set, syringe or three-way stopcock.

Product Code	Capacity	Tubing Length	Inner Diameter	Outer Diameter	Carton
	mL	cm	mm	mm	
SF*ET0222L22	0.25	25	1.1	2.6	50
SF*ET0727L22	0.75	75	1.1	2.6	50
SF*ET152EL22	1.5	150	1.1	2.6	50

Terumo Australia Pty Ltd

Macquarie Park NSW 2113 Australia
Mount Albert Auckland 1025 New Zealand
E: cs_australia@terumo.co.jp
Orders: au_orders@terumo.co.jp

T: 1800 837 866
T: 0800 66 77 57

SAFEED[®] is a registered trademark of Terumo Corporation.
© 2017 Terumo Australia Pty Ltd. GHP170529D728V2 | PBD00226

