



# NEW PRODUCT

## PML1 LED EXAMINATION LIGHT



The PML1 LED examination light is a great all-round diagnostic light. It provides a soft white light for better colour rendition and easier diagnosis and treatment.

With a small light head and convenient handle with integrated on/off switch it is easy to move without blocking the vision of the operator. The flexible goose neck allows positioning in a diverse range of angles for optimal illumination of the work area.

A stepless dimmer control allows the operator to set the precise illumination required for each use without being constrained to factory preset levels and the heat free LED provides bright white light for up to 50,000 hours of operation.

It is designed for use in Consulting Rooms, Recovery Areas, Emergency Rooms, Gynaecology Examinations and Minor Surgical Procedures.

### Features

- Powerful long-life LED with up to 50,000 hours of operation
- Low energy consumption
- IR and heat free cold white light
- Comfortable guiding handle with integrated On/Off switch
- Cooling through a heatsink so no fans or moving parts
- Superb colour rendering
- Mobile base and wall mount provided with every light

### Specifications:

Light Source	3W LED
Light Intensity	89,000 Lux @ 30cm 40,000 Lux @ 50cm
Light Field Diameter	80-150mm
Colour Temperature	4,500-5,000K
Lamp Head Diameter	75mm
Lamp Head Length	110mm
Base Diameter	530mm
Base Height	790mm
Gooseneck Length	760mm
Total Height	1,830mm

### Order Codes:

- PLM1WH WHITE
- PLM1BL BLACK



### Pacific Medical

Phone: 1300 076 332  
Email: [info@pacificmedical.com.au](mailto:info@pacificmedical.com.au)  
[www.pacificmedical.com.au](http://www.pacificmedical.com.au)



© Pacific Medical June 2021 - All rights reserved.



Team Medical [www.teammed.com.au](http://www.teammed.com.au) | [info@teammed.com.au](mailto:info@teammed.com.au) | phone 1300 22 44 50  
BP Medical (WA) | [www.bpmedical.com.au](http://www.bpmedical.com.au) | [info@bpmedical.com.au](mailto:info@bpmedical.com.au) | phone 1300 33 55 60