



# NEW PRODUCT

## PML2 LED EXAMINATION LIGHT



The PML2 LED examination light has 6 powerful LED's offering powerful illumination and great value. It provides white light for better colour rendition and easier diagnosis and treatment.

With a small light head with integrated handle it is easy to move without blocking the vision of the operator. The flexible goose neck allows positioning in a diverse range of angles for optimal illumination of the work area.

A stepless dimmer control allows the operator to set the precise illumination required for each use without being constrained to factory preset levels and the heat free LED provides bright white light for up to 50,000 hours of operation.

It is designed for use in Consulting Rooms, Recovery Areas, Emergency Rooms, Gynaecology Examinations & Minor Surgical Procedures.

### Features

- Powerful long-life LED with up to 50,000 hours of operation
- Low energy consumption
- IR and heat free cold white light
- Small light head with integrated handle
- Cooling through a heatsink so no fans or moving parts
- Superb colour rendering
- Mobile base and wall mount provided with every light

### Specifications:

Light Source	6W LED
Light Intensity	78,000 Lux @ 30cm 40,000 Lux @ 50cm 20,000 Lux @ 100cm
Light Field Diameter	70-200mm
Colour Temperature	4,500-5,500K
Lamp Head Diameter	80mm
Base Diameter	530mm
Base Height	790mm
Gooseneck Length	760mm
Total Height	1,815mm

### Order Codes:

- PLM2WH WHITE
- PLM2BL BLACK



### Pacific Medical

Phone: 1300 076 332  
Email: [info@pacificmedical.com.au](mailto:info@pacificmedical.com.au)  
[www.pacificmedical.com.au](http://www.pacificmedical.com.au)



© Pacific Medical June 2021 - All rights reserved.



Team Medical [www.teammed.com.au](http://www.teammed.com.au) | [info@teammed.com.au](mailto:info@teammed.com.au) | phone 1300 22 44 50  
BP Medical (WA) | [www.bpmedical.com.au](http://www.bpmedical.com.au) | [info@bpmedical.com.au](mailto:info@bpmedical.com.au) | phone 1300 33 55 60