

OSCIMED ERGO II

Cast cutter saw



ERGONOMIC, COMPACT, LIGHTWEIGHT AND EASY TO USE

OSCIMED ERGO II is a compact, lightweight plaster saw ergonomically shaped to enable precision work without any effort, even with only one hand. OSCIMED ERGO II is robustly built and equipped with a cooling fan for intensive use without overheating. The device is powered by a safe low voltage of 24 volts and complies with the latest medical standards.

Rated power	180 Watts
Length of the saw body	230 mm
Effective diameter	58 mm
Weight without power supply	895 gr.
Speed of motor rotation	14'500 rev/min
Noise level	67 dBA
Length of power supply cable	1,10 m
Length of the cable 24 VDC	3,00 m

OSCIMED ERGO Saw:

- with European plug: ref. OSC 340
- with British plug: ref. OSC 340-UK
- 120V/60 Hz Saw with U.S. plug: ref. OSC 340-US



Cooling fan for use without overheating.



A high-performance, easy-to-adapt dust vacuuming system.



Usable with one hand thanks to its very short (230 mm) and lightweight (860 g) shape.



Blade change without tools.