



LHH G2-TW

LHH G2-TW



FINISHES

- /// Tuneable white next-generation LED examination light
- /// SteriTouch® antimicrobial finish
- /// Wall, ceiling, table or trolley mounting

OVERVIEW

Tuneable white illumination has obvious benefits for examination lighting. The new LHH G2-TW heralds a new era for Luxo, offering a high output LED module that can be adjusted to one of three colour temperatures. Its highly flexible self-balancing arm in conjunction with excellent light output gives the LHH G2 TW diagnostic luminaires vast popularity among medical professionals. The frictionless 3d movement and ease of use make the LHH G2 TW the most flexible examination luminaire on the market? and the best in its class.

DESIGN

SteriTouch® antimicrobial white powder-coated aluminium and steel construction. 3D positioning with 3-axis movement and perfect balance. Protective covers on arms to prevent crush-points. Smooth membrane switch for ease of use and cleaning. Fully dimmable.

OPTIC

High output LEDs with refractor lenses for increased light output. LEDs are protected by tempered glass.

ILLUMINATION

25,000lx @ 1000mm.

65,000lx @ 500mm.

Light Field Ø160mm (d50) / Ø350mm (d10)

ELECTRICAL

Integral electronic gear.

LIGHT SOURCE

3500K / 4000K / 4500K CCT LED, Ra95 (R9:98 - 4000K).

VERSION	MOUNT ₁	WATTAGE (MAX.)	ILLUMINANCE @ 1M	D50 LIGHT FIELD DIAM.	D10 LIGHT FIELD DIAM.	COLOUR	CODE	DIMENSIONS
Image 1	Ceiling / Wall	20W	25,000lx	Ø160mm	Ø350mm	3500K/4000K /4500K CRI95 (R9:98 4000K)	LHG026542	Diagram 1

₁ Mounting accessories (e.g. Wall brackets, ceiling brackets, table mounts and trolley mounts) are not included - please specify required accessory to suit application.

VERSION



IMAGE 1

DIAGRAM




DIAGRAM 1

VERSION	DESCRIPTION	CODE	DIMENSIONS
Image 2	60cm Ceiling mounting bracket, suitable for a ceiling fixing point of 2.7m above floor	BRK002529LA	No Diagram
No Image	Extension rod accessory (must be added to a standard ceiling mount for ceilings between 2.7 and 3.1m affl)	CMBR1000	No Diagram
Image 3	Table clamp	BRK014904	No Diagram
Image 4	Wall bracket with locking device	BRK0017516	No Diagram
Image 5	Trolley mount	BRK025205	No Diagram

VERSION



IMAGE 2



IMAGE 3



IMAGE 4



IMAGE 5