

Standardise across facilities for greater efficiency and consistent data capture or customise for your facility's needs.

	WALL BOARD	GS 777 WALL TRANSFORMER	COAXIAL OPHTHALMOSCOPES	PANOPTIC™ OPHTHALMOSCOPES	DIAGNOSTIC OTOSCOPES	DIAGNOSTIC MACROVIEW™ OTOSCOPES	767 ANEROID	KLEENSPEC® DISPENSERS	LOCKING COLLARS	SURETEMP® PLUS THERMOMETER
										
Part #		Country	11720	11820	25020	23810	7670-01	52400-PF	76700	01690-300
Integrated Wall Systems										
77796-M6	•	77716	•			•	•	•		
77796-MX	•	77716	•			•	•	•		•
77796-2MP	•	77716		•		•	•	•		
77796-2MPX	•	77716		•		•	•	•		•
77796-MXNOBP	•	77716	•			•	•	•		•
77797-2X	•	77717	•		•		•	•		•
77797-M2	•	77717	•			•	•	•		
77797-M2X	•	77717	•			•	•	•		•
77797-2MP2	•	77717		•		•	•	•		
77797-2MP2X	•	77717		•		•	•	•		•
Wall Transformer Sets										
77516		77716						•		
77716-71M		77716	•			•			•	
77716-82M		77716		•		•			•	
77716 COMP		77716	•		•				•	
77717-72M2		77717	•			•				
77717-82M2		77717		•		•				
77735-2x		77718	•		•		• (KPA)	•		•
77747		77717	•		•				•	

\*Not all the configurations are available in every country. Please check with your local Welch Allyn representative.

The United States Green Building Council's Leadership in Energy and Environmental Design (LEED) is an internationally recognised green building certification system providing benchmarks for the design, construction and operation of high-performance green buildings.

Welch Allyn is proud to have earned a Gold rating in the LEED certification program.

To place an order, call our Customer Service Department at:

Australia: +61 29 638 3000

China: +86 512 68086078

Japan: +813 3219 0071

Malaysia: +603 788 43341

Singapore: +65 6419 8100

South Africa: +27 11 777 7555 or contact your local Welch Allyn representative.

© 2012 Welch Allyn SM4085AP Rev A www.welchallyn.com



WelchAllyn®  
Advancing Frontline Care™

Welch Allyn  
Green Series™ 777  
Integrated Wall System



ENERGY STAR  
PARTNER

WelchAllyn®

Advancing Frontline Care™

Welch Allyn Green Series™ 777 Integrated Wall System




## Welch Allyn Green Series™ 777 Integrated Wall System

Allows for reliable diagnostic patient data in a standardised and energy-efficient solution to improve your bottom line.

### Welch Allyn Green Series™ 777 Wall Transformer – configured to each facility's needs.

- Help reduce facility cost with a more energy-efficient product than the 767 Wall Transformer
- Help improve workflow and save energy with handles that automatically turn on/off
- Quickly and easily upgrade from the Welch Allyn 767 Wall Transformer to the Green Series 777 Wall Transformer
- No on/off switch reduces contact, to help reduce the risk of cross-contamination



### Welch Allyn Otoscopies & Ophthalmoscopes

- Standardise or choose the technology you need for each facility
- Upgrade to the Welch Allyn LED lamp for greater energy efficiency and lower cost of ownership



#### Welch Allyn MacroView™ Otoscope

- Almost twice the field of view and 30% greater magnification than traditional otoscopes
- Allows for simplified diagnosis of Otitis Media through superior visualisation of the tympanic membrane



#### Welch Allyn Pneumatic Otoscope

- Rotating lens with ultraseal enables easy pneumatic otoscopy
- Easier instrumentation to allow a faster and more comfortable exam for the patient



#### Welch Allyn PanOptic™ Ophthalmoscope

- Fast, easy entry into undilated pupils with 5X larger view of the fundus for a streamlined exam
- 25° field of view vs. the standard 5° field of view provides optimal visualisation



#### Welch Allyn Coaxial Ophthalmoscope

- 18 unique aperture/filter combinations for greater versatility
- 28 focusing lenses with a range of -25 to +40 diopters



Welch Allyn industry-leading blood pressure and thermometry devices provide you with fast, accurate and reliable patient data you can trust – to help improve patient care and clinical decisions.



#### Welch Allyn Wall Aneroid

- Recessed, laser-engraved dial face with high-contrast pointer provides certified accuracy to reduce errors and improve patient care
- Lifetime calibration warranty reduces cost of ownership



#### Welch Allyn Connex® ProBP™ Digital Blood Pressure Device

- Connectivity to most leading EHR systems, to improve patient documentation accuracy and efficiency
- Improve patient comfort and increase throughput with fast, 15-second NIBP readings
- Portable, so you can take BP readings just about anywhere



#### Welch Allyn Spot Vital Signs® Device

- Improve workflow by using a device that includes NIBP, pulse rate and MAP
- Optional SureTemp® thermometry and Masimo® or Nellcor® pulse oximetry



#### Welch Allyn Spot Vital Signs® LXi Device

- Improve patient comfort and increase throughput with 15-second NIBP readings taken as the cuff inflates
- Connectivity to most leading EHR systems, to improve patient vitals documentation accuracy and efficiency
- Improve workflow by using a device that includes NIBP, pulse rate, MAP, thermometry and SpO<sub>2</sub>



#### Welch Allyn SureTemp® Plus Electronic Thermometer

- Provides accurate oral temperatures in 4 to 6 seconds, axillary and rectal temperatures in 10 to 13 seconds and adult axillary temperatures in 12 to 15 seconds
- Reduce risk of cross-contamination with single-use probe covers
- Improve efficiency with convenient storage for 75 probe covers
- Takes approximately 6,000 readings on 3 AA batteries and provides last temperature recall