

**Baxter**

## **CP 150** ELECTROCARDIOGRAPH WITH OPTIONAL SPIROMETRY

---

Simple, Fast, Reliable





## HELP YOUR PRACTICE RUN MORE EFFICIENTLY WITH THE **CP 150** RESTING ECG




### SIMPLE

- Intuitive touchscreen display for more accuracy.
- Lead quality display identifies connection problems before completing a test.



### FAST

- Instant-On Feature helps eliminate typical start-up delays.
- Large display provides preview options for three, six or 12-lead ECGs.
- ECG Print Preview verifies test quality prior to printout.



of **CP 150** ECG users surveyed say the touchscreen display and simple user interface have helped them improve their practice efficiency.<sup>1</sup>



“

I find the **CP 150** [ECG] to be easy to use, reliable and safe for adult and pediatric patients.

”

- Brad Ronevich,  
Richmond Heights Medical Center RN

### SPIROMETRY OPTION

The **CP 150** ECG lets you perform two reimbursable diagnostic exams with one device. Help manage your COPD patients and maximize HEDIS scores with a spirometry solution that includes automatic interpretation and reversibility analysis. Plus, you can customize report formats and screen displays to meet your workflow needs.

- Perform FVC and pre and post BD testing with automatic interpretation and reversibility analysis.
- Customise report formats and screen displays to fit your needs and workflow.
- Ambient sensors can eliminate the need for the user to find and enter values for temperature, pressure and humidity.
- Manage patients with COPD to help you maximise HEDIS scores.



## ORDERING INFORMATION [SELECTED DEVICE CONFIGURATIONS\*]

PART NUMBER	DESCRIPTION
CP150A-1EN66	CP 150 Interpretive Electrocardiograph
CP150AS-1EN66	CP 150 Interpretive Electrocardiograph with Spirometry

## ACCESSORIES

PART NUMBER	DESCRIPTION
715006 –AUS	CP 150 Multifunction Electrode Adapter
108071	Tab Electrode Case (5000)
108770	Tab Electrode Box (500)
105353	CP 150 Recording Paper (1000 sheets)
719653	10-lead Patient Cable, AHA Banana, CP 150
105342	CP 150 Hospital Cart
105343	CP 150 Cable Arm & Shelf

## SPIROMETRY ACCESSORIES

PART NUMBER	DESCRIPTION
105660	Spirometry Upgrade for CP 150 (works with any standard CP150A-1ENB or CP 150 - 1ENB configuration)
720705	Disposable Flow Transducers, 25 pk, CP 150
720706	Disposable Flow Transducers, 100 pk, CP 150
410370	Spirometry Pressure Tube and Handle

## WARRANTY INFORMATION

PART NUMBER	DESCRIPTION
CP 150 Electrocardiograph	3 years
Patient Cable/Battery	90 days
Spirometer	1 year



Experience simple, streamlined ECG testing.  
Scan the QR code to learn more today.

For more information, please contact your Baxter representative or call Baxter Customer Care on 1800 445 576.

**ALWAYS READ THE LABEL AND FOLLOW THE DIRECTIONS FOR USE** CONTACT YOUR BAXTER REPRESENTATIVE FOR THE LATEST VERSIONS

\*For more device configurations, see CP 150 Part Number Configuration document.

### References

<sup>1</sup> Survey by Walker Information, Inc. of clinicians who have purchased the CP 150.

Hill-Rom reserves the right to make changes without notice in design, specifications and models.

The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

Baxter, CP 150 and Hillrom are trademarks of Baxter International Inc.

AUTHORISED AUSTRALIAN SPONSOR Welch Allyn Australia Pty. Limited. 1 Baxter Drive, Old Toongabbie NSW 2146.

[hillrom.com.au](http://hillrom.com.au)

**Baxter**

ANZ-FLC48-240001\_r1\_05/2024