



SAFETY DATA SHEET

Section 1: IDENTIFICATION

V-WIPES

Recommended Use: Cleaning/disinfection of surfaces and non-critical medical devices.

Product Code: 210561 (4 x 100 Canisters), 210579 (12 x 50 Flat Packs), 210634 (12 x 80 Flat Packs),

Whiteley Corporation Pty Ltd

(A.C.N. 000 906 678)

Postal Address: P. O. Box 1076 North Sydney NSW 2059

Telephone Number: (02) 9929 9155 Facsimile: (02) 9929 9077

Web: www.whiteley.com.au

Emergency Telephone Number: Poisons Information Centre (National) 131126

Distributed in New Zealand by:

Whiteley Corporation NZ Limited

Address: PO Box 22-519 Otahuhu, Auckland, New Zealand 1640

Telephone Number: 092737313, 0800 257 352

Emergency Telephone Number: New Zealand Poisons Centre 0800 POISON or 0800 764 766

Section 2: HAZARDS

GHS Classification

Not a dangerous substance or mixture according to GHS or HSNO criteria.

Signal Word None

Hazard Statements

None

Precautionary Statements

None

Response Statements

None

Storage Statements

None

Disposal Statements

None

Section 3: COMPOSITION INFORMATION

Ingredient	CAS No	Proportion
Ethanol	64-17-5	<10%
Ingredients deemed not to be hazardous	Not applicable	To 100%

Section 4: FIRST AID	
-----------------------------	--

Eye (Contact)	Hold eyelids apart and flush the eye continuously with running water.
Skin (Contact)	Remove contaminated clothing and flush skin and hair with running water.
Inhalation (Breathing)	Remove to fresh air.
Ingestion (Swallowing)	DO NOT induce vomiting. Give water to drink if conscious. For advice, contact a Poisons Information Centre (Phone 131126) or a doctor.
Advice to Doctor	Treat symptomatically for neutral detergent.
First Aid Facilities	Ensure an eye wash is available and ready for use.
Additional Information	No aggravated medical conditions are known to be caused by exposure to this product.

Section 5: FIREFIGHTING MEASURE	
--	--

Suitable Extinguishing Media	Use water fog (or if unavailable fine water spray), dry chemical, carbon dioxide or alcohol stable foam.
Hazards From Combustion Products	Carbon dioxide and carbon monoxide may be produced in the case of fire or during thermal decomposition.
Precautions For Fire Fighters and Special Protective Equipment	Keep containers cool by spraying with water.
Additional Information	Hazchem Code – Not applicable.

Section 6: ACCIDENTAL RELEASE MEASURES	
---	--

Emergency Procedure	SAA/SNZ HB76: Dangerous Goods – Initial Emergency Response Guide – Not applicable.
Spills / Clean up	Disposal to infectious waste if contaminated recommended. If not contaminated then disposal to normal waste acceptable. Refer to State/Territory Land Waste Management Authorities if applicable. Containers are recyclable and can be disposed of by a licensed waste contractor. Containers can be disposed of to general waste or rinsed thoroughly and recycled.

Section 7: HANDLING AND STORAGE	
--	--

Precautions for Safe Handling	Contact Whiteley Corporation sales representative for advice when using this product for any application other than that outlined on the label or technical bulletin. Any non-intended or non-authorized use of this product may result in personal injury or damage to equipment. Store product in original container. Wash thoroughly after handling product.
Conditions for Safe Storage	Store in a cool, dry, well ventilated area. Keep container tightly sealed. Store below 30°C.

Section 8: EXPOSURE CONTROL/PERSONAL PROTECTION	
--	--

National Exposure Standards – Source: National Exposure Standards for Atmospheric Contaminants in the Occupational Environment [NOHSC:1003] & New Zealand Workplace Exposure Standards 2002.

<u>Ingredient</u>	<u>CAS No</u>	<u>ES-TWA</u>	<u>ES-STEL</u>
Ethanol	64-17-5	1000ppm	-

Biological Limit Values No data available.

Engineering Controls Not required.

Personal Protective Equipment



Eye/face protection – Eye protection is not usually necessary. If risk of eye contact exists then safety glasses or chemical-resistant goggles are recommended.

Skin protection – Use of gloves is recommended to prevent skin contact.

Respiratory protection – Respirator is not usually necessary but if required use a half face filter respirator suitable for organic vapours.

Section 9: PHYSICAL AND CHEMICAL PROPERTIES

Appearance	White and pink apertured wipe	Odour	Citrus
Odour Threshold	No data available	pH (lotion)	6 – 8
Freezing Point	No data available	Initial Boiling Point	No data available
Boiling Range	No data available	Flash point	No data available
Evaporation Rate	No data available	Flammability	No data available
Upper flammability limit	No data available	Lower flammability limit	No data available
Vapour Pressure	No data available	Vapour Density	No data available
Relative Density	No data available	Solubility	No data available
Partition Coefficient: n-octanol/water	No data available	Autoignition temperature	No data available
Decomposition temperature	No data available	Viscosity	No data available

Section 10: STABILITY AND REACTIVITY

Chemical Stability	Product is stable and will not undergo any hazardous reactions under normal conditions of use and storage.
Conditions to avoid	Avoid high temperatures (store below 30°C).
Incompatible materials	None known. Do not mix with other chemicals.
Hazardous decomposition products	None known.
Hazardous reactions	None known.

Section 11: TOXICOLOGICAL INFORMATION

Acute toxicity	No data available.
Skin corrosion/ irritation	May cause irritation.
Serious eye damage/ irritation	May cause irritation and reddening.

Respiratory or skin sensitisation	No data available.
Germ Cell Mutagenicity	No data available.
Carcinogenicity	No components of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.
Reproductive toxicity	No data available.
Specific Target Organ Toxicity – Single Exposure	No data available.
Specific Target Organ Toxicity – Repeated Exposure	No data available.
Aspiration Hazard	No data available.

Section 12: ECOLOGICAL INFORMATION

Ecotoxicity	No data available.
Persistence and degradability	No data available.
Bioaccumulative potential	No data available.
Mobility in soil	No data available.
Other adverse effects	No data available.

Section 13: DISPOSAL CONSIDERATIONS
--

Disposal method	Disposal to infectious waste if contaminated recommended. If not contaminated then disposal to normal waste acceptable. Refer to State/Territory Land Waste Management Authorities if applicable. Containers are recyclable and can be disposed of by a licensed waste contractor. Containers can be disposed of to general waste or rinsed thoroughly and recycled.
Special precautions	None.

Section 14: TRANSPORT INFORMATION
--

Not classified as Dangerous Goods by the criteria of the Australian Dangerous Goods Code (ADG Code).

UN Number	Not applicable.
UN Proper Shipping Name	Not applicable.
Class and subsidiary risk	Not applicable.
Packing Group	Not applicable.
Special precautions for	Not applicable.
Hazchem Code	Not applicable.

Section 15: REGULATORY INFORMATION

Poisons Schedule (SUSMP): Not applicable.

HSNO Approval Code: Not applicable.

All ingredients are listed in the Australia Inventory of Chemical Substances (AICS).

This document has been produced in accordance with the requirements of the Globally Harmonised System of Classification and Labelling.

Section 16: OTHER INFORMATION

Date of preparation: 17th September 2021

© This document remains property of Whiteley Corporation Pty Ltd.
Alterations are not permitted without prior written authorisation from Whiteley Corporation Pty Ltd.