

HEALTH IS THE HEART OF OUR BUSINESS.

HOLTER-RR 24H ABPM

GENTLE AMBULATORY BLOOD PRESSURE MONITORING DIAGNOSTICS

Meaningful data: The Holter-RR guarantees reliable ambulatory blood pressure monitoring.

The Auto Feedback Logic (AFL) allows quiet and fast readings with low cuff pressure.

For precise interpretation threshold tables for children and adolescents are integrated.



MANAGE COMPLEX DATA EASILY

AMEDTEC ECGpro®



TECHNICAL DATA

Pressure measurement range	systolic 60 to 290 mmHg diastolic 30 to 195 mmHg
Accuracy	+/-3 mmHg in the range indicated
Static pressure range	0 to 300 mmHg
Pulse range	30 to 240 beats per minute
Method	oscillometric
Measurement intervals	0, 1, 2, 3, 4, 5, 6, 10, 12, 15, 20 or 30 measurements per hour
Monitoring protocols	2 programmable and 7 fixed protocols
Storage capacity	300 measurements
Operating temperatures	+10 °C bis +40 °C
Operating humidity	15 % to 90 %
Dimensions	128 mm x 75 mm x 30 mm
Weight	approx. 240 g incl. rechargeable batteries
Power supply	2 Mignon AA alkaline or rechargeable batteries
Interfaces	Cable: serial interface RS232, converter for USB

ACCESSORIES



pouch for Holter-RR



cuffs in size XS, S, M, L, XL



Team Medical Supplies

Phone 1300 22 44 50
Fax 1300 22 44 60
info@teammed.com.au
www.teammed.com.au



BP Medical (WA)

Phone 1300 33 55 60
Fax 1300 33 55 70
info@bpmedical.com.au
www.bpmedical.com.au

Manufacturer:

IEM GmbH
Gewerbepark Brand 42
52078 Aachen, Germany



Certificate:
0044 DIN EN ISO 13485

Errors and omissions excepted
M0731 018-0402-2_04

AMEDTEC Medizintechnik Aue GmbH
Schneeberger Straße 5 · 08280 Aue · Germany
Phone: +49 (0)3771 59 82 70
Fax: +49 (0)3771 59 82 790
info@amedtec.de · www.amedtec.de

