



Patient name: **Smith, John**
Patient Number: **00000999999999**
Date of Birth: **8/2/1945 (56 Y)**
Recording Date: **5/29/2002 12:18:00 PM**

Sex: **Male**
Weight: **170lbs**
Height: **70inch**
Address: **Democity, 12345678, 123 South Street**

Dr. D. Doe
21 North Street, Sample City
Phone: **555-12345** Fax: **555-123456**



Day interval

Valid measurements	34 of 34 = 100%	
Average	147/124/104mmHg, 76bpm	< 135/85mmHg
Pulse pressure	43.0mmHg	
≥ 140mmHg syst.	82.4%	< 25%
≥ 90mmHg diast.	100.0%	< 25%

Night interval

Valid measurements	13 of 13 = 100%	
Average	125/106/84mmHg, 75bpm	< 120/75mmHg
Pulse pressure	41.0mmHg	
≥ 125mmHg syst.	61.5%	< 25%
≥ 80mmHg diast.	84.6%	< 25%

Day/Night-decrease	15.0% / 19.2%	10% - 20%
	Normal Dipper	

Summary

Case history, diagnosis or comments go here.

Dr. D. Doe

Number	Time	Day / Night	Systolic (mmHg)	MAP (mmHg)	Diastolic (mmHg)	HR (bpm)
1	5/29/2002 12:18 PM	☀	144	117	108	72
2	5/29/2002 12:30 PM	☀	128	107	100	75
3	5/29/2002 1:00 PM	☀	134	115	102	74
4	5/29/2002 1:30 PM	☀	155	122	112	70
5	5/29/2002 2:00 PM	☀	142	118	103	72
6	5/29/2002 2:30 PM	☀	142	113	98	74
7	5/29/2002 3:00 PM	☀	155	133	112	69
8	5/29/2002 3:30 PM	☀	145	118	105	78
9	5/29/2002 4:00 PM	☀	142	123	99	72
10	5/29/2002 4:30 PM	☀	145	133	106	68
11	5/29/2002 5:00 PM	☀	138	111	94	72
12	5/29/2002 5:30 PM	☀	155	129	108	78
13	5/29/2002 6:00 PM	☀	172	143	131	75
14	5/29/2002 6:30 PM	☀	166	148	124	78
15	5/29/2002 7:00 PM	☀	156	135	114	78
16	5/29/2002 7:30 PM	☀	148	113	102	76
17	5/29/2002 8:00 PM	☀	152	122	108	77
18	5/29/2002 8:30 PM	☀	142	125	94	82
19	5/29/2002 9:00 PM	☀	144	124	98	74
20	5/29/2002 9:30 PM	☀	143	126	99	79
21	5/29/2002 10:00 PM	☀	138	122	92	73
22	5/29/2002 10:30 PM	☀	142	128	93	71
23	5/29/2002 11:00 PM	☀	133	107	94	69
24	5/29/2002 11:30 PM	☀	135	111	91	72
(25)	5/29/2002 11:52 PM	☾				
26	5/30/2002 12:00 AM	☾	128	115	84	74
27	5/30/2002 12:30 AM	☾	126	111	82	73
28	5/30/2002 1:00 AM	☾	127	114	83	80
29	5/30/2002 1:30 AM	☾	129	114	88	79
30	5/30/2002 2:00 AM	☾	125	103	84	72
31	5/30/2002 2:30 AM	☾	120	105	88	79
32	5/30/2002 3:00 AM	☾	128	111	90	74
33	5/30/2002 3:30 AM	☾	122	102	78	73
34	5/30/2002 4:00 AM	☾	121	101	82	75
35	5/30/2002 4:30 AM	☾	127	109	84	78
36	5/30/2002 5:00 AM	☾	123	92	80	68
37	5/30/2002 5:30 AM	☾	123	98	79	70
38	5/30/2002 6:00 AM	☾	131	100	88	76
(39)	5/30/2002 6:18 AM	☀				
40	5/30/2002 6:30 AM	☀	144	122	94	82
41	5/30/2002 7:00 AM	☀	153	126	103	82
42	5/30/2002 7:30 AM	☀	158	136	117	88
43	5/30/2002 8:00 AM	☀	163	145	128	83
44	5/30/2002 8:30 AM	☀	149	122	104	86
45	5/30/2002 9:00 AM	☀	144	130	102	87
46	5/30/2002 9:30 AM	☀	148	120	105	77
47	5/30/2002 10:00 AM	☀	152	123	110	79
48	5/30/2002 10:30 AM	☀	148	122	102	76
49	5/30/2002 11:00 AM	☀	144	124	98	81

Patient name:	Smith, John	Sex:	Male	Dr. D. Doe
Patient Number:	000009999999999	Weight:	170lbs	21 North Street, Sample City
Date of Birth:	8/2/1945 (56 Y)	Height:	70inch	Phone:555-12345
Recording Date:	5/29/2002 12:18:00 PM	Address:	Democity, 12345678, 123 South Street	Fax:555-123456

Obligatory

Valid measurements	47 of 47 = 100%	Norm values
Total average	141/119/99mmHg, 76bpm	> 90%
Total pulse pressure	42.0mmHg	< 130/80mmHg

Day interval

Valid measurements	34 of 34 = 100%	
Average	147/124/104mmHg, 76bpm	< 135/85mmHg
Pulse pressure	43.0mmHg	
≥ 140mmHg syst.	82.4%	< 25%
≥ 90mmHg diast.	100.0%	< 25%
Minimum syst.	128mmHg (5/29/2002 12:30 PM)	
Maximum syst.	172mmHg (5/29/2002 6:00 PM)	
Minimum diast.	91mmHg (5/29/2002 11:30 PM)	
Maximum diast.	131mmHg (5/29/2002 6:00 PM)	
Minimum HR	68bpm (5/29/2002 4:30 PM)	
Maximum HR	88bpm (5/30/2002 7:30 AM)	

Night interval

Valid measurements	13 of 13 = 100%	
Average	125/106/84mmHg, 75bpm	< 120/75mmHg
Pulse pressure	41.0mmHg	
≥ 125mmHg syst.	61.5%	< 25%
≥ 80mmHg diast.	84.6%	< 25%
Minimum syst.	120mmHg (5/30/2002 2:30 AM)	
Maximum syst.	131mmHg (5/30/2002 6:00 AM)	
Minimum diast.	78mmHg (5/30/2002 3:30 AM)	
Maximum diast.	90mmHg (5/30/2002 3:00 AM)	
Minimum HR	68bpm (5/30/2002 5:00 AM)	
Maximum HR	80bpm (5/30/2002 1:00 AM)	

Day/Night-decrease	15.0% / 19.2%	10% - 20%
	Normal Dipper	


Optional

Early morning average	149/124/99mmHg, 82bpm	
Standard deviation (Day)	9.5/10.0/9.9mmHg, 5.2bpm	< 17,0/13,0mmHg
Standard deviation (Night)	3.4/7.2/3.8mmHg, 3.6bpm	< 13,0/10,0mmHg

Patient name: **Smith, John**
Patient Number: **000009999999999**
Date of Birth: **8/2/1945 (56 Y)**
Recording Date: **5/29/2002 12:18:00 PM**

Sex: **Male**
Weight: **170lbs**
Height: **70inch**
Address: **Democity, 12345678, 123 South Street**


Dr. D. Doe
21 North Street, Sample City
Phone: **555-12345** Fax: **555-123456**



Patient name: **Smith, John**
Patient Number: **000009999999999**
Date of Birth: **8/2/1945 (56 Y)**
Recording Date: **5/29/2002 12:18:00 PM**

Sex: **Male**
Weight: **170lbs**
Height: **70inch**
Address: **Democity, 12345678, 123 South Street**


Dr. D. Doe
21 North Street, Sample City
Phone: **555-12345** Fax: **555-123456**



Summary

Case history, diagnosis or comments go here.

Dr. D. Doe



Summary
 Case history, diagnosis or comments go here.

Dr. D. Doe

Patient name: **Smith, John**
 Patient Number: **00000999999999**
 Date of Birth: **8/2/1945 (56 Y)**
 Recording Date: **5/29/2002 12:18:00 PM**

Sex: **Male**
 Weight: **170lbs**
 Height: **70inch**
 Address: **Democity, 12345678, 123 South Street**


Dr. D. Doe
21 North Street, Sample City
 Phone: **555-12345** Fax: **555-123456**

Number	Time	Day / Night	Systolic (mmHg)	MAP (mmHg)	Diastolic (mmHg)	HR (bpm)	Comments
1	5/29/2002 12:18 PM	☀	144	117	108	72	manuelle Messung
2	5/29/2002 12:30 PM	☀	128	107	100	75	
3	5/29/2002 1:00 PM	☀	134	115	102	74	
4	5/29/2002 1:30 PM	☀	155	122	112	70	
5	5/29/2002 2:00 PM	☀	142	118	103	72	
6	5/29/2002 2:30 PM	☀	142	113	98	74	
7	5/29/2002 3:00 PM	☀	155	133	112	69	
8	5/29/2002 3:30 PM	☀	145	118	105	78	
9	5/29/2002 4:00 PM	☀	142	123	99	72	
10	5/29/2002 4:30 PM	☀	145	133	106	68	
11	5/29/2002 5:00 PM	☀	138	111	94	72	
12	5/29/2002 5:30 PM	☀	155	129	108	78	
13	5/29/2002 6:00 PM	☀	172	143	131	75	
14	5/29/2002 6:30 PM	☀	166	148	124	78	
15	5/29/2002 7:00 PM	☀	156	135	114	78	
16	5/29/2002 7:30 PM	☀	148	113	102	76	
17	5/29/2002 8:00 PM	☀	152	122	108	77	
18	5/29/2002 8:30 PM	☀	142	125	94	82	
19	5/29/2002 9:00 PM	☀	144	124	98	74	
20	5/29/2002 9:30 PM	☀	143	126	99	79	
21	5/29/2002 10:00 PM	☀	138	122	92	73	
22	5/29/2002 10:30 PM	☀	142	128	93	71	
23	5/29/2002 11:00 PM	☀	133	107	94	69	
24	5/29/2002 11:30 PM	☀	135	111	91	72	
(25)	5/29/2002 11:52 PM	☾					Tag->Nacht-Taste
26	5/30/2002 12:00 AM	☾	128	115	84	74	
27	5/30/2002 12:30 AM	☾	126	111	82	73	
28	5/30/2002 1:00 AM	☾	127	114	83	80	
29	5/30/2002 1:30 AM	☾	129	114	88	79	
30	5/30/2002 2:00 AM	☾	125	103	84	72	
31	5/30/2002 2:30 AM	☾	120	105	88	79	
32	5/30/2002 3:00 AM	☾	128	111	90	74	
33	5/30/2002 3:30 AM	☾	122	102	78	73	
34	5/30/2002 4:00 AM	☾	121	101	82	75	
35	5/30/2002 4:30 AM	☾	127	109	84	78	
36	5/30/2002 5:00 AM	☾	123	92	80	68	
37	5/30/2002 5:30 AM	☾	123	98	79	70	
38	5/30/2002 6:00 AM	☾	131	100	88	76	
(39)	5/30/2002 6:18 AM	☀					Tag<-Nacht-Taste
40	5/30/2002 6:30 AM	☀	144	122	94	82	
41	5/30/2002 7:00 AM	☀	153	126	103	82	
42	5/30/2002 7:30 AM	☀	158	136	117	88	
43	5/30/2002 8:00 AM	☀	163	145	128	83	
44	5/30/2002 8:30 AM	☀	149	122	104	86	
45	5/30/2002 9:00 AM	☀	144	130	102	87	
46	5/30/2002 9:30 AM	☀	148	120	105	77	
47	5/30/2002 10:00 AM	☀	152	123	110	79	
48	5/30/2002 10:30 AM	☀	148	122	102	76	
49	5/30/2002 11:00 AM	☀	144	124	98	81	

Patient name: **Smith, John**
Patient Number: **000009999999999**
Date of Birth: **8/2/1945 (56 Y)**
Recording Date: **5/29/2002 12:18:00 PM**

Sex: **Male**
Weight: **170lbs**
Height: **70inch**
Address: **Democity, 12345678, 123 South Street**

Dr. D. Doe
21 North Street, Sample City
Phone: **555-12345** Fax: **555-123456**




Patient name: **Smith, John**
Patient Number: **000009999999999**
Date of Birth: **8/2/1945 (56 Y)**
Recording Date: **5/29/2002 12:18:00 PM**


Sex: **Male**
Weight: **170lbs**
Height: **70inch**
Address: **Democity, 12345678, 123 South Street**

Dr. D. Doe
21 North Street, Sample City
Phone: 555-12345 Fax: 555-123456


Overall



Day



Night



Patient name: **Smith, John**
Patient Number: **000009999999999**
Date of Birth: **8/2/1945 (56 Y)**
Recording Date: **5/29/2002 12:18:00 PM**

Sex: **Male**
Weight: **170lbs**
Height: **70inch**
Address: **Democity, 12345678, 123 South Street**

Dr. D. Doe
21 North Street, Sample City
Phone: **555-12345** Fax: **555-123456**

	SBP Mean (mmHg)	SBP SD (mmHg)	MAP Mean (mmHg)	MAP SD (mmHg)	DBP Mean (mmHg)	DBP SD (mmHg)	HR Mean (bpm)	HR SD (bpm)	Readings
12:00 AM - 1:00 AM	127	1.4	113	2.8	83	1.4	74	0.7	2
1:00 AM - 2:00 AM	128	1.4	114	0.0	86	3.5	80	0.7	2
2:00 AM - 3:00 AM	123	3.5	104	1.4	86	2.8	76	4.9	2
3:00 AM - 4:00 AM	125	4.2	107	6.4	84	8.5	74	0.7	2
4:00 AM - 5:00 AM	124	4.2	105	5.7	83	1.4	77	2.1	2
5:00 AM - 6:00 AM	123	0.0	95	4.2	80	0.7	69	1.4	2
6:00 AM - 7:00 AM	138	9.2	111	15.6	91	4.2	79	4.2	2
7:00 AM - 8:00 AM	156	3.5	131	7.1	110	9.9	85	4.2	2
8:00 AM - 9:00 AM	156	9.9	134	16.3	116	17.0	85	2.1	2
9:00 AM - 10:00 AM	146	2.8	125	7.1	104	2.1	82	7.1	2
10:00 AM - 11:00 AM	150	2.8	123	0.7	106	5.7	78	2.1	2
11:00 AM - 12:00 PM	144	0.0	124	0.0	98	0.0	81	0.0	1
12:00 PM - 1:00 PM	136	11.3	112	7.1	104	5.7	74	2.1	2
1:00 PM - 2:00 PM	145	14.8	119	4.9	107	7.1	72	2.8	2
2:00 PM - 3:00 PM	142	0.0	116	3.5	101	3.5	73	1.4	2
3:00 PM - 4:00 PM	150	7.1	126	10.6	109	4.9	74	6.4	2
4:00 PM - 5:00 PM	144	2.1	128	7.1	103	4.9	70	2.8	2
5:00 PM - 6:00 PM	147	12.0	120	12.7	101	9.9	75	4.2	2
6:00 PM - 7:00 PM	169	4.2	146	3.5	128	4.9	77	2.1	2
7:00 PM - 8:00 PM	152	5.7	124	15.6	108	8.5	77	1.4	2
8:00 PM - 9:00 PM	147	7.1	124	2.1	101	9.9	80	3.5	2
9:00 PM - 10:00 PM	144	0.7	125	1.4	99	0.7	77	3.5	2
10:00 PM - 11:00 PM	140	2.8	125	4.2	93	0.7	72	1.4	2
11:00 PM - 12:00 AM	134	1.4	109	2.8	93	2.1	71	2.1	2