

## Key Differences Between 600M and 800M Audiometers:

The 800M has various improvements:

- **Accessories:** Now includes high-phonics isolation headphones instead of standard headphones.
- **Software Enhancements:** New welcome page with patient visualisation, including an added column showing the most recent test.
- **Patient Management:**
  - Customisable patient ID field added
  - Operator-based patient database: each operator only sees patients they created (this function can be turned off)
- **Test History:**
  - Ability to compare up to three different audiometric tests--ideal for tracking hearing changes over time.