

TM-2440 24-hour Ambulatory Blood Pressure Monitor



A&D’s TM-2440 is an entry-level model designed for ABPM measurement and comes with AFib+™ as standard.

- AFib+ Irregular Heartbeat detection (IHB) and AFib screening
- Automatically adjusts: Proper Cuff Pressure, Exhaust, End of Measurement
- Compact - only 66mm x 24.5mm x 95mm (WxHxD)
- Lightweight - only 120g (without batteries)
- Can be integrated into an Electronic Patient Record System (EPR)
- New pulse correction algorithm for faster measurement & noise rejection
- Battery powered (2 alkaline or Ni-Cd “AA”)
- Output of data to analysis software (included with device)
- Self measurement display, OLED – Ambulatory monitoring display
- USB cable data output to computer or printer



TM-2440 ABPM Kit includes:

- A&D TM-2440 ABPM
- Medium Cuff (20-3cm)
- Carry Pouch
- A&D BPM Data Analysis Software

\$1560

ANDTM2440EX

TM-2441 24-hour Ambulatory Blood Pressure Monitor



The TM-2441 can process five types of self-measurement while recording environmental data (temperature, pressure, activity) through a multi-sensor. It is an “All-in-One” model that covers any environment from the doctor’s office to home.

Includes all features of the TM-2440 model plus:

- PatientView 360: The new TM-2441 provides the Clinician with a complete and comprehensive multi-sensor monitoring solution, allowing the tracking and recording of additional patient-centric information whilst they are wearing the device.
- By recording Activity, Temperature and Air Pressure, the new TM-2441 provides a more inciteful view of the patient beyond just blood pressure recording and heart rate information – helping build a better picture for better clinical decisions.
- New pulse correction algorithm for faster measurement & noise rejection.
- Compact - only 66mm x 24.5mm x 95mm (WxHxD)
- Lightweight - only 135g (without batteries)
- Self measurement option and display, OLED – Ambulatory monitoring display
- Bluetooth (BLE) or USB cable data output to computer or printer
- LCD: Self-BP measurement display
- Auto, measurement mode slide switch



TM-2440 ABPM Kit includes:

- A&D TM-2440 ABPM
- Medium Cuff (20-3cm)
- Carry Pouch
- A&D BPM Data Analysis Software

\$1895

ANDTM2441

For more information on A&D ABPMs download the product brochure from our website



Cuffs for TM-2440 and TM-2441 24- Hour ABPM

A&D Cuff Small 15-22cm	ANDTMCF202A	\$172.55
A&D Cuff Medium 20-31cm	ANDTMCF302A	\$172.55
A&D Cuff Large 28-28cm	ANDTMCF402A	\$210.95
A&D Cuff Extra-Large 36-50cm	ANDTMCF502A	\$231.55



Welch Allyn®

ABPM 7100  
Ambulatory Blood Pressure Monitor

Help avoid white coat hypertension, track sleep readings and tailor drug therapy regimens with simple, 24-hour ambulatory blood pressure monitoring.

Powerful analysis software lets you examine readings more effectively so you can diagnose potential patient complications. Enhance blood pressure management for your patients while making out-of-office monitoring efficient for your practice.

Get fast and accurate results on demand, taking up to 300 measurements in 24 hour period whether it be for tracking sleep blood pressure readings or testing new drug therapy regimens.

ABPM 7100 Kit includes:

- Welch Allyn ABPM 7100
- Adult Cuff (24-32cm)
- Large Adult Cuff (32-38cm)
- Carry Pouch
- Software 7100

\$3495

WAABPM7100S

Visit our website to watch the product video and download the fact sheet



1YR warranty



MEET OUR TEAM: ANGUS

SALES - NSW

Angus has hung up his 8-year warehouse, customer service and marketing hat and moved into sales. We're sure customers will value his extensive product knowledge and his passion for helping. It's been a big year for Angus; he also got married in February and recently welcomed Indie the Dachshund to the family. When Angus isn't playing baseball or watching F1, he likes chilling with mates, gaming and going to the gym.

What is your favourite food? Pizza



MEET OUR TEAM: JACK

SALES - NSW

Jack joined the sales team just over a year ago. With plenty of business experience and a strong ambition to get results for his customers, Jack has immersed himself well in the role. In his downtime, he enjoys watching and playing AFL, going to the gym, doing yoga and spending time with his daughter Frankie. He is also big on snow sports hoping for his next holiday to be in Japan or New Zealand.

Do you play any instruments? The guitar  
What is your biggest pet peeve? Not being organised  
What is something we may not know? I have a twin sister

