

Buffalo ViroVac® 220 Smoke Evacuation System



ViroVac surgical smoke plume evacuator provides ideal flow for general practice, skin cancer clinics and day surgeries. It combines customisable (variable) filter life, a compact shelf design, 3-port filter with filter life tracking technology, while exceeding expectations for value and efficacy.

- **Compact Design**
Allows for easy placement on shelf or counter and accommodates stackability
- **Patented 3-Port Filter Design**
Versatile 3-port filter design features covered, hidden male ports and accommodates a variety of tubing sizes providing versatility in procedure and accessory setup
- **Variable Filter Life (up to 35 hours)**
Delivers unsurpassed value and reduces per procedure costs
- **Filter Tracking Technology**
Automatically tracks and records filter life through a patented technology
- **Simple LED Display**
Two button design (motor speed and on/standby) makes this system extremely easy to operate
- **Quiet Operation**



For more information download the product brochure from our website

\$2785

CMLV220

Buffalo EzLink® Automatic Activation Device

Automatic activation device, can be used to remotely turn the evacuator on and off, saving filter life.

\$99.00

CMLEZLINK01



Buffalo ViroVac® 3 Port Filter

Four layers of proven filtration gross particulate, ULPA media, activated carbon, and post filter.

3 port 6.4mm, 9.5mm and 22mm



\$418.40

CMLV353

Buffalo SnapEvac™ Electrosurgical Pencil Adaptor

Ergonomically designed smoke evacuation electrosurgical pencil adaptor. Transforms most hyfrecator hand-switch pencils and electrosurgical pencils into a surgical smoke evacuation device that captures smoke plume at the point of creation.

- **Snap and Evac** - Just snap in handpiece of your choice
- **Efficient Workflow** - Eliminates need for OR personnel to hold tubing
- **Transparent Capture Port** - Improves visibility of the surgical site
- **Cost-Effective** - Economical option for surgical smoke evacuation
- **Lightweight Tubing**



Box/20

\$470.05

CMLSNAPEVAC20

Buffalo Smoke Evacuation Tubing



22mm x 3.1m Non-Sterile

CMLVT1052410

Box/10

\$275.85

22mm x 1.8m Non-Sterile

CMLVT786

Box/6

\$126.90




What is surgical plume and what are the risks?

Surgical plume is a by-product from the use of electrosurgical instruments when cutting, vaporising or coagulating tissue during surgical procedures. It is generated when thermal damage of tissue causes a release of smoke. The smoke is often not visible and has an unpleasant odour.

Surgical plume can contain a mixture of potentially hazardous components including ultra-fine particulates, bacteria, viruses, cellular debris, noxious and toxic aerosols, gases, vapours and fumes. Human papillomavirus (HPV) and the human immunodeficiency virus (HIV) have also been found in the plume, although it has not been confirmed if that includes infective virus.

Whilst there is insufficient scientific research on adverse health effects from exposure to surgical plume, there is evidence that smoke from surgical plume is a potential health hazard, at least similar to smoking cigarettes, and research emphasises the basic principle that exposure to smoke of any kind should be minimised.

Reference: WorkSafe Victoria, updated on 22/06/2021, accessed 17/05/2022 <<https://www.worksafe.vic.gov.au/managing-surgical-plume-exposure-healthcare>>



Hyfrecator® 2000

The Hyfrecator has been the #1 product-of-choice by physicians in office-based electrosurgery for almost 70 years

Features only available in the Hyfrecator® system:

- High or Low Output Ports; allows you convenient switching back and forth between settings depending on your procedure. Save your most common procedural settings in both high and low mode. Precise output at low and high settings.
- Dual Microprocessors; one microprocessor controls the output and functionality of the system while the second terminates activation if a safety hazard is detected. This makes the Hyfreator one of the safest office-based electrosurgery units in the world.

1-year warranty

CONTENTS

- Power up/down handswitching pencil with 3m cable
- Wall mount kit
- 5 sharp and 5 blunt disposable tips
- Operator's manual, warranty card

The "Common Dermatologic Procedures" video is available online.
<http://www.vimeo.com/14252858> Password: Electrosurgery



\$2195
CML7900230

Electrolase disposable electrodes sharp tips

Sharp Non-Sterile	CML710012BX	Box/100	\$131.10
Sharp Sterile	CML71008BX	Box/50	\$131.10

Electrolase disposable electrodes blunt tips

Blunt Non-Sterile	CML710112BX	Box/100	\$131.10
Blunt Sterile	CML71018BX	Box/50	\$131.10

Disposable pencil sheaths

Pencil Sheath Non-sterile	CML779618BX	Box/100	\$64.70
Protects the pencil body from external contamination during procedures. 20cm			
Handpiece Sheath Sterile	CML779619CS	Box/25	\$230.65
Protects the handpiece body from external contamination during procedures. 92cm			




Sharp Tips - Ideal for pinpoint coagulation



Blunt Tips - Ideal for broad based coagulation



Reusable Hyfrecator Accessories



Bipolar Starter Pack: Most common starter combination = a + d + l

Bovie DERM 942 High Frequency Desiccator

The DERM 942 and its handpiece are designed to be a reliable and durable desiccator.

- 40 watt high frequency desiccator
- Adjustable in 1/10th watt increments in the first 10 watts, in monopolar and bipolar
- More accurate power with Universal Power Control 100-240 VAC
- No patient plate required
- Made in the USA
- 4-year warranty

Package Includes

- A942 Unit & Power Cord
- 3-Button Autoclavable Handpiece
- 5 non-sterile sharp dermal tips
- 2 sterile sharp dermal tips
- 5 non-sterile blunt dermal tips
- 2 sterile blunt dermal tips
- 2 non-sterile DERM-Elite™ Premium blunt dermal tips
- 2 sterile DERM-Elite Premium blunt dermal tips
- 2 sterile disposable handpiece sheath
- 2 non-sterile disposable handpiece sheath
- Wall mount kit



\$3125

DALA942

Dimensions: 18.8 x 22.8 x 10.5cm

Aaron Disposable Dermal Tips



Sharp, non-sterile	DALA804	Box/100	\$142.60
Blunt, non-sterile	DALA806	Box/100	\$142.60

Aaron Disposable Pencil Sheaths



Non-sterile	DALA910	Box/100	\$99
-------------	---------	---------	------

Disposable High Temp Fine Tip Cautery

- A highly portable and easy to use cauterising tool.
- To operate, depress the finger switch. Current flows through the burner tip.
- Fully disposable including the outer plastic handle, the batteries inside and the tip
- Individually sterile packed, 10 per box
- Includes 2 x AA batteries

Box/10	DALAA01	\$343.45
--------	---------	----------



Silver Nitrate Applicators

With well-established antiseptic and antimicrobial properties, silver nitrate from an applicator or pencil is applied topically usually in aqueous solution and used as a haemostatic, for pain relief in certain conditions and for the removal of small quantities of unwanted tissue.

- Simple, Effective, Economic
- 75% Silver Nitrate
- Single use 15cm Bendable Applicator
- 100 Applicators in a resealable tube
- Removal of non-genital warts (including verruca's)
- Removal of unwanted tissue
- Cauterisation



SUTSNA7484	
Pkt/100	\$90.95

