

CryoClear®
The all new, clearly better cryotherapy solution.
Multi-Use Disposable Cryosurgical Device

CryoClear uses Carbon Dioxide which freezes at -79°C to a depth of no more than 2mm into the skin, delivering a very safe and clinically effective treatment to superficial skin indications.

CryoClear's high efficacy - over 95% - combined with less destruction of healthy tissue makes it perfect for treating age spots, sun spots, and skin tags. Whether used as a stand-alone treatment, or as a combination treatment to chemical peels, facials, micro-needling or dermaplaining, CryoClear offers a better, faster cryotherapy solution with significant benefits for clients and estheticians.

Up to 300 secs of treatment time per device.

Indications for Use

- Acrochorda - Skin Tag
- Actinic Keratosis - Sun Spot
- Sola Lentigo - Liver or Age Spot

\$645
ECO15001001



Ideal for skin tags and sun spots.



CryoOmega®
Multi-Use Disposable Cryosurgical Device

Clinically effective direct spray of liquid nitrous oxide (-89°C) in a convenient and economical device.

- Delivers a continuous, direct spray of nitrous oxide (-89°C) to treatment site
- Lower operating temperature compared to canister-based systems, provides clinically effective freezing with pinpoint precision
- Simple, hassle-free operation, simply activate the internal cartridge once
- 30 to 45 five-second sprays per device
- Up to 300 secs of treatment time per device

\$485
ECO15001001



CryoOmega® is an ideal choice in portable cryosurgery for physician offices looking for convenience and great first- outcomes.



HistoFreezer® Portable Cryosurgical system

A professional-grade product that offers medical professionals an easy and effective method to treat nine different types of skin lesions, including:

- Acrochordon (Skin Tags)
- Actinic Keratosis (Facial & Non-Facial)
- Condyloma Acuminata (Genital Warts)
- Lentigo (Facial & Non-Facial)
- Molluscum Contagiosum (MCV)
- Seborrheic Keratosis
- Verruca Plana (Flat Warts)
- Verruca Plantaris (Plantar Warts)
- Verruca Vulgaris (Common Warts)

ECO KIT Contents:
2 x 80ml Canisters of HistoFreezer Cryogen
24 x 2mm Applicators
36 x 5mm Applicators

\$144
ECO10010412



Cryo Professional Portable Cryosurgical system

The Cryo Professional can be used to treat various types of benign skin lesions on patients over 4 years old. Examples of the types of skin lesions treatable with the Cryo Professional include:

- Genital warts
- Molluscum Contagiosum
- Fibroma, Acrochordon (skin tags)
- Verruca (Plantaris, Vulgaris, Plana)
- Seborrheic Keratosis (SKs)
- Aktinica Keratosis (AKs), both in the face and on other body parts
- Lentigos, on the face and other body parts

\$132
CRIUTM0171

Contents:
1 x 170ml aerosol dispenser mix
30 x 2mm foam sticks
30 x 5mm foam sticks



CRYOSUCCESS Liquidfreezing

CryoSuccess with 1 and 4 mm tip

The unique method with direct contact between the liquid and the skin allowing a constant penetrating temperature of -89°C

- Medical Grade Gas
- Low treatment costs
- Simple and secure handling
- Autoclavable tips
- Extremely precise application, with a constant penetrating temperature of -89° C
- Approximately 300 seconds of treatment time per cartridge
- Cartridge can be removed and reattached without wasting gas
- Gas Cartridges have no shelf life or evaporation
- Excellent cosmetic effects; almost no scar formation
- Non bleeding treatment, no anaesthesia required
- Practically painless, therefore suitable for children and sensitive patients

\$950

STMCSUC01

Kit includes:

- 1 x CryoSuccess
- 2 x 23.5g N2O Cartridge
- 1 x 1mm Tip
- 1x 4mm Tip



Visit our website to watch the product video

Cryosuccess 1mm Tip

- For fine work
- Autoclavable

\$335

STMSTS01



Cryosuccess 4mm Tip

- For larger areas
- Autoclavable

\$390

STMDTS04



Large 23.5g Refill Cartridge

\$104.90

STMCA2310



Walmay Cryospray

- Made in Australia
- 0.35L Flask capacity
- Fixed stainless steel probe with high tensile strength
- Regulator can be positioned on the right of left to suit the operator
- Set of four cone-shaped spray tips for accurate visual aiming during treatment
- One curved spray tip for body orifice surface treatment

WALWM227

\$895



Premier Nitrospray Plus

Premier Nitrospray Plus Cryosurgical Unit constructed of a double-walled vacuum stainless steel container (canister) to hold liquid nitrogen. The outflow track extension port, located on the top of the instrument has a luer-lock connector to easily engage/disengage spray tips and other attachments as required.

Cryosurgical Unit

- Easy to use, comfortable to hold
- Unobstructed view of treatment site
- 6-8 hours static holding time
- Luer-lock spray tip connector
- 3-year warranty
- Made in the USA

Includes:

- Cryoflask
- Set of 5 luer-lock spray tips
- Spray tip protector.

Holding Times:

- Static 6 – 8 hours
- Intermittent Use 3 – 6 hours of patient treatment time



Nitrospray Lite 300ml

PACLNSLITEPLUS

Nitrospray Plus 475ml

PACLNSPLUS

\$899

Nitrospray Accessories



Tip Set
Set of tips: Sizes 16-20 Gauge
\$89
PACLNSTSET



Tip Extension Bent
\$35
PACLNSTEB



Tip Extension Rigid
\$99
PACLNSNER



Spray Tip Protector
\$39
PACLNSTP



Luer Lock Tip
\$19.95
16G PACLNST16G
17G PACLNST17G
18G PACLNST18G
19G PACLNST19G
20G PACLNST20G

Cryotherapy Starter Kit

Nitrospray Cryospray	\$899.00
Dewar 20 ltr	\$1,525.00
Extractor 20 ltr	\$181.00
Ladle 20 ltr	\$133.00
Measuring rod	\$29.05
Polar gloves	\$42.00
Total	\$2,809.05

Total Package **\$2459**



Dewar

Liquid Nitrogen Container

- Made of high-strength aluminium alloy
- Adopting advanced powder coating technology
- Each unit include leather bag and ladle
- Designed for large capacity storage and low liquid nitrogen consumption
- Safety lock lid
- 5 year vacuum warranty
- 2 year container warranty

Model	LNC-10-50	LNC-20-50	LNC-30-50
Capacity	10L	20L	30L
Weight empty (kg)	6.0	9.5	13
Caliber (mm)	50	50	50
Outer diameter (mm)	300	394	462
Height (mm)	552	672	706
Static evaporation loss (L/d)	0.12	0.12	0.12
Static holding time (days)	86	168	254
Code	PACLNC10L	PACLNC20L	PACLNC30L



10 ltr **PACLNC10L** **\$1357**
20 ltr **PACLNC20L** **\$1525**
30 ltr **PACLNC30L** **\$1747**



each dewar includes ladle

Gloves Liquid Nitrogen

- Gloves for cryogenic applications
- Length of glove is approx. 35cm
- Can bare temperatures from -180°C to -250°C
- External waterproof material, inner material has a high quality cold-proof sponge and lining



PACLNCGLV **\$42.00**

ActivArmr® 23-700 Glove

- Gloves that shed water, resist salts and chemicals
- Rough-finish PVC coating provides gripping power on wet and slippery surfaces



AN23700
\$32.60

Extractor 20/30 Ltr

Liquid Nitrogen extractor for filling cryospray units.

PACLNCXT **\$181**



We recommend having a ladle as backup when purchasing an extractor.

Ladles

10 Ltr
PACLNL10L **\$129**
46cm for 10 ltr ladle

20 Ltr
PACLNL20L **\$133**
56cm for 20/30 ltr ladle



Measuring rod

The measuring rod will clearly show the LN2 level of a dewar in inches and centimetres to 92cm/36".

TAYR0338C11 **\$29.05**



MEET OUR TEAM: NICK

With over 25 years of experience in the medical industry our East Melbourne Sales Rep Nick is a true asset to our team. A people person, Nick loves engaging with customers and sharing his expertise and detailed product knowledge.

For fun, he enjoys fishing, beer brewing and cooking, often choosing recipes that go well with finds from his wine cellar.

What sport do you love watching? AFL (Essendon's Mighty Bombers – established in 1872). Our VIC team is keeping up the friendly rivalry!

What is your biggest pet peeve? Teabags in the sink.

Have you been at a show recently? The Point at Portsea - like Big Day Out for old people.

