

Setting up a new practice, upgrading or expanding?

We're the right team for the job. With the largest GP Sales Team and our extensive experience, we look after your team from beginning to end. No set-up is too big or too small for Team and BP Medical.

Need a Quote Today?

Give us a call or use our online set-up calculator for a quick estimate.



HVF128 TENS Therapy Device

Engage TENS technology for drug-free treatment.

- Block the pain messages to the brain
- Trigger the release of endorphins (natural pain killers)
- Improve blood circulation (as result of repeated muscle contraction and muscle relaxation)

Adapt your treatment according to your preference and pain symptoms.

- Hyper Mode - For relief of stiffness with gradual increase of intensity during treatment
- Soft Mode - For relief with soft stimulation
- Repeat Mode - Repeats your choice of stimulation
- Point and Wide Mode - Treats a small or wide area respectively

Washable & Durable Long-Life Pads, reusable up to 150 times



OMRON

\$229\*  
DAVHV128

Your patients can buy all these items via our online shop  
www.teammed.com.au

Pro Tens Machine

A regular TENS with adjustable pulse width, pulse rate, intensity, modes and an accurate timer of 15, 30 and 60 minutes, will stop output at presenting time.


- Channel: Dual, isolated between channels  
Pulse Amplitude: 0-80mA peak to peak into 500 ohm load  
Pulse Rate: Adjustable, from 2Hz ~ 150Hz  
Pulse Width: Adjustable, from 30uS ~ 260uS  
Mode: B, N, M Timer: Selectable, 15, 30, 60 minutes (or 15, 30, constant)  
Wave Form: Asymmetrical Bi-Phasic Square Pulse  
Voltage: 0-40V peak into 500ohm load



ALLPROTENS \$54.90\* Supplied with Carry Case, 2 x TENS Leads, Operating Manual, 1 x 9V Battery, 1 x Pack of Self-Adhesive Electrodes allcare

Tens Electrodes

- 5 x 5cm square ALLF35050 Pkt/4 \$12.45\*  
5 x 9cm rectangle ALLF35090 Pkt/4 \$13.95\*  
5cm round ALLF350R Pkt/4 \$11.95\*



Perfect MamaTENS

Perfect MamaTENS uses the latest technology in TENS to offer drug-free, more natural pain relief during labour.

It is designed to be portable and simple to use and to provide an alternative to more intrusive methods of pain relief during childbirth

- Portable– can be used at home when labour starts and then taken into the hospital
- Easy-to-use controls– allowing the Mother to remain in control of her own pain relief
- Dual channel stimulation
- Specifically designed programs– 3 pre-set programs designed to help counter labour pain at every stage
- 4 large electrode pads– for maximum nerve coverage
- Up to 60 steps of intensity– for gentle increases in stimulation
- Boost Button– to provide additional pain relief during contractions
- Active Power Ramp Mode– provides “brief intense” pain relief during the final contractions
- Safety feature– returns intensity to ‘0’ if there is a break in the circuit connection

DAVKPMT \$194.95

Your patients can buy all these items through our online shop [www.teammed.com.au](http://www.teammed.com.au)

Elise Pelvic Floor Exerciser

The Elise Pelvic Floor Exerciser is an easy-to-use, flexible pelvic floor exerciser with a range of clinically tested programs – providing relief from stress, urge and mixed incontinence. Pelvic floor exercisers are one of the most effective forms of continence management providing a discreet, affordable and long-term solution for incontinence sufferers. The Elise comes with a Vaginal Probe but can be used with an Anal Probe (sold separately) to treat Faecal Incontinence and Erectile Dysfunction

- One Touch Memory – remembers the last program, strength and setting used
- LED Backlit Screen – easy-to-use, modern, colour screen for clearer visibility
- Comfortable Stimulation – gentle, comfortable stimulation with fine-tune adjustment settings for different levels of intensity
- Rechargeable Li-ion Battery – lightweight & compact with an external charger
- Safety Duration Override – unit will turn off after 20 minutes to ensure the pelvic floor muscle is not overworked

Advantages of using the Elise:

- It can reduce leakage – not simply contain it. Most women see significant results in as little as 3 weeks.
- It is drug-free with no side effects.
- It is easy to use with preset programmes to suit your needs.
- It is discreet and can be used at home during your own time.
- It may help to avoid surgery

DAVKELISE \$199



Welcare Stay-Dry Bedwetting Alarm

BWA100 uses behavioural techniques to train your child out of bedwetting (Nocturnal Enuresis).

- Bedwetting alarms have been found to be one of the most effective methods of training a child out of bedwetting
- When the child wets the bed, a moisture-sensitive alarm will sound and/or vibrate
- Helps the child to recognise and respond to the sensation of a full bladder while asleep
- User Friendly, Light & Compact
- Adjustable Arm Band
- Sound & Vibration Mode
- Vibration Only Mode
- LED Light Indicator
- USB Rechargeable Battery with Charging Cable
- Progress Chart & Stickers

\$86.75  
DAVWBA100



Hwato Acupuncture needles

Hwato Acupuncture Needles are China's finest quality traditional needles. Established 140 years ago, the Hwato factory is the world's oldest, largest and most respected manufacturer of hand-made needles (metal handle needles).

20 x 30mm tube	ACUHT2030
20 x 40mm tube	ACUHT2040
22 x 13mm tube	ACUHT2213
22 x 30mm tube	ACUHT2230
25 x 30mm tube	ACUHT2530
25 x 40mm tube	ACUHT2540
25 x 50mm tube	ACUHT2550
25 x 75mm tube	ACUHT2575
30 x 30mm tube	ACUHT3030
30 x 40mm tube	ACUHT3040
30 x 50mm tube	ACUHT3050
30 x 60mm tube	ACUHT3060

Box/100	\$14.70
5+	\$13.50
10+	\$12.15



Seirin Acupuncture Needles (J-Type Plastic Handle)

The SEIRIN J Type is made with a new technology that allows for smooth painless insertion using less pressure. Features an easy grip, colour coded, lightweight 2.0 x 20mm plastic handle with a soft-touch top.

20 x 30mm tube	ACUSJ2030
20 x 40mm tube	ACUSJ2040
25 x 30mm tube	ACUSJ2530
25 x 40mm tube	ACUSJ2540
25 x 50mm tube	ACUSJ2550
30 x 30mm tube	ACUSJ3030
30 x 40mm tube	ACUSJ3040

Box/100	\$35.20
5+	\$32.00
10+	\$30.75



Quick Ultrasound Gel

- Acoustically correct; for the broad range of frequencies used
- Completely aqueous; will not stain clothing or damage transducers
- Water-soluble and odourless
- Hypoallergenic, bacteriostatic, non-sensitising and non-irritating



Quick-Eco Gel Blue

blue colour

250ml	BIO9980020	\$9.10
5Ltr	BIO9980021	\$35.00



Quick-Cristal Gel Clear

no colour, totally transparent

250ml	BIO9980025	\$9.10
5Ltr	BIO9980022	\$34.50



Aquasonic 100  
Ultrasound Gel

Most widely used gel for diagnostic and therapeutic ultrasound

250ml	BIBUS0108	\$16.35
5Ltr	BIBUS0150	\$86.30



Tensive® Conductive  
Adhesive Gel

Safe, non-flammable, odourless, conductive adhesive gel. Recommended for TENS, EMG, ECG and other electromedical procedures.

50G Tube	AUSDITG50	\$26.25
----------	-----------	---------



Can't find the right product? Need more information?  
Visit our website [www.teammed.com.au](http://www.teammed.com.au) or call our team. We will look after you.



MEET OUR TEAM: DEE  
CUSTOMER SERVICE - QLD

When it comes to customer service in the medical supplies industry, it would be hard to find someone with more experience than Dee. Dee joined the team 4 years ago and has been in the industry for a total of 18 years. Her favourite part of working at Team Medical is *"The people, Team Medical is a great place to work with plenty of support."*

Like her favourite animal the tiger, Dee is tenacious yet calm by nature.

Next Holiday? Melbourne

Pet Peeve? Toilet paper hanging wrong way (paper should face outward)

Favourite Food? Chinese



MEET OUR TEAM: JULIE  
CUSTOMER SERVICE - WA

Julie has been part of our BP and Team Medical team for 15 years, helping customers with orders and products. When not helping customers, Julie loves watching AFL football (West Coast Eagles) and MAFS.

She recently came back from a trip to Canada. *"It is incredibly beautiful. The Rocky Mountaineer was definitely a highlight!"*

Favourite Food? Fish and chips

A fact about you that not many people know? I did salsa dancing for a couple of years





S+M Instant Ice Pack

230 x 130mm

AXP10130331  
Each \$3.05  
12+ \$2.20



Instant Ice Pack

230 x 150mm

LRINT11069  
Each \$2.00 40+ \$1.85



S+M Hot/Cold Pack

Medium (110 x 220mm)  
AXP10130311

\$5.55

Large (135 x 260mm)  
AXP10130312

\$6.10  
Heat: Microwavable or hot water



Bodichek Premium  
Hot/Cold Pack  
(with towel bag)

Medium (125 x 270mm)  
AXP62630353

\$11.30

Large (175 x 270mm)  
AXP62630352

\$19.70

Heat: Microwavable or hot water



Bodichek Reusable  
Hot/Cold Gel Pack  
(clear plastic)

200 x 115mm

Heat: hot water only

AXP62630311  
\$6.95



Hot/Cold Gel Pack Reusable

PVC  
100 x 150mm



LRINTGEL3057  
\$3.45  
Heat: Microwavable or hot water

Hotteeze Hand Warmers

- Keeps hands and fingers warm
- Self-heating - air activated
- Non-stick
- Designed to hold in your hand, keep in your pockets, or tuck into your gloves
- Up to 10 hours of heat
- Max temp 63° C

Pad size – 100mm x 70mm  
DAVHHW10

Pkt/10 \$11.75



Hotteeze for Feet

- Keeps your toes nice and warm
- Self-heating - air activated
- Non-stick
- Designed to fit into men's and women's shoes
- Up to 6 hours of heat
- Max temp 39° C

Pad size – 96mm x 70mm  
DAVHFP5

Pkt/5 \$11.75



Hotteeze Heat Pads

- Keeps hands and fingers warm
- Self-heating - air activated
- Non-stick
- Designed to stick over clothes or to a bodigrip tube
- Up to 10 hours of heat
- Max temp 63° C

Pad size – 132mm x 100mm  
DAVHP10

Pkt/10 \$27.35



ALL-Band 1.5m Exercise Band

All-Band Exercise Band is used in rehabilitation and exercise programs. Available in a handy 1.5 metre length single pack or 25m rolls.

Resistance		1.5m x 14.5cm		25m x 14.5cm	
Light	Yellow	ALL04002401	\$9.00	ALLACBAND25YE	\$43.30
Medium	Red	ALL04002402	\$10.70	ALLACBAND25R	\$47.25
Firm	Green	ALL04002403	\$11.85	ALLACBAND25G	\$59.95
Extra Firm	Blue	ALL04002404	\$10.15	ALLACBANDBL	\$58.50
X-Extra Firm	Black	ALL04002405	\$10.15	ALLACLBANDBK	\$65.80




Exercise and Posture Ball

A great tool to enhance:  
Posture, Performance, Flexibility, Strength and Muscle tone

- A durable, flexible and firm ball
- Feature a non slip surface with ribbing for extra security
- An anti-burst ball rated at 550kg


Red	55cm	ALL03001512	\$33.55
Orange	60cm	ALL03001527	\$32.60
Green	65cm	ALL03001533	\$34.40
Blue	75cm	ALL03001554	\$39.10
Pearl	85cm	ALL03001568	\$40.95
Red	90cm	ALL03001569	\$52.45



Foam Rollers

- Made of high density eva foam
- Provides firm support for all body types
- Excellent for balance and motor skill development
- Also use for self massage, de-activating trigger points
- Use for soft tissue work, stretching, massage, stability exercises, core strengthening, injury prevention
- Strong enough to be used by large adults

Small Round	30L x 15Dcm	ALL700461R	\$21.60
Large Half Round	90L x 7.5Dcm	ALL700463H	\$32.40
Small Half Round	30L x 7.5Dcm	ALL700461H	\$18.20



Wheelchair  
Patient mover

An essential in every clinic

- 110kg user weight capacity
- Seat 46cm W x 41cm D
- Frame: Single cross brace, powder coating
- Front wheel: 8", PVC tires
- Rear wheel: 12" PU tires
- Armrest: Flip back with desk length
- Footrest: Swing away footrest

TMF18WCM \$349\*



Wheelchair - Standard

A versatile and practical wheelchair

- 110kg user weight capacity
- Seat 45cm W x 40cm D
- Frame: Single cross brace, powder coating
- Front wheel: 8"x/6", rubber tires
- Rear wheel: 24" spoke wheel, pneumatic tires
- Armrest: Flip back with desk length
- Footrest: Swing away footrest

TMF18WCS \$329\*



Crutches

Forearm Medium	890-1200mm	AUS4010
Forearm Large	890-1370mm	AUS4020

\$56.50\*

Underarm Small	1400-1600mm	AUS4050
Underarm Medium	1570-1780mm	AUS4051
Underarm Large	1570-1980mm	AUS4052

\$59.50\*



Wheelchair - Bariatric

High Capacity wheelchair

- 160kg user weight capacity
- Seat 55cm W x 40cm D
- Frame: Single cross brace, powder coating
- Front wheel: 8"x2", PU tires
- Rear wheel: 24" spoke wheel, PU tires
- Armrest: Flip back with desk length
- Footrest: Swing away footrest
- Seat and backrest: Nylon with sponge

TMF22WCB \$459\*

