



Fingertip Pulse Oximeter

- Splash-proof and drop resistant
- OLED display with SpO2 and pulse rate
- Spot check and continuous measure model
- Automatic power on/off
- Alarm with low battery indicator
- Accurate and sensitive

PACPMFPO

\$80



Fingertip Pulse Oximeter - PC60B1

A very strong, reliable and durable Pulse Oximeter. Accurately measures SpO2 value and Pulse Rate value.

- Dual colour OLED display of SpO2, PR, PI, waveform and pulse bar
- Innovative 4-direction display
- Water resistant and drop resistant
- Lightweight and portable
- Automatic power on/off
- Audible & visible alarm function
- More than 30 hours working time
- Includes neck lanyard, batteries and carry pouch
- 12 month warranty

MDQPC60B1BL Blue
 MDQPC60B1O Orange
 MDQPC60B1Y Yellow **\$80**



1YR warranty

TOP SELLER

Paediatric Fingertip Pulse Oximeter - PC-60D2

- OLED display of SpO2, PR, PI, Pulse Bar and waveform
- Special design for paediatric
- Reliable and accurate performance against low perfusion
- 4-direction display with Plethysmogram
- Audible & visible alarm function
- Includes neck lanyard, batteries and carry pouch
- 12 month warranty


MDQPC60D2P **\$100**



Finger Pulse Oximeter with Paediatric probe (PC-60E-plus)

- SpO2 value and pulse rate value
- Bright colour LED and multi-directional display
- Automatic display direction adjustment
- Over limit indication
- External paediatric probe
- Neonatal probe optional
- Includes neck lanyard
- 12 month warranty

\$320
 MDQPC60EPLUS



Optional probes

- MDQ15044048 Neonate reusable Y sensor with two disposable wraps
- MDQ15040014 Paediatric long term (reusable) rubber probe (included with PC-60E)



Handheld Pulse Oximeter - Model PC-66B



- High resolution LCD to display SpO2, Pulse Rate and plethysmogram
- Includes adult and paediatric probe
- Optional infant and neonate probes
- Stored data can be uploaded to PC for further review
- SpO2 trend curves through data manager software
- Up to 384 hours of oximetry data storage
- Records data once every 1/2/4/8 seconds
- Over limit indication
- Includes desk stand
- 12 month warranty



Light, portable and very user-friendly



\$425

MDQPC66B

Handheld Pulse Oximeter - SP-20



Ideal for Adult, Paediatric and Neonate

- 3.5" auto rotation display for SpO2, PI, Plethysmogram, PR
- Patented MOVEOXY™ SpO2 technology with artefact removal and anti-motion performance
- Supplied with Adult and Paediatric Probe (neonate optional)
- Up to 20 hours continuous monitoring with built-in lithium rechargeable battery
- Adjustable alarm limits, pitch tone level, key tone and brightness
- Up to 1000 hours SpO2 data storage
- Stored data can be uploaded to PC for further review
- Supplied with charging cradle, USB power cable and instructions



\$475


MDQSP20



Optional Probes for Creative Medical Oximeters



	PC66B	SP20	SP20 Version
1. Adult probe – clip D type	MDQ15040009	MDQ15040022 MDQ15040091	Models Pre June 2020 Models Post June 2020
2. Paediatric Probe – clip D type	MDQ15040015	MDQ15040055	
3. Infant Probe –Y type	MDQ15040084	MDQ15040050	
4. Neonate velcro wrap	MDQ15040072	MDQ15040074 MDQ15040092	Models Pre June 2020 Models Post June 2020
5. Oximeter extension cable	MDQ29090004	not required	



Happy Team Medical Customer

“Words to describe Team Medical? Eager to please, efficient, knowledgeable, supportive, negotiable, systems, approachable.”

Long-standing NSW-based Customer



Simple fast accurate vital signs at a great price

Vital Signs Monitor: SpO2, NIBP, PR

Model PC-900PRO

This small multi-parameter patient monitor is light-weight and user friendly. It monitors the vital physiological parameters including Non-Invasive Blood Pressure (NIBP), Pulse Oxygen saturation (SpO2) and Pulse Rate (PR).

- 3.5" LCD display, easy to read data, alarm events and menus
- Standard with SpO2, PR and NIBP
- Storage for 12000 groups NIBP and 2000 groups SpO2 measurement
- Auto/Manual/Continuous BP mode
- Inbuilt rechargeable lithium battery with 6-hour continuous monitoring
- Nurse call function is available
- USB connectivity for data uploading and software upgrading
- Supplied with Adult & Large Adult NIBP Cuff
- 3-year warranty

\$1400

MDQPC900P



Mobile trolley for vital signs monitors

Suitable for PC-900PRO and PC-3000

MDQCGX

\$535




Optional Probes & Cuffs for Creative Medical Vital Signs monitors

Probes & cuffs for PC-900PRO, PC-3000, PC9000-C Vital Signs Monitors		
1.	Adult probe – clip D type	MDQ15044074
2.	Paediatric Clip	MDQ15040021
3.	Neonatal reusable Y type probe	MDQ15040019

Probes & Cuffs for K-Series Vital Signs Monitors		
4.	Adult probe – clip D type	MDQ15040077
5.	Paediatric Clip	MDQ15040097
6.	Paediatric Clip Soft	MDQ15040098



MANAGEMENT -NSW



Meet our Team: Matt

Matt has seen Team Medical grow from the ground up. Matt became our first Warehouse Team Member when he decided to join the company 'just part-time' to help out his aunt - Marijke. Thankfully, Matt has never looked back; today, he is our Procurement Director, an essential part of our day-to-day operations, and still a much-loved part of the Team.

Outside of work, Matt is also a Team Player, playing Football (Soccer) with his brother Tom (General Manager) and lifelong friend Jake (Operations Manager); once a keen tennis player, Matt's now trying his luck with paddle tennis. We can't forget to mention his beloved dog Boston, a Rhodesian Ridgeback bringing joy and trouble, who also likes to pop up in the occasional video meeting.

- Other hobbies?** Travelling taking drone photos
- Favourite Food?** Hot Chips
- What is one of your pet peeves?** Bad Drivers



Baxter

Welch Allyn

Spot Device & Monitor

Spots Vital Signs 4400 Device

BP, SpO2, Temperature
Simple, Accurate, Connected

The Welch Allyn Spot Vital Signs® 4400 Device gives you a simple, all-in-one touchscreen solution for measuring patient vital signs.

- SureBP technology can record blood pressure in just 15 seconds
- BP averaging helps ensure accurate hypertension detection
- SureTemp Plus oral/axillary thermometry
- 7" color touchscreen for a simplified workflow
- EMR connectivity helps your care team avoid documentation errors
- Includes one adult blood pressure cuff

WA44WT6

\$3490

Please include your option of stand

Desk Stand

WA4400DST


\$355



Mobile Stand

WA4400MBS

\$795



Connex Spot Monitor

Make smarter, more informed clinical decisions.

Discover comprehensive monitoring—from basic vitals capture to early warning scoring and interval monitoring—in a simple, secure, connected solution.

- Enhance productivity and patient outcomes by capturing a full set of vital signs in less than a minute
- Easy-to-use, vivid touchscreen display
- SureBP technology can record blood pressure in just 15 seconds
- Help ensure accurate hypertension detection with BP averaging

Connex Spot Vital Signs Monitor
with SureBP®, Nonin SpO2 and Braun ThermoScan® PRO 6000 Thermometry

WA71WE6

\$3665

Classic Mobile Stand

WA7000MS3

\$790



Mobile WorkStand

WA7000MWS

\$935

Baxter

Multi-Parameter Patient Monitor: ECG, RESP, SpO2, NIBP, PR, TEMP



Model PC-3000

Features:

- 7" high resolution color TFT display, 5-channel waveform display
- Standard parameters: ECG, RESP, SpO2, NIBP, PR, TEMP
- Arrhythmia analysis and ST segment analysis
- Up to 12000 groups parameter data & 72 hour ECG waveform storage and recall
- 1000 hour graphic and tabular trends for all parameters
- 7-channel ECG waveforms display simultaneously
- Protection against interference from defibrillator and electrosurgical knife
- Adult / Paediatric / Neonatal measurement modes
- Visual and audible alarms
- Built-in rechargeable battery
- 2-year warranty



2YR
warranty



Mobile trolley for vital signs monitors

Suitable for PC-900PRO and PC-3000

\$1720

MDQPC3000

MDQCGX **\$535**

AiView V10/V12 Patient Monitor

- Standard configuration: 5 lead ECG, Respiration, SpO2, PR, NIBP, Temperature
- Optional 2-IBP, EtCO2, Nellcor SpO2, SunTech NIBP, 12-lead ECG, WIFI
- Full-screen design, capacitive projection touch screen, support gesture operation
- Updated UI interface, preset H1 menu, optimized user experience, IPX2 waterproof
- 3/5 electrodes ECG collection, ST wave filter mode, myocardial ischemia assessment, optional embedded AI chip, multi-leads simultaneous analysis.

Allows for multiple application scenarios. Experience a new monitoring concept



V10 - 11.6 inch TFT display

MDQAIIEWV10

\$2695*



V12 - 13.3 inch TFT display

MDQAIIEWV12

\$3235*

Explore

Take a closer look at all the features of the AI View patient monitor.



For more information download the product brochure from our website



Keep up to date with us on Social Media

For stories on our team, updates on events, special offers and new products specials follow us on Facebook, LinkedIn YouTube and Instagram.

