

NATIONWIDE FIELD SUPPORT

Providing high quality service Australia wide



Annual Programmed Maintenance

To ensure efficient and reliable operation of your vaccine fridge it is essential that your fridge undergoes regular maintenance. Avem Quirks offer maintenance plans for medical refrigeration equipment to keep your facility running smoothly. All maintenance is performed as per the Strive for 5 recommendations.

Preventative maintenance includes:

- Check fridge seals
- Checking and cleaning of compressor and motor
- Testing and cleaning of drains
- Check of temperature setting - minimum/maximum temperature and data logger
- Cleaning and checking of evaporation tray
- 12 months access to the AQ Cloud monitoring portal
- Calibration of AQ box probes
- Test and Tag of unit
- Certificates and documents for your records




Book your service by calling 1800 658 120 or visiting www.avemquirks.com.au

Strive for 5

For advice and information on vaccine storage management refer to the Vaccine Storage Guidelines; Search online "Strive for 5"

Strive for 5 - Strive for 5 degrees Celsius (°C) – that is, the point midway between +2°C and +8°C which is the temperature range recommended for vaccine storage.





Medisafe Plus Vaccine Fridges

For 10 Years the Medisafe Plus range has been the vaccine fridge of choice for Medical Centres. You won't find a fridge better equipped for the storage of your vaccines.

The Medisafe Plus medical fridges are purpose-built for vaccine storage and follow the "Strive for 5" National Vaccine standard.

All Medisafe fridges come equipped with an AQ Box with Medisafe Plus cloud monitoring. The cloud monitoring acts as an automated data logger which allows you to view and print your temperature chart with min and max data in real time on any web enabled smartphone, tablet, or computer.

- Fully commercial construction
- Superior quality: no domestic fridge components used
- Highly accurate electronic controls (incorporates dual temperature probes and thermostat)
- Digital Temperature control system pre-configured to hold temperature between 2°C to 8°C (approved QCPP specification)
- Equipped with temperature alarms, AQ Box and electronic thermostat
- Adjustable shelves
- 2 Year Warranty



Additional shelves available



Model	Capacity	Dimensions	Power	Code	Price
HB80	80 Litre	50W x 56D x 83H cm	3.08 kWh/24h	QUIFKG080	\$2595
311 Model	306 Litre	59.5W x 59.5D x 163H cm	3.69 kWh/24h	QUIFKG311	\$3395
371 Model	381 Litre	59.5W x 59.5D x 193H cm	4.17 kWh/24h	QUIFKG371	\$3695
650 Model	610 Litre	74.5W x 74.5D x 220H cm	5.90 kWh/24h	QUIAKG650	\$6995
1000 Model	980 Litre	113W x 74.5D x 220H cm	8.50 kWh/24h	QUITME1000	\$8995





Invest in an UPS system:
Save time
Save money
Provide peace of mind

UPS and Battery Back Up Pack
(Uninterrupted Power Supply)

The Medisafe Plus UPS and Battery Pack is specifically designed for Medisafe Plus Medical refrigeration providing extended power back-up during a power outage.

Rated at 2kVA, this UPS provides sufficient capacity to cope with the draw required to power a 2 door vaccine fridge. The UPS will also regulate power surges, circuit shorts and spikes while providing a consistent power output of 230v at all times.

UPS and Battery pack provides up to 8 hours run time when attached to a medium sized refrigerated cabinet (average 300 L). Units can be combined with an additional battery pack unit (2+) to increase run times of up to 7-8 additional hours per unit.

- Each UPS and Battery Pack comes with a market leading 3 year warranty
- Protection against power outages
- Protection against power variances
- UPS and battery pack is vertically stackable

Specifications	
Dimensions (W x D x H)	UPS: 192 x 318 x 428mm Battery: 192 x 318 x 428mm
Weight	26Kg
Capacity	2kVA
Power Requirement	10amp GPA
Input Voltage	110VAC-288VAC
Input Frequency	40Hz-70Hz
Output Voltage	230V
Phase	Single Phase

	Code	Price
UPS + Battery Back Up Pack	QUI1P2000TSB	\$3095
Additional Battery Pack	QUIESBP3RV	\$1595



*Run times are based on a new unit and will reduce with usage.



AQ Box

Wireless temperature data logging with
automatic cloud-based reporting

- Urgent notification via text message and email of any temperature faults, alerts continue until the issue is resolved
- Equipped with two temperature probes providing accurate measurements of both air and simulated product temperature for more accurate monitoring
- The AQ Box portal allows access to easily view, download or print reports with date and time information, anytime, anywhere on any device with an internet connection
- Built-in 8 hour battery back-up ensures your AQ box will continue to monitor and provide alerts when a power failure occurs

The AQ box comes standard with all Medisafe Plus Vaccine Fridges or can be fitted to your existing vaccine storage unit.

12 Months Subscription included
with Medisafe Plus Vaccine Fridge





AQ Medical Vaccine Fridge

AQ Medical Vaccine Fridges are a robust, cost effective solution for all your vaccine refrigeration needs

- Self closing door
- Temperature keeps from 2°c - 8°c accurately and evenly
- High-low temperature audible and visual alarm
- Factory installed lock ensures security
- Electronic temperature controller and display
- 2 year warranty
- National service support



Model	Capacity	Dimensions	Power Supply	Code	Price
AQM 110	145 Litre	59.5W x 59.5D x 82H cm	10amp	QUIAQM110	\$1650
AQM 280	306 Litre	58W x 58.4D x 162.5H cm	10amp	QUIAQM280	\$1895
AQM 350	381 Litre	58W x 58.4D x 190H cm	10amp	QUIAQM350	\$3080



ICS PACIFIC Pharma 1000 Vaccine Fridge 145L

Under counter or bench top vaccine fridge. Fully compliant with Quality Care Pharmacy Program (QCPP) Cold Chain Standards.

- Easy-use High / Low / Reset temperature display
 - Digital temperature control & display
 - Temperature range: +3°C to + 7°C
 - No freeze fully protected internal rear wall
 - Audible and visual alarm with external alarm facility
 - Self-closing lockable door with internal LED lighting
 - Virtually maintenance-free
 - Fully ventilated cooling system
 - Automatic defrost
 - Data logger access valve
- Insulation - 50mm
 - 59.5W x 64D x 83Hcm
 - 4 adjustable shelves plus raised base shelf
 - Effective capacity - 145 litres
 - Energy consumption - 0.7 kWh/24h
 - Optional door alarm available
 - Made in Denmark

ICS1000GD

\$2295



Emergency Vaccine Esky 35L

For emergency storage of vaccines and cold chain products.

- More insulation and thicker walls than most fibreglass eskies for longer lasting coolness
- Made from roto moulded food grade U.V stabilised polyethylene - will not chip, crack or break
- Equipped with advanced refrigeration grade Iso and Polyol insulation
- Compact esky size for storage in your practice
- Dimensions: 46.8(L) x 38.6(W) x 38.9cm(H)
- Made in Australia
- 3 Year Warranty

TEAMT35EVE

\$245



CleverLogger

Starter Kit - Gateway + Wireless logger

CleverLogger is designed for non-technical users. It measures the temperature and saves the results to the Internet every few minutes.

It is browser-based, so there is nothing to install or update. Clever Logger creates a live log of temperatures that’s accessible from anywhere.

It can create PDF or Excel reports and email them to you automatically every day/week/month.

If something goes wrong, Clever Logger will send you an immediate alert.



ONSCLK01 \$317.85

CLEVER
LOGGER

LogTag Data Loggers

LogTag

LogTag with Internal Sensor	A	ONSLOGTAG	\$79.00
LogTag with Internal Sensor and Built-In USB Plug	B	ONSUTRIX16	\$152.00
LogTag with Internal Sensor and Temperature Display	C	ONSLOGDISP	\$165.00
Logtag with External Probe ST100S 30mm tip ø5mm 1.5m	D	ONSLOGEXS15	\$162.40
Logtag with External Probe ST100J 40mm tip ø3.2mm 1.5m		ONSLOGEXJ15	
LogTag Replacement - No Probe (logger only; probe required not included)	E	ONSLOGEXNP	\$118.40
Logtag with Display and External Probe	F	ONDISPEXSET	\$227.50
Logtag LogTag Replacement with Display - No Probe (logger only; probe required not included)	G	ONSLOGDISPEXNP	\$189.00
External Probe Only 40mm ø3.2mm for LogTag	H	ONST100J	\$66.00
LogTag software/reader kit (cradle and software; logger not included) You will need this reader to use the LogTag (Except USB LogTag).	I	ONSLTIKIT	\$131.10



Max Min
Thermometer



The unit incorporates two temperature sensors, a remote water resistant probe with one metre PVC lead for the appliance temperature and an internal sensor for the room temperature. The external probe can be mounted onto the fridge internal wall using the suction pad or mounting bracket supplied.

- Range of -49.9 to 69.9°C
- Large LCD display
- Max/min memory function to record the highest and lowest temperatures
- Programmable high/low audible alarm

DIGMEDIT
\$40.95

Correct Storage and Transport is Paramount.

Because, even one broken link in the chain can mean your vaccines aren't viable.

At Team Medical, we adhere to the 'Strive for 5' National Vaccine Storage Guidelines. Below is a snapshot of what we do to ensure your vaccines are kept at the correct temperature from the moment we receive the stock to when the stock is delivered.



Medi-Log™ II Temperature Data Logger

Automatic setting data logger for vaccine storage

- Meets all the latest Strive for Five guidelines and WHO requirements
- Fully self contained, sets its own parameters automatically for va
- Simple one button operation for all functions
- Free software
- Includes Glycol Vial , 2 Meter External probe and USB Cable
- 16,000 reading memory
- Time interval between readings 1 minute to 99 minutes
- Visual on LCD Screen and Buzzer alarm
- Download data using free software suits Windows computers (USB connection cable - included)

HLPMELOGII

\$149





Quality Guaranteed

We pride ourselves on being a market leader in the importation and distribution of high quality medical equipment at a competitive price. We have put the work in to get the best products for your practice.

All our products meet the strictest Australian Standards.

Pacific Medical, premium equipment,
no substitute for quality.



WWW.PACIFICMEDICAL.COM.AU

Free Standing Exam Couch

- Safe working load 280kg
- Adjustable headrest
- Easy lift and lock mechanism for easy assembly
- Powder coated frame
- Heavy duty PVC vinyl upholstery
- Flat packed for self assembly
- 195L x 71W x 82H cm

Also available with breathing hole

Grey [TMF100MTGR](#)
Black [TMF100MTBL](#)
Navy Blue [TMF100MTNB](#)

\$515



Free Standing Exam Couch with 3 Drawers

Popular fixed height table now available with drawers for extra storage.
Ideal for Doctor's Rooms.

- Loading weight: 280kg
- Adjustable headrest
- Easy lift and lock mechanism for easy assembly
- Powder coated frame
- Flat packed for self assembly
- 3 Drawers: 434mm wide, 499mm deep, 100mm high
- Dimensions: 195L x 71W x 82H cm

Black [TMF100MTBL3D](#)
Grey [TMF100MTGR3D](#)
Navy Blue [TMF100MTNB3D](#)

\$858





Classic 3 Section Couch

- Electric Hi Lo operation by foot switch
- Gas lift head section as standard
- Heavy duty PVC vinyl upholstery that is fire, stain, mildew, oil, water and wear resistant
- Sleek design with rounded frame legs for safety
- Lifts to maximum height in 30 seconds
- Gas strut release for head and foot section from -40° to +65°
- Head section with face hole and plug
- Padding thickness 5cm
- Adjustable feet for uneven surfaces
- Retractable dual locking castors on all 4 corners
- 225kg safe working load
- Warranty 5 years frame, 3 year electrics

Specifications:

- Total Length: 190cm
- Length (head-mid-leg | 62-62-62cm)
- Width: 71cm
- Height: adjustable from 46.5cm to 92.5cm

- TMF100NB33 Navy Blue
- TMF100BL33 Black
- TMF100GR33 Grey
- TMF100CG33 Charcoal Grey



Visit our website to watch the product video

Classic 2 Section Couch

Same features and specifications as Classic 3 Section above except:

- Adjustable head section through dual gas operated struts from -20° to + 75°
- Total Length: 185cm (head-leg | 63-11cm)
- Height: adjustable from 47cm to 92cm



- TMF10002NB Navy Blue
- TMF10002BL Black

\$1595



Visit our website to watch the product video

3 Section Bariatric Couch

- Dual motor. Electric Hi Lo and Electric head section through foot switch.
- Minimum height to maximum height in 30 seconds
- Adjustable head section through motor and foot switch -40° to +65°
- Adjustable leg sections featuring double gas struts, angle adjustable from -40° to +65°
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 75mm castors for safety and mobility
- Retractable dual locking castors on all 4 corners
- Padding thickness 5cm
- Heavy duty PVC vinyl upholstery
- 280kg safe working load
- Warranty 5 years frame, 3 year electrics

Specifications:

- Total Length: 190cm
- Length (head-mid-leg | 62-62-62cm)
- Width: 71cm
- Height: adjustable from 46.5cm to 92.5cm

TMF100NB44

\$2195



5YR
warranty
frame

3YR
warranty
electrics



TOP FEATURES:

- Electronic head rest
- Electronic up and down
- Two motors
- 280 kg capacity

3 Section Couch (short head - 66cm width)

Same features and specifications as
Classic 3 Section on previous page except:

- Total Length: 195cm (head-mid-leg | 40-54-93cm)
- Width: 66 cm
- Height: adjustable from 47.5cm to 93cm
- Most suitable for physio/massage

TMFET3NB

Navy Blue

TMFET3BL

Black

\$1785



5YR
warranty
frame

3YR
warranty
electrics



Paper Roll Holder

- Powder coated
- Fits electric hi-lo tables and free-standing tables



TMFHOLD

\$119



MEET OUR TEAM: BELINDA
CUSTOMER SERVICE - NSW

Belinda joined the team just shy of twelve months ago in the new role of Customer Service Supervisor. Bringing with her a wealth of experience, she sure is enjoying the journey so far. "I love the feeling of being in a big family who are helping to make a difference by providing for the medical industry and donating via B1G1. My favourite thing is knowing that there is constant support from colleagues who genuinely care."

Funnily enough, Belinda once hated coffee and now can't live without it. Her last live concert was Diesel – AMAZING, and she's seen him live at least a dozen times.

What is a pet peeve of yours? Bad grammar

Next travel destination? Banff, Canada

What sport do you love watching? Cricket



5 Section Couch

- Electric Hi Lo operation through foot switch.
- Minimum height to maximum height in 30 seconds
- Adjustable backrest featuring single gas strut, angle adjustable from 0° to + 70°
- Adjustable head section featuring single gas struts, angle adjustable from -35° to +45°
- Adjustable middle section, angle adjustable from 0° to + 25°
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 75mm castors for safety and mobility
- Retractable dual locking castors on all 4 corners
- Padding thickness 5cm
- Heavy duty PVC vinyl upholstery
- 225kg safe working load
- Warranty 5 years frame, 3 year electrics

Specifications:

- Total Length: 195cm
- Length (head-mid-leg | 43-51-96cm)
- Width: 66cm
- Height: adjustable from 50cm to 96cm



TMFET5NB

\$2495



5YR
warranty
frame

3YR
warranty
electrics



Visit our website to
watch the product video



Podiatry / Multi purpose chair

- Electric Hi Lo operation by foot switch
- Electric seat angle adjustment by foot switch, angle adjustable from 0° to + 25°
- Adjustable backrest through dual gas operated struts to +80°
- Minimum height to maximum height in 30 seconds
- Adjustable foot section featuring single gas strut
- Head section has face hole and headrest plug
- Comfortable flip up armrests
- Retractable dual locking castors on all 4 corners
- Padding thickness 5cm
- Heavy duty PVC vinyl upholstery
- 75mm retractable castors for safety and mobility
- 225kg safe working load
- Warranty 5 years frame, 3 year electrics

Specifications:

- Total Length: 190cm without leg extension
- Length (back-seat-leg-foot) | 91-56-38-21cm)
- Width: 60cm
- Height: adjustable from 56.5cm to 98cm



TMF100POD

\$2695



5YR
warranty
frame

3YR
warranty
electrics



Visit our website to
watch the product video



BP Medical (WA) | www.teammed.com.au
PHONE 1300 33 55 60 | info@bpmedical.com.au

All prices are subject to 10% GST unless marked with * (GST Free)



Team Medical | www.teammed.com.au
PHONE 1300 22 44 50 | info@teammed.com.au

[173]

Effortless, Affordable.

ACERO All Electric Examination Couch




The couches are fitted with 100mm fully braked castors that enables transportation where necessary.

Designed and constructed with twin support frames, the Acero Examination Couch delivers excellent stability and weight distribution.

The unique base design gives clear access for patient lifters and over bed tables.

*1-year warranty castors, 2-year warranty electrics



Generous Low Height

Allows for easy ingress by patients and the couch can be effortlessly raised to a comfortable height for the care giver.



All Electric Operation

The all electric operation not only provides ease of use, but is ergonomic and OH&S friendly.



2 And 3 Section Models

Ideal for general examination, assessments and clinical treatment.

2 Section

AX46002SECT

- Electric operated variable height
- Electric operated head section
- Padded, upholstered top, with face cut out
- 100mm castors; fully lockable
- Safe working load: 200kg

Height Range: 490 – 900mm

Head Section: Positive tilt to 85°; 595mm

Size: 1900mm long x 700mm wide

CABAX46002SECT


\$2445



Optional:

Foot Control: AX3021

Side Rails: AXSS70



3 Section

AX46103SECT

- Electric operated variable height
- Electric operated head section
- Electric operated foot section
- Padded, upholstered top, with face cut out
- 100mm castors; fully lockable
- Safe working load: 200kg

Height Range: 490 – 900mm

Head Section: Positive tilt to 85°; 595mm

Foot Section: 720mm

Mid Section: 530mm

Size: 1900mm long x 700mm wide

CABAX46103SECT


\$2825



Optional:

Foot Control: AX3017

Side Rails: AXSS70



+61 (3) 9922 0300
info@axishealth.com

axishealth.com.au



Gynaecology Table with foot pads, stirrups and leg extensions



- Electric seat angle adjustment through foot switch, angle adjustable from 0° to + 20°
- Adjustable back rest section featuring dual gas struts, angle adjustable from 0° to +65°
- Fully adjustable / removable patient leg supports
- Minimum height to maximum height in 30 seconds
- Includes a plug and stainless steel debris tray for examination
- Retractable dual locking castors on all 4 corners
- Padding thickness 5cm
- Heavy duty PVC vinyl upholstery
- 75mm retractable castors for safety and mobility
- 225kg safe working load
- Warranty 5 years frame, 3 year electrics

Specifications:

- Total Length: 185cm
- Length (backrest-seat | 91-46.5cm)
- Width: 660cm
- Height: adjustable from 64.5cm to 100cm

5YR
warranty
frame

3YR
warranty
electrics



Visit our website to
watch the product video

TMF100GYN Navy Blue
TMF100GYNBL Black

\$2995



Premium Gynaecology Table with foot pads

NEW



This Gynaecology Table provides electric positioning for gynaecology procedures, as well as a full length table surface for most general procedures.

- Electric Hi Lo operation through foot switch
- Electric seat angle adjustment through foot switch, angle adjustable from 0° to + 18°
- Adjustable back rest section featuring dual gas struts, angle adjustable from 0° to +71°
- Fully adjustable / removable patient leg supports
- Minimum height to maximum height in 30 seconds
- Includes a stainless steel debris tray for examination
- Retractable dual locking castors on all 4 corners
- Padding thickness 5cm
- Heavy duty PVC vinyl upholstery
- 75mm retractable castors for safety and mobility
- 225kg safe working load
- Warranty 5 years frame, 3 year electrics

Specifications:

- Total Length: 181cm
- Length (headrest-backrest-seat | 61-76-40cm)
- Width: 57cm
- Height: adjustable from 57 - 118cm



TMFGCPRNB Navy Blue

\$4895



Chiropractic Table

This Chiropractic Drop Table is a versatile cost effective table providing flexibility for your practice.

- Electric Hi Lo operation through foot switch.
- Minimum height to maximum height in 30 seconds
- Adjustable leg extension
- 50mm castors for safety and mobility
- 2 fixed wheels and 2 adjustable legs
- Padding thickness 8cm
- Heavy duty PVC vinyl upholstery
- 225kg safe working load
- Warranty 5 years frame, 3 year electrics

Designed with 4 manual drops:

- Cervical
- Thoracic
- Lumbar
- Pelvic



TMFCTNB **\$3895**



Specifications:

- Total Length: 180cm (without leg extension)
- Length (head-mid-leg | 40-43-69.5cm)
- Width: 66cm
- Height: adjustable from 50cm to 96cm

Ergo-Step™ Stool

Light-weight, stackable step stools that provide an alternative to traditional heavy metal step stools.

- Weight capacity of ~220kg
- Anti-skid feet will keep it firmly planted
- Size 13H x 35W X 45Lcm
- Interlocking feature creates a safe work platform
- Distinctive in its bright orange color and easily recognisable
- If stacking, you will need a bottom and top configuration



Excellent stability



Top model **AN1170T**
Bottom model **AN1170B**
(clicks into top)

\$149

Step Up Stools

- Stainless steel frame
- Non-slip ribbed rubber platform
- Reinforced rubber tips
- Manufactured with 1" steel tubing

Single Step

- 45.5 L x 24.5 W x 20 H cm

TMFSSTEP **\$149**

Double Step

- 54 L x 40 W x 34.5 H cm
- First step 165mm; Second step 34.5cm

TMFDSTEP **\$169**



Beauty Day Spa Massage Couch

This electric hi lo table is perfect for use in Beauty, Days Spas and Massage.

- Electric Hi Lo operation through foot switch
- Adjustable angle of back rests 0° to +85°
- Adjustable angle of leg rests -75° to +5°
- Luxurious soft and comfortable head section
- Adjustable feet for uneven surfaces
- Retractable easy move wheel system for safety and mobility

Specifications:

- Padding thickness 7.5cm
- Foam Density 32g/m2
- Table length 186cm
- Loading weight 225kg



TMFBTW
TMFBTNB

White
Navy Blue

\$1695




Contoured Massage Couch

- Electric Hi Lo operation through foot switch.
- Minimum height to maximum height in 30 seconds
- Adjustable head section with manual operation
- Angle adjustable from -75° to +30°
- Arm rest and cradle
- Adjustable feet for uneven surfaces
- 75mm castors for safety and mobility
- Retractable dual locking castors on all 4 corners
- Padding thickness 5cm
- Heavy duty PVC vinyl upholstery
- 225kg safe working load
- Warranty 5 years frame, 3 year electrics

Specifications:

- Total Length: 195cm
- Length (head-body | 32-165cm)
- Width: 70cm
- Height: adjustable from 48cm to 93.5cm



TMFCMNB

\$2035



Portable Massage Couch

- Aluminum frame and legs
- Manufactured without electric welding to join the aluminum
- 220kg capacity
- 2-year warranty on frame
- Open Size: 185 x 66cm x 65-80(adjustable) cm (LxWxH)
- Folded Size: 93 x 67 x 17cm (LxWxH)
- Adjustable Height 65-80cm



TMFPMTNB
TMFPMTBL

Navy Blue
Black

\$379



Portable Massage Chair

- Strongly made steel frame
- Removable chest pad, seat pad, knee support, face rest
- Folds up easily
- 180kg capacity
- 2-year warranty on frame
- Open Size: 54 x 76 x 120-127 (adjustable)cm (LxWxH)
- Folded Size: 54 x 76 x 118cm (LxWxH)



TMFPMC

\$325



Surgeon Stools Standard

- Five star chrome base
- Adjustable height swivel
- Height adjustment 49 to 63cm
- 6cm foam thickness
- Cushion diameter 35cm



Navy Blue TMF100GSNB
Black TMF100GSBL
Grey TMF100GSGR **\$179**

Surgeon Stools Premium with foot ring


- Five star chrome base
- 61cm foot ring
- Adjustable height swivel
- Height adjustment 49 to 63cm
- 6cm foam thickness
- Cushion diameter 35cm

Navy Blue TMF100MA08NB
Black TMF100MA08BL
Grey TMF100MA08GR **\$249**



Saddle stools Standard

- Five star chrome base
- Adjustable height swivel
- Height adjustment 49 to 63cm
- 6cm foam thickness
- Cushion diameter 35cm



Navy Blue TMF100SSNB
Black TMF100SSBL
Grey TMF100SSSG

\$229

Saddle stools Premium with foot ring

- Five star chrome base
- Adjustable height swivel
- Height adjustment 49 to 63cm
- 6cm foam thickness
- Cushion diameter 35cm

Navy Blue TMF100MA07NB
Black TMF100MA07BL
Grey TMF100MA07GR

\$259



Blood Collection Chair

- Extra wide
- Comfortable, soft, supportive high density foam
- Height: 108cm
- Width: 55cm incl arm rests 91.5cm
- Length: 65cm (incl foot rest 84cm)
- Fold away foot rest
- Weight capacity 175kg

\$695

Grey TMFBCC
Black TMFBCCBL



Does your practice meet the national standards to provide a safe, designated children’s area in your waiting room?

Scan the QR code to view our range of Wall panels and Play Centres for kids.



Disposable Curtains

Significant cost savings and improved infection control

Haines® Disposable Curtains are designed to replace your existing reusable curtains in healthcare facilities such as hospitals, day surgeries, medical centres and other healthcare facilities that require the use of curtains. They offer a cost effective alternative to purchasing and washing reusable curtains and with their proven ability to inhibit bacterial growth, they can improve infection control at your facility.

- Manufactured from 100% polypropylene, non-woven fabric – 120gsm
- Fire retardant (AS 2755.2-1985) – Failed to Burn
- Impregnated with anti-bacterial and anti-mildew chemicals including nanometer silver (curtain and mesh)
- Tie Back – sewn on so won't fall off
- Curtains always look fresh, pleated and new.
- Identification Label - each curtain is supplied with an easy to read identification label this is to monitor length of use and provides an easy way to monitor required changes of curtains

Previous studies have demonstrated that pathogens such as Vancomycin resistant enterococci (VRE), methicillin resistant staphylococcus aureus (MRSA) and Clostridium difficile are still present on non-disposable curtains even after one week post laundering (Ohl,M. et al 2011). Even with steam cleaning, the curtains still remained contaminated (White, L.F; Dancer, S.J; Robertson, C; 2007).



Standard Range



2.5m x 2m drop

summer blue	HAIDC25BLU	
cream	HAIDC25C	
latte	HAIDC25L	
grey	HAIDC25GY	Box/15
lilac	HAIDC25LILAC	\$352.65

4.5m x 2m drop

summer blue	HAIDC45B	
cream	HAIDC45C	
latte	HAIDC45L	
grey	HAIDC45GY	Box/8
lilac	HAIDC45LILAC	\$271.25

7.5m x 2m drop

summer blue	HAIDC75B	
cream	HAIDC75C	
latte	HAIDC75L	
grey	HAIDC75GY	Box/5
lilac	HAIDC75LILAC	\$227.85



Print Range



2.5m x 2m drop

spring mood	HAIDC25MOOD	
aussie animals	HAIDC25PAEDAA	Box/15
safari	HAIDC25PAED	\$420.90

4.5m x 2m drop

spring mood	HAIDC45MOOD	
aussie animals	HAIDC45PAEDAA	Box/8
safari	HAIDC45PAED	\$363.60

7.5m x 2m drop

spring mood	HAIDC75MOOD	
aussie animals	HAIDC75PAEDAA	Box/5
safari	HAIDC75PAED	\$368.90

Roller Hooks

For curtains that are difficult to slide along the track, use Haines Roller Hooks for easy sliding.


Box/50 \$55.45 HAIDCROLL50



"S" Type Hooks

Use to connect Disposable Curtains to existing curtain hooks. Ideal for window applications.

Box/50 \$38.00 HAIHOOKS



Curtain Length

Measure lengths of curtains required

Curtain A: 4.5m length for short sections dividing cubicles

Curtain B: 7.5m length for long sections i.e. curved track at front of cubicle

(2.5m length also available for small cubicles, privacy curtains, across windows - fits sizes 0.5m to 1.8m width areas)

A QUICK Q&A WITH OUR ACCOUNTS TEAM



MEGAN

An interesting fact about you?
I'm ambidextrous (using right and left hands equally well)



LINDY

An interesting fact about you?
My first job was at the Zoo



DARINA

What is your favourite takeaway food?
Lebanese kebab



FIONA

What are your hobbies?
Hiking, knitting & day trips on the motorbike



HELEN

What are your hobbies?
Paused triathlon to train my Vizsla dog to be my running buddy



KELLY

Do you have any special talents?
I can say words backwards and remember number plates



JODY

What sport do you love watching?
V8 Supercars



SUSAN

An interesting fact about you?
I'm a beekeeper

FAST AND RELIABLE DELIVERY THROUGHOUT AUSTRALIA

Distributing orders from four warehouses - Sydney, Brisbane, Melbourne and Perth.

Join the team that delivers.



MEET OUR TEAM: ALEC

Alec started out in the warehouse 7 years ago. it was a natural progression that Alec joined our IT Team after finishing his uni degree. *"Team Medical is a great place to work, with a real family culture that values nurturing talent. I feel lucky that I can work in a job I enjoy, for a company that supports my career growth and where I can gain experience working on many different projects."*

He has his hands on all things computer and system related and has been a critical piece of many of our recent digital innovations.

Alec is not to be messed with at table tennis or Mario Kart during lunch breaks! Alec follows Parramatta Eels in the NRL and loves to play tennis.

What is an app you use almost daily? Wordscapes

What series have you been watching recently? Miracle Workers

What is one of your pet peeves? When people don't know their passwords

IT - NSW



Hospital Curtain Track and Curtains

- Ready to go shapes to suit almost every situation
- Easy to install. Use your own contractor or DIY
- Track can be cut to size on site for perfect fit
- Includes fixing brackets & hooks
- Coverage 4 metres long x 1.95 metres drop



PLEASE NOTE:

J shape is not a reversed L shape!
The rails are on one side only. So choose carefully.
Not sure? Send us a copy of your plans and we'll help you out.

≤ 2500

≤ 1300

L-SHAPE

"L" Shape
2500 x 1300
CRFLSHAPE
\$349

≤ 2500

≤ 1300

J-SHAPE

"J" Shape
2500 x 1300
CRFJSHAPE
\$349

STRAIGHT TRACK
≤ 3800

Straight track
3800 long
CRF3800STR
\$280.50

Need a different size?
Listed above are the most common sizes.
Different sizes and shapes are available.
Call customer service for more info.



Curtains

4m x 1.95m drop

This fabric is reversible and shown railroaded. 100% Polyester FR.

\$340



Pearl
CRFPEARL

Dune
CRFDUNE

Monaco Green
CRFGREEN

Storm Blue
CRFBLUE

Folding Mobile Screen



- Each panel size: 90cmW x 166cmH
- Panels can fold fully out into a straight line
- Powder coated steel frame
- Privacy screen starts low to the ground
- Curtains charcoal grey



DAL1360

\$650



Ferno Mobile Privacy Screen

A full height, lightweight, stable and secure screen

Each panel is made of lightweight vinyl that is heat sealed and easy to clean. The rigid frame is constructed from anodised aluminium for easy manoeuvrability. The legs and castors are removable so the unit can be packed flat for storage or freight.

Folded: 66cm
Extended: 192cm
Height: 158cm



FERCRP3M

\$968.30




Quality Single Stainless Steel Instrument trolleys


- High grade 304 stainless steel throughout
- Soft closing magnetic drawers
- German made ball bearing drawer runners
- 125mm deep drawers
- 95mm smooth rolling castors for easy movability
- Locks on front castors
- Top and bottom side rails on all 3 sides
- Double safety folds
- Smooth welding



★★★★★ PREMIUM PRODUCT



No Drawer
50L x 50W x 90H cm
TMF200001
\$446



1 Drawer
50L x 50W x 90H cm
TMF200002
\$583
1 Drawer with Lock
50L x 50W x 90H cm
TMFSST01LK
\$636



2 Drawer
50L x 50W x 90H cm
TMF200003
\$776
2 Drawer with Lock on Top Drawer
50L x 50W x 90H cm
TMFSST02LK
\$803



3 Drawer
50L x 50W x 90H cm
TMF200005
\$987
3 Drawer with Lock on Top Drawer
50L x 50W x 90H cm
TMFSST03LK
\$1016



4 Drawer
50L x 50W x 90H cm
TMF200004
\$1172
4 Drawer with Lock on Top Drawer
50L x 50W x 90H cm
TMFSST04LK
\$1225



Quality Double Stainless Steel Instrument trolleys

- High grade 304 stainless steel throughout
- Soft closing magnetic drawers
- German made ball bearing drawer runners
- 125mm deep drawers
- 95mm smooth rolling castors for easy movability
- Locks on front castors
- Top and bottom side rails on all 3 sides
- Double safety folds
- Smooth welding



2 Drawer side by side

80L x 50W x 90H cm

TMF200007

\$890




4 Drawer - 2 over 2

80L x 50W x 90H cm

TMF200008

\$1131



No Drawer

80L x 50W x 90H cm

TMF200006

\$501



1 Drawer (Full width)

80L x 50W x 90H cm

TMF200009

\$836



2 Drawer (Full width)

80L x 50W x 90H cm

TMF200010

\$995



Powder Coated Trolleys


- Soft magnetic closing drawers
- Smooth running castors
- Lockable front wheels
- 125mm deep shelves
- Top and bottom rails on 3 sides
- Outstanding finish



No Drawer

50L x 50W x 90H cm
TMF200200


\$433



1 Drawer

50L x 50W x 90H cm
TMF200201

\$485



2 Drawer

50L x 50W x 90H cm
TMF200202

\$561



3 Drawer

50L x 50W x 90H cm
TMF200203

\$662



4 Drawer

50L x 50W x 90H cm
TMF200204

\$729



Anaesthetic Trolleys

- High quality stainless steel
- 125mm heavy duty castors



4 Drawer

60L x 49W x 97H cm
• 2 Drawers - 125mm deep
• 2 Drawers - 250mm deep

TMFSSAT04

\$1425



6 Drawer

90L x 49W x 97H cm
• 1st drawer 75mm
• 2nd & 3rd drawer 150mm
• 4th drawer 75mm
• 5th & 6th drawer 100mm,
• 2 shelves with rails

TMFSSAT06

\$1952



Resuscitation Trolleys

- High quality stainless steel
- 125mm heavy duty castors



5 Drawers

75L x 49W x 97H cm
• 4 half drawers - 125mm deep
• 1 Drawer full width - 125mm deep

TMFSSRT05

\$2095



7 Drawers

75L x 49W x 97H cm
• 6 half drawers - 125mm deep
• 1 Drawer full width - 250mm deep

TMFSSRT07

\$1643



Emergency Trolley

This rugged, portable crash cart is constructed with an ABS modular surface, making it easy to clean and resistant to corrosion. It also features stainless steel pillars, aluminum composite back/side boards and aluminum alloy detachable drawers.

Drawer loading capacity is 50kg, overall loading capacity is 200kg

- Central or seals lock for medicine security
- Trolley has central locking system through side key
- 100mm noise resistant castors with brakes
- Luxurious push handle
- Soft closing drawers
- Dimensions: 650 x 480 x 980mm (LxWxH)
- Drawer sizes: 3 x 110mm; 1 x 135mm; 1 x 195mm

All accessories included:

- Defibrillator shelf
- Extension table
- 304 Stainless steel IV pole
- Cardiac CPR board
- Waste bin
- External power holder & plug board
- Oxygen tank holder
- 304 Stainless steel wire basket
- ABS divider
- Non slip mat



TMFETRED

\$1876



Treatment Trolley

Multi-function treatment trolley with a range of accessories.

- Soft Closing Drawers
- Castors with brakes
- Stainless steel pillars
- ABS modular surface
- Luxurious push handle
- Central or seals lock
- Drawer size 110mm
- Top Loading weight of 75kg
- Dimensions: 650 x 480 x 900mm (LxWxH)

All accessories included:

- Non slip mat
- Extension Table
- Waste Bin
- 304 Stainless steel wire basket



1 drawer TMFTTR01

\$928

2 drawer TMFTTR02

\$1051

3 drawer TMFTTR03

\$1121



BP Medical (WA) | www.teammed.com.au
PHONE 1300 33 55 60 | info@bpmedical.com.au

All prices are subject to 10% GST unless marked with * (GST Free)



Team Medical | www.teammed.com.au
PHONE 1300 22 44 50 | info@teammed.com.au

Computer Trolley

with printer shelf

- Display arm
- Disinfectant Holder
- Lock & Key for laptop box (separate keys for each side)
- Laptop sits inside laptop box underneath worktop and connects to external screen
- Gas Lift Handle for height adjustment
- Slide keyboard holder
- Pull-out mouse tray
- Mouse holder
- Printer tray

Dimensions: 54D x 70W cm
Height to top: 119 -142H cm (A)
Height to worktop: 81-116H cm (B)
Base (printer tray): 51 x 68 cm
Castors: 4 x 100mm twin wheel castors with brakes
Working Top: 53L x 48W cm
Height adjustability: Manual Gas Spring adjustment



Top

Worktop

(B) 81 -116H
(A) 119 -142H

PACMWST

\$1995



Medication Trolley

Medicine trolley with four drawers, ideal for storage of practice medications.

- Soft Closing Drawers
- Castors with brakes
- Stainless steel pillars
- ABS modular surface
- Luxurious push handle
- Central or seals lock
- 3 drawers x 110mm
- 1 drawer x 140mm
- Top Loading weight of 75Kg
- Dimensions: 650 x 480 x 900mm (LxWxH)

All accessories included:

- Non slip mat
- Extension Table
- Waste Bin
- 304 Stainless steel wire basket



TMFMEDT

\$1372



Mayo Instrument Table

- Welded stainless steel construction
- Telescopic height range: 90 - 135cm
- Removable tray: 56 x 41 x 2cm deep
- 75mm diameter twin wheel castors



VOSV70 \$561



Mayo Instrument Table

- Stainless steel construction
- Removable tray: 60 x 40 x 2cm deep
- Telescopic height adjustment from 85cm to 140cm
- 65mm castors



TMFMIT \$449



IV Stand Stainless Steel

- Height range 1350 - 2050mm
- Stainless steel construction
- 5 Leg weighted base
- Smooth rolling and lockable castors



2 Prong

TMFIVP2H

\$251.00



4 Prong

TMFIVP4H

\$294.00



IV Stand with Plastic Base

- Moulded plastic base with stainless steel pole
- Height: 1350mm to 2200mm
- Base diameter: 600mmW
- Five 50mm twin wheel castors

2 Prong

DAL1564

\$324.50

4 Prong

DAL15644P

\$259.00



Mobile Base

with stainless steel base



RBKML131
\$179

- Comes with both 13 & 16mm adapter
- 5 castors (2 brake castors)
- Flat packed for self assembly



Linen Skip

- Stainless steel



Single no lid

38 x 54 x 85cm

TMFSLT

\$315

Double no lid

78 x 54 x 85cm

TMFDLT

\$469

Single with foot operated lid

38 x 54 x 85cm

TMFSLTL

\$418

Double with foot operated lid

78 x 54 x 85cm

TMFDLTL

\$624



MEET OUR TEAM: JAKE
MANAGEMENT - QLD

As our National Operations Manager, Jake leads our four warehouses with his unique blend of enthusiasm, passion and positivity. A big believer in getting things done right the first time in an efficient manner, Jake is always on the hunt for new ways to improve our fast, smooth and safe operations.

Always up for a challenge, Jake has recently project managed the seamless transition to our new Team Medical home in Sydney. Ensuring we have the right infrastructure and continued scalability to grow alongside our growing customer base.

Jake attributes our success to our company culture being built on strong family values and teamwork. *“Every Team Member is valued. Our Team is invested in their jobs and are willing to go the extra mile.”*

Always on the go, you’ll often find Jake fishing, bushwalking, playing soccer, watching a live gig or cheering on his beloved Sydney Roosters in his spare time.

What is your favourite takeaway food? Indian – Chicken Vegetable Korma.

What is your next holiday destination? No trips planned yet, I’ve recently returned from 5 weeks in Central America.

What is something people may not know about you?
When I was a kid, I had pet blue tongue lizards.



www.teammed.com.au

Everything you
need at your
fingertips



With a teammed.com.au login you can access over 7,000 products, place orders, create order templates, pay invoices, track your orders and so much more.

For feel-good ordering register online today.



WE ARE A BUSINESS
FOR GOOD.



In August, Team Medical rolled over our first-year anniversary of giving with B1G1. For the year we funded the planting of over **4,338 trees across the globe**. This included **756 tangerine trees** for families in Malawi to generate an income, and **2793 trees in India** to increase forest cover. While supporting these communities, we also offset the kilometres driven by our sales reps, as they travel to meet our customers.

Learn more about our giving impacts at www.teammed.com.au/teamgiving



Storage Trolleys

Wire chromed modular system

- Ready-to-assemble MOSYS system
- Chrome plated steel
- Multiple configurations are possible
- Ideal for storage, transport
- 4 x 100mm castors (two with brakes)

Model P181940415

- 1750mm H
- 4 shelves
- 460mm W x 1500mm L
- Includes castors

MDQP1940415 **\$450**

Model P181940615

- 1750mm H
- 4 shelves
- 610mm W x 1500mm L
- Includes castors

MDQP1940615 **\$575**

Additional accessories:

A wide range of accessories and components are available.

- Shelf dividers
- Edge fences
- Push handle
- Label holders
- Wire baskets and more



francehospital®
trolleys for healthcare



WALLY

Hanging waste collector with knee opening system

The bin comes with an inner removable bin with flaps for fixing the waste bags

Wall fixing screws and plugs are supplied as standard together with mounting instructions .

Wally 8 Ltr

22 x 22 x 36 H cm

Grey/Grey **MLQWAL1**

\$68

Maxi Wally 20 Ltr

35 x 25 x 52 H cm

Grey/Grey **MLQMWALGY**

Grey/Blue **MLQMWALBLG**

\$126

francehopital®
trolleys for healthcare



8 Litre Model



20 Litre Model



Waste bag trolley



- Stainless steel construction
- Holds 55cm wide plastic bags
- Solid platform to support the bag
- Foot pedal operated ABS plastic lid
- 5cm diameter twin wheel castors
- Overall dimensions 41 x 37 x 84cm high



VOSV286

\$407

Waste bags for Vost trolley

60L Yellow

CNW006

Box/200 **\$146.00**

60L Blue

GAR006

Box/200 **\$139.95**

Pedal Bin - Matte Steel Fingerprint Proof

Motion Control Whisper Lid

- Solid metal lid
- "Motion Control" whisper lid, designed to ensure smooth and silent operation
- Fingerprint proof
- Clever hinge design - lid stays open if desired
- Strong plastic inner bucket - removable and easy to clean
- Robust pedal mechanism and high quality materials - corrosion resistant
- Non-skid base - pedal bin stands stable, even on wet or polished floors
- Plastic floor ring - less chance of damage to the floor
- 10 year guarantee



brabantia®
solid company



30 Litre

66H x 39.5W x 29.5Dcm

SHE06112

\$160.75



BP Medical (WA) | www.teammed.com.au
PHONE 1300 33 55 60 | info@bpmedical.com.au

All prices are subject to 10% GST unless marked with * (GST Free)



Team Medical | www.teammed.com.au
PHONE 1300 22 44 50 | info@teammed.com.au

X-Ray Viewers Standard

- High luminance
- Lightweight
- Easy change fluorescent globes
- Easy grip film attachment on all sides of the viewer



Single Bay	44L x 12W x 59H cm	TMFSBXST	\$164
Double Bay	80L x 12W x 59H cm	TMFDBXST	\$278
Triple Bay	116L x 12W x 59H cm	TMFTBXST	\$418

X-Ray Viewers Slimline

- LCD display
- 2.4cm thickness
- 100,000 hours working time
- No globes required
- Lightweight
- High luminance
- Sleek design



Single Bay	47L x 2.4W x 55H cm	TMFSBXSL	\$356
Double Bay	95L x 2.4W x 55H cm	TMFDBXSL	\$562
Triple Bay	123L x 4.1W x 55H cm	TMFTBXSL	\$922

Drug Safe

- Constructed from 10mm steel plate to Health Department specifications
- Finished in durable grey enamel paint
- Six lever pick resistant key lock
- Four 13mm mounting holes are provided in the back plate of the safe
- Removable shelves and concealed hinges



- 3 Shelves
- 320mmW x 150mmD x 450mmH
 - Weight: 46kg

DAL1620 \$1045

- 4 Shelves
- 450mmW x 175mmD x 600mmH
 - Weight: 70kg

DAL1622 \$1171.50



Dexion Maxi Bins

The inherent strength of the Maxi bin comes from the robust design using the best quality polypropylene materials. Maxi bins are most commonly hung on louvred panels, but their unique interlocking design enables them to be stacked securely without the bin splaying under load. This makes them perfect for incorporating within your shelving system or on a bench top.

- Bins can be stacked or hung on louvred panels
- Smooth internal features with curved edges makes bins easy to clean

Maxi Bins (W x D x H)

P05	100mm x 60mm x 120mm	COP80102	\$4.20
P10	100mm x 85mm x 178mm	COP80122	\$7.05
P20	125mm x 135mm x 225mm	COP20BLUE	\$9.35
P30	165mm x 205mm x 280mm	COP30BLUE	\$17.90
P40	165mm x 410mm x 280mm	COP40BLUE	\$28.15


Dividers

P30 Maxi Bin Divider	(Bin P30 can hold 1 divider)	COP30DIV	\$6.40
P40 Maxi Bin Divider	(Bin P40 can hold 3 dividers)	COP40DIV	\$7.25

Wall mounted Louvre Panels

Single sided panels are designed to be wall mounted. Perfect for small parts picking. The bin sits securely on the panel, but is easily removed when required

Louvre panel 45W x 47Hcm	COPLP18	\$34.85
Louvre panel 45W x 91Hcm	COPLP36	\$67.20



Clocks



Model CS1067WH
225mm diameter analogue clock

\$32.40
DAVJCS1067



Model 6200
300mm diameter analogue clock

\$39.45
DAVJC6200

Stainless Steel Oxygen Trolley
D Size

- 304 stainless steel construction
- Safety chain for securing oxygen cylinder

TMFOTD
\$315.00



Oxygen Bottle Trolley C size

- Powdercoated
- Height of handle: 1m
- Bottle holder: 140mm dia. x 300mmH


AXH537
\$325.00



C-Size Oxygen Cylinder Holder

- Hooks over trolley rail
- Full stainless steel construction
- 125mmW x 110mmD x 350mmH

DAL1480
\$165





Team Medical is the leading distributor of Melag Autoclaves

Why Melag?

Read the philosophy of a family business

MELAG is an owner operated family business that has consistently focused its efforts on products for practice hygiene since the company was first founded in 1951. This focus has permitted us, as a medium-sized company, to continue to develop our products using a highly specialised development team so that they are market-ready and of a leading international standard.

We have constantly strived to enhance quality while delighting in innovation and combined this with the highest safety objectives, all of which have resulted in our position as world leader in the area of hygiene, having supplied over 450,000 devices to customers. Production is carried out exclusively by our 300 employees in Berlin in a premises of approximately 20,000 m².

We have never seriously considered relocating production to low-wage countries, although this practice is widespread elsewhere. Other firms have often reported that while this relocation leads to savings, quality is often compromised. We have also resisted the trend to expand our range of products that are not part of our core competency: surgical hygiene. We are convinced that consistent specialisation is a decisive criterion for the quality of our products. Our current mantra now and in the future shall remain:

"competence in hygiene" and most importantly: Quality - made in Germany

Things to consider when buying an Autoclave

1. What type of items are you sterilising?

Hollow, porous, cannulated items – you'll need a Class B Autoclave

Solid instruments – you'll need a Class S Autoclave

2. How many sterilising cycles do you perform per week?

This will be a determining factor in your choice of autoclave (i.e. cycle time, size, load of instruments).

3. Where are you planning to put your autoclave?

Your available bench space and ventilation are important factors


4. What does the calibration/installation and validation process entail?

The installation process confirms the autoclave is positioned with sufficient bench space and ventilation. Dedicated power outlets are checked and a sterilisation cycle is run to ensure no interference from other equipment in the practice.

Validation determines your steriliser meets the criteria and performance capabilities.

Validation should be completed annually or after any commissioning by technicians.

Installation and Validation of Instrument Processing equipment is a requirement of Australian and New Zealand standards AS/NZS 4815 & 4187.






MELAG

competence in hygiene



Quality – made in Germany

MELAG
competence in hygiene



MELAtherm®10

The washer-disinfector

With active drying and short operating times

The MELAtherm 10 is a washer disinfector complying with the EN 15883 European Standard. It is ideally used in processing instruments by an automated process.

With compact external dimensions, the MELAtherm 10 offers full features without any additional charge:

- Feeding process agents with the dosing unit
- Batch documentation with the CF card or network
- Active drying in accordance with the Standard
- Available in single phase or 3 phase models

MELAtherm 10 automatically monitors the rotation of the rinse arms, the rinse pressure and the filter sieves. Errors in the program mode or during operation are thus prevented. MELAtherm 10 can be fitted in floor units without any problem, the standard width being only 60 cm.

Should the MELAtherm 10 be free-standing on top of a floor unit, MELAG offers a stainless steel cover plate. The optional floor unit, also available, facilitates ergonomic working height and creates additional storage space. MELAtherm 10 is equipped with a basic basket, with or without injector rail, at your option. Up to 11 hollow body instruments can be processed with the injector rail, e.g. suction tubes or dental handpieces and contra angles etc.

The water-softening unit integrated in the MELAtherm 10 provides consistent results in the treatment of water, irrespective of the quality of the tap water. Those practices seeking to minimise stain formation on their instruments can connect the MELAtherm 10 to an external water treatment unit, e.g. MELAdem® 53/53 C, which supplies the device with de-ionized water (DI water).

\$11,895

BIO10021



LARGE ITEM FREIGHT
CHARGES WILL APPLY



For more information
download the product
brochure from our website



 Quality – made in Germany

Now released the MELAtherm® Evolution

The washer-disinfector for those busy practices
that need to process many instruments in one cycle.

Contact Us for more information.

