

# AMEDTEC ECG Pro®

## PC-Based ECG Machine


ECG with high flexibility, Simple handling, Excellent measurement technology & Reliable pacemaker detection. Perfect for paperless workflow.

**The 12-lead resting ECG with optimal workflow. Fast, precise and simple.**

- One-key operation - Simply press the start button to start the ECG recording
- Switch between pre-configured ECG views while recording
- Many configuration options - compile individual diagnostics programs for different patient types, i.e. geriatrics and paediatrics
- Freeze and scroll back option during examination while live ECG continues on bottom of screen
- Predefined text blocks help you compile the final report
- Wide selection of print formats with ability to easily add other formats
- Banana Plug leads for each electrode connection
- Electrodes check; faulty application is easily identified by graphic indicator
- German quality
- 3-year warranty

### Contents:

- 1 x ECG Cables
- 1 x Pkt/10 Alligator Clips
- 1 x ECG Pro Software and Interpretation
- Test Certificate
- Laminated cheat sheets



## USB Model

CRI100434 **\$2695\***

## Bluetooth Model

CRI100437 **\$3235\***




**Switch between pre-configured ECG views while recording**  
**Freeze and scroll back during examination while live ECG continues on bottom of screen**



★★★★★ PREMIUM PRODUCT

**3YR**  
warranty



## Meet Our Team: Thiago

Having grown up in Brazil, you could be forgiven for assuming Team Medical's latest addition to the IT team simply lived and breathed for all things soccer. But Thiago had other plans, pursuing a path in computers instead. He worked first in IT infrastructure and then sharpened his skills, enabling him to move into programming and development.

Outside of computers, Thiago loves to lead an active life. He is something of a gym junkie and enjoys going on long hikes and drives to explore new places. Although he insists he is too old now, he also used to spend much of his time skating, surfing and even whipping down snow-crested mountains on his snowboard. In his youth he was also quite a versatile musician and even spent some time learning the craft as a DJ producer. Who has time for soccer, right?





**An advanced paperless ECG solution**  
**Clean trace & comprehensive reporting**  
**Medical Director compatible**

### PC-Based 12 Lead ECG

- Easy interface to Medical Director and other EMR (Electronic Medical Records) and HIS (Hospital Information Systems)
- Measurement and interpretation included
- Comprehensive display and printout for caliper QRS and ECG measurements
- Advanced digital filters for real-time and stored ECGs to create exemplary ECG traces
- Plain paper print-out to local or network printers
- Enabled for post processing of saved data
- Export ECG test results into formats: PDF, HL7, JPG, XML, SCP, DICOM
- 12 months warranty



### NORAV 1200M+I3

BIO1200MI3

\$2585\*

Includes I3 Software; It has advanced algorithms for diagnosis. Best suited for GP's who specialise in Cardiology or Cardiologists.

### NORAV 1200M Bluetooth

BIO1200BT

\$2695\*

Wireless Bluetooth connectivity



### Neo ECG S120 ECG Tablet

Ideal to use between rooms or for home visits

- 7" ultra-high resolution colour touch screen
- Size: 19.7cm(L) x 11,25cm(W) x 2.61cm(H)
- Easy to operate
- Portable design, compact in size
- Synchronous acquisition and display of 9/12-lead waveform, as well as heart rate detection
- Automatic pacing detection and marking, RR analysis, HRV, medicine test, ECG event mode
- Freeze the ECG waveform on the screen
- Output files in multiple formats including PDF, BMP, HL7, DICOM, SCP
- Store, preview, review, edit, export, upload, print and search patient data
- Wireless transmission of ECG data via Wi-Fi and mobile networks
- Built-in Wi-Fi and USB Ports for transmission of ECG data and printing
- Supports laser printer
- Glasgow Analysis Algorithm

**SCAN ME**  
 Visit our website to watch the product video and download the fact sheets

MDQ31300273

\$1890\*



### Neo ECG T180


Also available as a 10" tablet with base and build-in thermal printer.  
 Tablet Size: 27.26cm(L) x 18.12cm(W) x 2.58cm(H)

MDQ31300302  
 \$3015\*

## Diagnostic Cardiology Suite 12-Lead PC-Based ECG

Designed with your workflow in mind, the Diagnostic Cardiology Suite ECG provides results you can trust in one easy-to understand format.

- Intuitive design for simplified workflow
- Provides high-fidelity data to support better decision-making
- "Best 10" technology can capture ECG data with the least amount of noise to help reduce the need for repeat tests
- One device for spirometry and PC-based ECG software allows for a common workflow and setup
- Flexibility to use as a standalone ECG or spirometer, or both devices together
- Complies with ACC, AHA and HRS recommendations for adult and paediatric ECG acquisition



Laptop not included

### PC-Based ECG

#### Wireless Acquisition Model (WAM)

(most popular)

**\$3070\*** [WACCRXXWAXX](#)

### PC-Based ECG

[WACCRXXAAXX](#)

**\$2910\***

### Optional Spirometry Suite

[WADCSSCX](#)

**\$3015\***

Includes Calibration syringe

This suite can be purchased on its own; ECG module is not required



## CP 150™ Resting ECG

Experience simple, streamlined ECG testing with a solution that's designed to help your practice run more efficiently. Designed for your primary care practice, the CP 150 ECG helps streamline your workflow with clinical decision support and full-size print outs. The 7" colour touchscreen display lets you easily enter patient data.

- Intuitive touchscreen display - can improve accuracy, minimize training and save time
- Lead quality display - help avoid repeat ECGs by identifying lead connection problems before completing a test
- Instant-on feature - help eliminate typical start-up delay
- 12-Lead ECG preview large display provides three, six or 12-lead ECG preview options
- ECG Print Preview can save time and reduce paper costs by verifying the ECG test quality prior to printout
- Optional spirometry module lets you perform two tests with one device

### CP 150™ Spirometry Option


- Perform FVC and pre- and post-BD testing with automatic interpretation and reversibility analysis
- Customise report formats and screen displays to fit your needs and workflow
- Integrated ambient sensors can eliminate the need for the user to find and enter values for temperature, pressure and humidity

### CP 150 ECG

[WACP150A66](#)

RRP \$4446

**\$4690\***



### CP 150 ECG with Spirometry

[WACP150AS1EN66](#)

RRP \$6384

**\$6795\***




## EDAN SE-601C ECG

with PDF reporting

This compact 6-channel 12 lead ECG is perfect for both regular tests and emergencies. The portable, lightweight design makes it convenient for mobile use and offices with limited space.

- Large 5.7 inch colour screen
- Heart rate variability analysis
- Pacemaker detection which meets the requirements of ANSI/AAMI EC11
- Internal memory for up to 100 ECGs
- Open communication with data expert in PDF through USB drive or LAN (WiFi optional)
- Alphanumeric keyboard
- Built-in rechargeable Li-ion battery AC/DC power supply
- Automatic interpretation for adults and paediatrics
- Internal thermal printer (110mm)
- 3-year warranty



\$2475\*

BIOEDSE601



Trolley Options:

Rolling Stand MT201 for SE601

Rolling Stand

EDAMT201 \$700

## EDAN SE-1200 EXPRESS BASIC ECG

with PDF reporting

Combines compact, lightweight design and brilliant performance with a built-in high-resolution thermal printer

- Alphanumeric keyboard and one touch operation
- Standard 12 lead ECG
- 8 inch LCD screen
- High resolution blue monochrome LCD screen
- Records 12 ECG traces on A4 size paper
- Automatic measurement and interpretation tested with authoritative CSE database
- Built-in rechargeable Lithium battery
- Internal memory for up to 800 ECGs
- Built-in thermal printer A4
- Open Communication with data export in PDF (USB drive)
- 3-year warranty



\$2910\*

BIOEDSE1200EB



A4 format | 200 ECG storage | Full interpretation |  
High resolution colour touch screen | PDF report | 3-year warranty



## EDAN SE-1200 EXPRESS ECG

**With touch screen and one touch operation, SE-1200 Express brings you a great user experience with easy and fast operations.**

- Easy to use high-resolution 8" colour screen with touch function
- Alphanumeric keyboard and one touch operation
- Advanced Algorithms used for interpretation and measurement (included)
- PDF report in easy steps
- Waveform freeze, review and edit
- 120 seconds ECG waveform reviewable
- Diagnose results preview before printing
- Internal storage 800 ECGs
- Heart rate variability analysis
- Two-step exercise test with periodic recording
- Heart rate variability (HRV) analysis
- Built-in rechargeable Li-ion battery
- Standard 12 lead recording with A4 size printing
- 3-year warranty



**3YR**  
warranty

**\$3665\***

BIOEDSE1200EX

**Ask for your demo and be impressed!**

## Meet Our Team: Tom

Tom's journey at Team Medical started in 2013 when he decided to join the family business alongside his brother - Matt, Aunt - Marijke and Dad - Pieter. Today, Tom leads our management team as our General Manager.

He also looks after digital and customer experience strategies. He has a wealth of experience in how digital can help make not just our customers' jobs easy but how it can help make the day-to-day operations of our business run more effectively.

Tom values the company culture which supports him in pioneering change within the industry: *"We're innovators, always looking to make our customers' jobs easier. I love my role because I feel I can make a real change for our people and our customers."*

A true sculptor of our team culture with his ongoing positivity and enthusiasm – he literally runs up the stairs every morning and is always keen for some fun at lunch in the recreation room.

He is a huge soccer fan and proud supporter of Arsenal. On winter weekends, he plays soccer with Matt and coaches his sister's team.

**What is your favourite takeaway food?**

Mexican - Tacos

**What TV show have you been streaming recently?**

True Detective



MANAGEMENT - NSW



### Which Type?

Ambu differentiates its range of electrodes based on the type of backing material (breathable or occlusive), the type of gel (wet gel or solid gel), style of connector (stud, wired, slot, tab, radio-translucent etc.), adhesive type (high-tack, normal tack, sensitive) and size.

### Short Term or Long Term


When choosing an electrode for purchase, we recommend taking these features into account. For example, the clinician should consider the wear-time of the electrode and determine whether they would be better off selecting a short-term or long-term style based on the clinical area of use, e.g. simple monitoring (short-term) or holter monitoring (long-term).

### Wet Gel or Solid Gel

Similarly, Ambu recommends considering the respective benefits of wet gel (best suited for diagnostic use) and solid gel (better suited for general monitoring use) to determine the best fit for an organisation: a diagnostic style wet gel can be used for monitoring but a solid gel monitoring electrode may not give an optimal diagnostic result.

### Best Trace, First Trace

Correct electrode choice enables clinicians to achieve 'best trace, first trace' with less time spent on trouble-shooting artifact issues, fewer repeat tests and less waste of clinical time. This provides cost benefit efficiencies that may not necessarily be realised simply by selecting the cheapest electrode available.



## WhiteSensor

The flexible, reliable and price optimised range of Ambu® electrodes

## WhiteSensor

### Functionality at lower costs

Responding to the growing demands of standard solutions to cover a broader range of applications, WhiteSensor has been developed as a range of electrodes that unites price optimisation with the well-known Ambu® quality and functionality.

#### Ambu® White Sensor WS

- Highly conductive solid gel
- Strong adhesion
- Flexible foam backing
- Standard centre stud connection



BIOAMWS00S

Pkt/50 \$14.85

12+ \$12.90

#### White Sensor 0415M

- Solid gel tab electrode
- Latex free
- PVC free (in contact with skin)
- Fish tail tab design reducing electrode lift-offs
- 34L x 23W mm



BIO0415M

Pkt/100 \$14.85

10+ \$13.20

#### White Sensor 4540

- Highly conductive wet gel
- Comfortable foam backing
- 4mm fitting connector
- 32L x 34W mm



CON4540M

Pkt/60 \$29.00

5+ \$24.25

### Looking for the best deal on quality equipment? It's got to be Team Medical

Visit our online sale page for top prices on all your equipment essentials.

Head to: [www.teammed.com.au/onlinesale](http://www.teammed.com.au/onlinesale).



**YOU GET THE BEST DEAL WITH TEAM MEDICAL**



# BlueSensor

The world's most time and cost efficient electrode

## Exact monitoring - clear picture

The BlueSensor electrodes give a fast and reliable trace without the need for skin preparation, maximising comfort for the patient. Fast, exact monitoring that ensures a clear picture of the patients' heartbeat is essential to accurate diagnoses. BlueSensor electrodes will give you exactly that – an absolute quality experience.

### Ambu Blue Sensor M



M-00-A; 4mm fitting  
M-00-F; 3mm fitting

M-00-S; Stud

The small Ambu Blue Sensor M features highly conductive wet gel and superior adhesion to ensure optimal signal quality during short to medium term applications. Thanks to the special occlusive backing material, the electrode will not absorb liquids.

- Highly conductive wet gel
- Superior adhesion
- Occlusive backing material
- Ideal for resting ECG

size 40.8L x 34W mm

M-00-A 4mm	BIOAMMOOA	Pkt/50	\$26.90	20+	\$24.00
M-00-F 3mm	BIOAMMOOF	Ctn/1000	\$432.25		
M-00-S stud	BIOAMMOOS	Pkt/50	\$25.25	20+	\$20.45

### Ambu Blue Sensor L

- Highly conductive wet gel
- Combination of instant and long-term adhesives for excellent contact
- Breathable microporous material
- Comfortable to wear, ideal for long-term monitoring



L-00-S stud  
BIOAMLOOS

Pkt/25 \$16.70  
20+ \$14.50

size 68.2L x 55W mm

### Ambu Blue Sensor SP

- Highly conductive wet gel
- Superior adhesion
- Comfortable foam backing; gentle to skin
- Unique offset concept
- Ideal for short-term stress testing



SP-00-S stud  
BIOAMSP0OS

Pkt/50 \$27.60  
20+ \$23.80

size 42.5L x 38W mm

### Ambu Blue Sensor T

- Highly conductive wet gel
- Superior adhesion
- Oblong shape permits the placement of six closely spaced electrodes
- Comfortable foam backing; gentle to skin
- Large measuring area




T-00-S stud  
BIOAMT0OS

Pkt/25 \$16.70  
20+ \$14.45

size 57L x 35W mm

### Ambu Blue Sensor SU



SU-00-A; 4mm fitting  
SU-00-F; 3mm fitting

SU-00-C; Clip  
size 49L x 33W mm

### The Ambu Blue Sensor SU is the world's first tab electrode featuring a wet gel.

The combination of wet gel and superior quick adhesion provides optimal signal quality during 12-lead diagnostic resting applications. Ten SU electrodes are packed on a single liner so there is no need to peel off individual backing material. Thanks to the special occlusive backing material, the electrode will not absorb liquids.


- Highly conductive wet gel
- Superior quick adhesion
- Ten electrodes on a single liner



SU-00-A 4mm	BIOAMSU0OA	Pkt/60	\$22.75	10+	\$19.60
SU-00-F 3mm	BIOAMSU0OF	Pkt/60	\$18.15	10+	\$16.15
SU-00-C tab	BIOAMSU0OC	Pkt/100	\$18.80	10+	\$16.20

## How to place your Electrodes?

Run a finger around the outside of each electrode without touching its centre thus preventing the gel from spreading onto the adhesive area and ensuring that the electrode stays in place




## 3M Red Dot™ Electrodes

For over 40 years, clinicians have trusted 3M™ Red Dot™ Electrodes for reliable ECG monitoring solutions, tailored to their specific requirements.

### Red Dot™ Foam Monitoring Electrode with Sticky Gel

- Multipurpose foam monitoring electrode that is impermeable to fluids.
- Utilises an advanced, high performance conductive adhesive technology for long, secure wear
- Ideal for long term monitoring (can stay in place for 5 days)
- The high adhesion foam backing makes it ideally suited for stress testing
- Size 4cm x 3.5cm

Pkt/50 \$20.05  
 10+ \$18.15  
 3M2560



### Red Dot™ Soft Cloth Monitoring Electrode

- Multipurpose monitoring electrode with soft and breathable cloth backing
- Gentle on frail and elderly skin
- Solid gel electrode
- Size 6cm diameter

Pkt/50 \$22.30  
 10+ \$20.55  
 3M2238



### Red Dot™ Micropore™ Monitoring Electrodes

A round shaped monitoring electrode with soft and breathable micropore backing that conforms to the patient's skin.

- Solid gel electrode
- Suggested wear time up to 4 days
- Size 6cm diameter

Pkt/50 \$22.30  
 10+ \$20.55  
 3M2239



### Red Dot™ Resting ECG Electrode 2330

Designed to meet the unique needs of Diagnostic ECG monitoring. A semi-stiff tab style electrode with a water based conductive adhesive that is suitable for resting diagnostic testing.

- 2cm x 2cm Resting ECG Electrode
- Use with Alligator/Universal Clips
- Gentle to skin

Pkt/100 \$20.00  
 10+ \$19.25  
 3M2330



### Red Dot™ Trace Prep

To abrade the skin prior to applying electrode to improve trace quality

18mm x 5m  
 3M2236 \$20.35



### Red Dot™ Resting ECG Electrode 2360


Same features as #2330 above. 2360 electrode has a stronger adhesive than the 2330 electrode. Use for added adhesion on challenging sites

Pkt/100 \$20.80  
 10+ \$19.15  
 3M2360



### ECG universal clip

- Locks on to tab electrode
- Snap and lock overstud
- Fits 3 and 4mm pins
- Holds electrodes flat
- High quality tracing



Pkt/12 \$113.90 BIOECGCLIP



## ECG Electrodes


### KENDALL™ CA610 DIAGNOSTIC ECG TAB

Offers a distinctive shape and a broader tab. The clean vinyl substrate permits easy detection of any potential skin reaction.

23 x 32mm

Pkt/100 **\$15.10**  
10+ **\$12.90**

[COV31447793](#)



### KENDALL™ GOLD 5500 DIAGNOSTIC ECG TAB

The 5500 tab electrode features our highest tack adhesive gel providing maximum adhesion.

22mm x 31mm

Pkt/100 **\$20.45**  
20+ **\$17.45**

[COV30807732](#)



### Medi-Trace™ 450 Foam Series

Thinner conformable foam substrate. Rectangle, stud, general monitoring electrode, 35mm x 43mm

Pkt/50 **\$35.45**

[CAR22450](#)



### Medi-Trace™ 200 Foam Series

Round, stud, general monitoring electrode, Ø 38mm

Pkt/100 **\$29.95**  
10+ **\$27.45**

[COV31050522](#)



#### EDAN

SE601	110mm Z-fold	Pkt/140	<a href="#">ECGPAP125</a>	<b>\$21.00</b>
SE1200	210mm Z-fold	Pkt/150	<a href="#">ECGPAP126</a>	<b>\$30.45</b>
SE1200	210mm	Roll	<a href="#">ECGPAP127</a>	<b>\$14.60</b>

#### Cardioline AR series

AR600	60mm	Roll	<a href="#">ECGPAP160</a>	<b>\$9.10</b>
AR1200	120mm Z-fold	Pkt/300	<a href="#">ECGPAP130</a>	<b>\$17.20</b>
AR2100	210mm	Roll	<a href="#">ECGPAP140</a>	<b>\$33.50</b>
AR2100	210mm Z-fold	Pkt/200	<a href="#">ECGPAP150</a>	<b>\$53.90</b>

#### Cardioline DELTA series

DELTA 1	60mm	Roll	<a href="#">ECGPAP180</a>	<b>\$6.40</b>
DELTA 3	112mm Z-fold	Pkt/300	<a href="#">ECGPAP190</a>	<b>\$14.00</b>
DELTA 60	210mm	Each	<a href="#">ECGPAP200</a>	<b>\$35.60</b>

#### Fukuda

FUKUDA	63mm x 30mm	Each	<a href="#">ECGPAP260</a>	<b>\$7.50</b>
FUKUDA	145mm x 75mm	Roll/5	<a href="#">ECGPAP250</a>	<b>\$20.00</b>

#### Mortara

ELI50	50 x 70mm	Each	<a href="#">ECGPAP280</a>	<b>\$10.10</b>
ELI100	108mm x 25mm	Roll	<a href="#">ECGPAP270</a>	<b>\$15.40</b>
ELI250	210 x 300mm	Pkt/250	<a href="#">ECGPAP300</a>	<b>\$31.75</b>

#### Schiller

AT-1	90mm	Pkt/400	<a href="#">ECGPAP290</a>	<b>\$16.15</b>
AT-2	210mm	Pkt/200	<a href="#">ECGPAP291</a>	<b>\$34.85</b>

#### Welch Allyn

CP100/CP200	210 x 280mm (generic) Z-fold	Pkt/200	<a href="#">ECGPAP292</a>	<b>\$40.40</b>
CP100/CP200	210 x 280mm (WA Original) Z-fold	Pkt/200	<a href="#">WA105353P</a>	<b>\$74.15</b>
CP50	114mm Z-fold	Box/4	<a href="#">ECGPAP325</a>	<b>\$125.80</b>
CP50	114mm	Pkt/4	<a href="#">WA406020</a>	<b>\$86.45</b>

## ECG PAPER

