

HEINE XP® Emergency Set with XP Disposable Blades

Specially developed for emergency use featuring the HEINE XP® disposable laryngoscope blades. These sets are complete with the 5 most frequently used blade types and sizes for infants and adults and the world's best-selling F.O. blade with integrated fiber optic bundle. The set now includes a disposable Handle Shell for the F.O Handle. See below.

\$330 HEF25710815 RRP \$326

HEINE QUALITY
MADE IN GERMANY



JOIN THE TEAM MEDICAL FAMILY TO ACCESS YOUR SHOP, YOUR PRICES, YOUR ORDER HISTORY, AND EXCLUSIVE OFFERS!



Exceptional Customer Service



Access to over 15,000 products



Industry leading online ordering system



Highly experienced Sales Team



Multiple regional warehouses



Fast & reliable delivery



Competitive pricing on quality products

OPEN AN ACCOUNT TODAY - HEAD TO WWW.TEAMMED.COM.AU/REGISTER

Disposable Laryngoscope set

This disposable plastic laryngoscope set is manufactured with high quality polycarbonate. It consists of a AA handle, three Macintosh blades, 3 Miller blades and a plastic storage case.

A strong single use item designed for field use. Replacement items are available for individual sale.



\$100.00* ADTPN1151

Laryngoscope set

Blades and Handle in protective hard case

McIntosh Set:

- Blades 1 to 4
- 32mm handle
- Stainless steel

\$345.00

BYDAD0186
Standard illumination

\$545.00

BYDAD0201
Fibre Optic illumination



Endotracheal tubes Cuffed

Sz 4.0mm	MMDAN021002
Sz 4.5mm	MMDAN021003
Sz 5.0mm	MMDAN021004
Sz 5.5mm	MMDAN021005
Sz 6.0mm	MMDAN021006
Sz 6.5mm	MMDAN021007
Sz 7.0mm	MMDAN021008
Sz 7.5mm	MMDAN021009
Sz 8.0mm	MMDAN021010
Sz 8.5mm	MMDAN021011
Sz 9.0mm	MMDAN021012



Each **\$7.15*** 10+ **\$5.70***

Endotracheal tubes Uncuffed

Sz 2.5mm	SMI100141025
Sz 3.0mm	SMI100141030
Sz 3.5mm	SMI100141035
Sz 4.0mm	SMI100141040
Sz 4.5mm	SMI100141045
Sz 5.0mm	SMI100141050
Sz 5.5mm	SMI100141055
Sz 6.0mm	SMI100141060
Sz 6.5mm	SMI100141065
Sz 7.0mm	SMI100141070
Sz 7.5mm	SMI100141075



Each **\$13.65*** 10+ **\$13.05***

Yeescop

m|devices

The yeescop is a one-piece disposable laryngoscope incorporating handle and blade for single patient use.

- Function of the light can be tested in the sealed original pack
- Single use - no cross infection
- Powerful illumination



Mac 3 **MMDAN010000**
Mac 4 **MMDAN010001**

Each **\$28.00*** 10+ **\$26.20***

Silicone Disposable Laryngeal Airway Mask

		Patient Weight	
Sz 1	Black	(0-5kg)	MMDAN030000
Sz 1.5	White	(5-10kg)	MMDAN030001
Sz 2	Red	(10-20kg)	MMDAN030002
Sz 2.5	Orange	(20-30kg)	MMDAN030003

		Patient Weight	
Sz 3	Green	(30-50kg)	MMDAN030004
Sz 4	Yellow	(50-70kg)	MMDAN030005
Sz 5	Pink	(70-100kg)	MMDAN030006

Each **\$20.70*** 10+ **\$19.40***

Each **\$20.70*** 10+ **\$19.40***



m|devices

Guedel Airway Kit

The Guedel Airway Kit consists of 8 Guedel airways (one of each size) supplied in a plastic case. Ideal for medical kits and paramedic cases.

\$26.75*

ADTPN0523E



Guedel Airways Non-Sterile

Sz 00	50mm	blue	MIASLGA00	Each	\$1.55*
Sz 0	60mm	black	MIASLGA0	10+	\$1.30*
Sz 1	70mm	clear	MIASLGA1		
Sz 2	80mm	green	MIASLGA2N		
Sz 3	90mm	yellow	MIASLGA3N		
Sz 4	100mm	red	MIASLGA4N		
Sz 5	110mm	orange	MIASLGA5N		



Guedel Airways Sterile

Sz 00	40mm	pink	COV2877838	\$6.90*	10+	\$5.90*
Sz 0	50mm	blue	COV2877636	\$4.75*	10+	\$4.10*
Sz 1	60mm	black	COV2877635	\$4.25*	10+	\$3.60*
Sz 2	70mm	clear	COV2877634	\$4.25*	10+	\$3.60*
Sz 3	80mm	green	COV2877633	\$4.25*	10+	\$3.60*
Sz 4	90mm	yellow	COV2877632	\$4.25*	10+	\$3.60*
Sz 5	100mm	red	COV2877631	\$4.75*	10+	\$4.10*
Sz 6	110mm	orange	COV2877843	\$5.00*	10+	\$4.30*


i-gel® Supraglottic Airway

The first major development since the laryngeal mask airway, i-gel has changed the face of airway management and is now widely used in anaesthesia and resuscitation across the globe.

Made from a medical grade thermoplastic elastomer, i-gel has been designed to create a non-inflatable, anatomical seal of the pharyngeal, laryngeal and perilaryngeal structures whilst avoiding compression trauma.

- i-gel has a soft, gel-like, non-inflatable cuff, designed to provide an anatomical impression fit over the laryngeal inlet.
- The shape, softness and contours accurately mirror the perilaryngeal anatomy - an innovative concept meaning no cuff inflation is required.
- Working in harmony with the patient's anatomy, compression and displacement trauma are significantly reduced or eliminated.

i-gel®



Size 1	Neonate (2-5kg)	INT8201000
Size 1.5	Infant (5-12kg)	INT8215000
Size 2	Small Child (10-25kg)	INT8202000
Size 2.5	Large Child (25-35kg)	INT8225000
Size 3	Small Adult (30-60kg)	INT8203000
Size 4	Med Adult (50-90kg)	INT8204000
Size 5	Large Adult (90+kg)	INT8205000



Each **\$30.45**
10+ **\$29.10**

i-gel O2 Resus Pack

Includes an i-gel O2 airway with green hook ring, sachet of lubricant, airway support strap and a 12FG suction tube

Small Adult (30-60kg) - with sz 3 i-gel airway	INT8703000
Med Adult (50-90kg)- with sz 4 i-gel airway	INT8704000
Large Adult (90+kg)- with sz 5 i-gel airway	INT8705000

Each **\$45.15**



Oxygen Accessories

m|devices

Oxygen mask

adult [MMD061002NS](#) \$2.25* 50+ \$2.50*

Oxygen mask + 2.1m Tube

adult [MMD061003](#) \$4.80* 50+ \$3.90*

child [MMD060003](#) \$4.45* 50+ \$3.55*

Oxygen Tubing 2.1m

[MMD130000](#) \$2.50* 50+ \$2.15*

Nasal oxygen cannula with 2.1m tubing

adult [MMD070002](#) \$4.70* 50+ \$3.65*

child [MMD070001](#) \$4.90* 50+ \$3.85*

infant [MMD070000](#) \$3.95* 50+ \$3.10*



Hi-Oxygen Concentration Mask

- Soft, transparent, odour free vinyl
- Aluminium fixer and elastic strip to ensure the best fit and user comfort
- One way valve and reservoir to increase the efficiency of oxygen concentration

adult [LRESA0301](#)
child [LRESA0302](#)

\$6.35*



Keyring Resuscitation Masks

Be prepared for an emergency situation anytime, anywhere with this convenient one-way breathing rescue mask.

\$8.25
[FERRMAGK](#)



CPR Pocket Mask

Includes mask, storage case and filter. The PN-1141 is ideal for resus and paramedics packs.

\$15.00*
[ADTPN1141](#)



Disposable resuscitator/anaesthesia mask

Soft PVC Cushioned Mask (single use)

m|devices



Sz 1	Neonate	MMDAN140001NS
Sz 2	Infant	MMDAN140002NS
Sz 3	Child	MMDAN140003NS
Sz 4	Small Adult	MMDAN140004NS
Sz 5	Adult	MMDAN140005NS


\$7.50*

Disposable Resuscitators

Includes bag, mask, patient & intake valves, reservoir bag and 2 metres oxygen tubing.

Infant [ADTRE24610](#)
Child [ADTRE24510](#)
Adult [ADTBE2501](#)

\$36.75*



Reusable Resuscitators

The fully autoclavable Ferno Reusable Resuscitators offer great value and exceptional quality. They are manufactured to the highest standards using 100% latex-free, transparent silicone.

Includes:

- Bag
 - Mask
 - Patient and intake valves
 - Reservoir bag
 - 2m of O₂ drive-line
- | | |
|--------|---------------------------|
| Infant | FERRBMFIR |
| Child | FERRBMFCR |
| Adult | FERRBMFAR |

\$225.00*



Adult Silicone Resuscitation Kit

This Latex-free Resuscitation Kit is a great addition to your clinic. All devices are in a handy tool box for easy storage.

Contents:

- BE2101 Resuscitator
- size 4 & 5 mask
- Size 2, 3 & 4 Guedel airways
- 2500ml reservoir bag
- Oxygen tubing
- Durable hard carry case



ADTBE2100KIT

\$195.00*

Oxygen Therapy Kit with Resus Bag

Suitable for emergency resuscitation and oxygen therapy

- 1 x Oxygen therapy bag
- 1 x Oxygen therapy regulator
- 1 x Single use adult manual resuscitator set with 60cm H2O pop off
- 1 x CPR pocket mask
- 1 x Guedel airway kit
- 2 x Oxygen mask adult disposable with 7ft tube
- 2 x Oxygen mask child disposable with 7ft tube
- 2 x Nasal canula adult disposable
- 2 x Nasal canula child disposable

ADTOTK12026

\$695.00*
(bag included)



Spiracle BRAVO2 Oxygen Regulator

- Brass regulator with 11 variable flow rates from 0.25 to 25 LPM
- Safe and effective
- Made from quality crafted internal brass components
- Nickel plating on the regulator provides an attractive corrosion resistant finish
- Fits on all cylinders
- Meets Australian Standards



\$440.00* FEROXRS833

Piston Therapy Oxygen Regulator

The Piston Oxygen Therapy Regulator is a reliable oxygen therapy delivery device. It offers a Heavy Duty and compact design allowing it to be used in most therapy applications. Full Brass body make this regulator a durable and reliable solution to any environment where oxygen therapy is needed.

- Full brass body and internals
- Barb outlet
- 11 therapy flow settings (0.5, 1, 2, 3, 4, 5, 6, 8, 10, 12, 15, LPM)
- Weighs approximately 700g

\$188.75 ADTME870



Oxygen therapy regulator

The Series-O regulator has been the industry leader in Australia for years. The simple but effective design and ease of use make this single stage regulator a very low maintenance yet high performing product that can be used in a variety of environments.

- Triple filter system increases the regulator performance by preventing foreign particles from getting in contact with the high pressure gas or the regulator seat
- No tools are required to routinely disconnect the regulator from the cylinder valve
- A safety relief valve is built in to the regulator to protect staff and patients in the event of high pressure failures
- Gas specific inlet fittings and outlet fittings to prevent connecting to the wrong gas type
- Pin indexed yoke inlet
- 400kPa max pressure outlet
- 20,000kPa @ 15°C max rated inlet pressure



ADT518804 \$540.00*

Simpflow Oxygen Flowmeter 15LPM

- Provides an accuracy of +/-5% (of full scale reading) and is pressure compensated to ensure the flowmeter maintains its accuracy independent of varying gas supply
- Flow limited output orifice and fail safe shroud protects the patient from excessive flows in the event of excessive gas pressures
- Gas specific inlet fitting to prevent connection to wrong gas type
- Colour coded inlet fitting, outlet fitting and control knobs for quick recognition of gas used by the flowmeter
- The flowmeter design provides clear readability with large bold print lines and numbers providing 180° visibility



\$205*

ADTX011020XT

Twin-O-Vac multi-purpose device

The Twin-O-Vac is a low cost versatile unit that is ideal when extra or portable sources of oxygen or suction are required. Each Twin-O-Vac uses one oxygen outlet (i.e. an oxygen pipeline wall point or cylinder regulator) to provide variable suction plus two separate sources of oxygen.

All three outlets can be used simultaneously on up to three different patients if necessary. The Twin-O-Vac can be connected or linked to an assortment of equipment including resuscitators and flowmeters.




\$1030*

ADTTM117G

Do you know what to use and when?

The below chart will show you how and when to set up your oxygen kit.



Broselow™ Paediatric Emergency Tape

The Broselow Tape is recognised in most medical textbooks and publications as a standard for the emergency treatment of children.

The tape is calibrated for paediatric patients between 46cm and 146.5cm in length, which corresponds to a weight range of 3kg to 34kg. The tape is divided into nine colour zones with pre-calculated drug doses, fluid volumes, and equipment sizes for each zone.

DAL2132

\$172.10

